

# BIOFIRE® Joint Infection Panel Targets

## GRAM-POSITIVE BACTERIA

*Anaerococcus prevotii/vaginalis*  
*Clostridium perfringens*  
*Cutibacterium avidum/granulosum*  
*Enterococcus faecalis*  
*Enterococcus faecium*  
*Finnegoldia magna*  
*Parvimonas micra*  
*Peptoniphilus*  
*Peptostreptococcus anaerobius*  
*Staphylococcus aureus*  
*Staphylococcus lugdunensis*  
*Streptococcus* spp.  
*Streptococcus agalactiae*  
*Streptococcus pneumoniae*  
*Streptococcus pyogenes*

## GRAM-NEGATIVE BACTERIA

*Bacteroides fragilis*  
*Citrobacter*  
*Enterobacter cloacae* complex  
*Escherichia coli*  
*Haemophilus influenzae*  
*Kingella kingae*  
*Klebsiella aerogenes*  
*Klebsiella pneumoniae* group  
*Morganella morganii*  
*Neisseria gonorrhoeae*  
*Proteus* spp.  
*Pseudomonas aeruginosa*  
*Salmonella* spp.  
*Serratia marcescens*

## YEAST

*Candida* spp.  
*Candida albicans*

## ANTIMICROBIAL RESISTANCE GENES

### Carbapenemases

IMP  
KPC  
NDM  
OXA-48-like  
VIM

### ESBL

CTX-M

### Methicillin Resistance

*mecA/C* and MREJ

### Vancomycin Resistance

*vanA/B*

39  
TARGETS  
~1hr

FDA cleared | CE<sub>2797</sub>

## Sample Requirements

0.2mL of synovial fluid

## Overall Performance

- 91.7% Sensitivity<sup>20</sup>
- 99.8% Specificity<sup>20</sup>

# Guidelines

- European Bone and Joint Infection Society. Guideline for the management of septic arthritis in native joints (SANJO). Ravn C., et al, J. Bone Joint Infect, 8, 29–37
- 2025 International Consensus Meeting: Diagnostic Techniques: Molecular Tests. Martinazzi, Brandon J. *et al*, J Arthroplasty
- American Academy of Orthopaedic Surgeons Evidence-Based Clinical Practice Guideline for Diagnosis and Prevention of Periprosthetic Joint Infections <https://www.aaos.org/pjicpg>

# References

1. Premkumar A, *et al.* (2021) J Arthroplasty. 36(5) 1484-1489.e3
2. Singh JA, Yu S. (2017) PLoS One. 12(8): e0182577
3. Carpenter CR, *et al.* (2011) Acad Emerg Med. 18(8): 781–796
4. Tande JT, Patel R, (2014) Clin Microbiol. 27(2): 302–345
5. Zimmerli W, *et al.* (2004) N Engl J Med. (16):1645-54
6. Parvizi J (2018) J Arthroplasty. 33(5):1309-131
7. Ghirardelli S, *et al.* (2024) J Exp Orthop. 11: e12097
8. Berinson B, *et al.* (2025) Microbiol Spectr. 10.1128/spectrum.01263-25
9. Gaillard T, *et al.* (2024) Clin Microbiol Infect. 7:905-910
10. Pascual S, *et al.* (2024) J. Bone Joint Infect. 9: 87–97
11. Kesavan P, *et al.* (2024) Infect Dis Clin Pract. 32(6)
12. Agosta M, *et al.* (2025) Diagn. 15: 566
13. Copley L, *et al.* (2025) J Clin Microbiol. 63:e00278-25
14. Lin TL, *et al.* (2025) Microbiol Spectr. 13:e00286-25
15. Esteban J, *et al.* (2023) 61(11):e0035723
16. McNally M, *et al.* (2021) Bone Joint J. 103-B(1):18-25
17. Bauer TW, *et al.* (2019) J Arthroplasty. 34:S351–S359
18. Parvizi J, *et al.* J Arthroplasty. 2018; 33(5):1309-1314.e2
19. Osmon DR, *et al.* Clin Infect Dis. 2013 56:e1–e25.
20. Overall performance based on prospective clinical study for the BIOFIRE® Joint Infection Panel, data on file, BIOFIRE Diagnostics.

Product availability varies by country. Consult your bioMérieux representative.

# Contact Us

bioMérieux S.A.  
69280 Marcy l'Etoile  
France  
Tel.: +33 (0) 4 78 87 20 00  
Fax: +33 (0) 4 78 87 20 90  
[biomerieux.com](http://biomerieux.com)

Learn more about the BIOFIRE range of commercially-available panels for syndromic infectious disease diagnostics.





# Clinical Impact of the BIOFIRE® Joint Infection (JI) Panel

39

TARGETS

~1<sup>hr</sup>

# What's the Problem?

Joint infections cause a tremendous burden for patients and society.<sup>1,2</sup>

**Septic arthritis is a medical emergency requiring prompt diagnosis and treatment.**

Delayed diagnosis is associated with permanent disability and increased mortality, which can be as high as 15%.<sup>3</sup>

**Prosthetic joint infections (PJIs) are costly to treat and on the rise.**<sup>1</sup>

The cost to treat a PJI is 3 to 6 times more expensive than the initial arthroplasty.<sup>4</sup> When missed or undertreated, PJIs can lead to unnecessary surgical revisions causing poor function or disability, considerably impacting quality of life.<sup>5</sup>

---

## Diagnosis of Joint Infections is Complicated

- Joint infection diagnostics lack standardization of specimen type and preparation, test media, and methods.<sup>4</sup>
- Culture negative PJIs occur in up to 35% of infections.<sup>6</sup>
- Joint infections are associated with difficult fastidious organisms, anaerobes, biofilm-forming organisms, and polymicrobial specimens.<sup>4</sup>
- Complex society-developed diagnostic criteria vary considerably in diagnostic agreement.<sup>6</sup>



# The Right Test, The First Time

BIOFIRE's syndromic approach combines several potential targets into one fast test, helping clinicians identify pathogens that produce non-specific symptoms like red, hot, swollen joint(s) in a clinically actionable timeframe.

## Syndromic Testing




## Faster Than Traditional Methods

Traditional methods require multiple tests and can take up to two weeks to provide a pathogen identification result. The BIOFIRE® Joint Infection (JI) Panel can reduce time to pathogen ID by up to 4.4 days.<sup>7</sup>

## Improved Diagnostic Yield

The BIOFIRE JI Panel demonstrated increased diagnostic yield for on-panel targets,<sup>7-13</sup> including difficult-to-grow anaerobes & fastidious organisms,<sup>10-11</sup> and patients with antibiotic exposure.<sup>9,12</sup> The BIOFIRE JI Panel increased diagnostic yield by 16.7% for on-panel targets.<sup>8</sup>



**16.7%**  
Increased diagnostic  
yield for on-panel  
targets<sup>8</sup>

## Polymicrobial Detections

The BIOFIRE JI Panel prospective clinical trial demonstrated polymicrobial detections. Of the 242 positive specimens detected by the BIOFIRE JI Panel, 16 involved co-detections.<sup>15</sup>

## Pathogen Guided-Patient Management

Pathogen identification is a central component of septic arthritis and PJI treatment guidelines.<sup>16-19</sup> The BIOFIRE JI Panel provides clinicians with the results they need to make pathogen-guided decisions for the management of patients with native joint septic arthritis (SA) and prosthetic joint infections (PJI).



**Fast, accurate pathogen identification can aid clinicians in antibiotic therapy decisions, including appropriate escalation, de-escalation, and administration route of antibiotics.<sup>8</sup>**



**Antibiotic therapy was modified in 100% of patients with positive JI Panel results and previous antibiotic exposure.<sup>12,\*</sup>**



**The JI Panel reduced the time on empiric antibiotics by 66.4 hours.<sup>14</sup>**

\*N=9 patients with previous antibiotic exposure and positive BIOFIRE JI Panel results.

## Informed Surgical Decision Making

With the BIOFIRE JI Panel, 19.2% of PJI patients avoided 2-stage revisions and underwent single-stage or 1.5-stage revisions due to BIOFIRE JI Panel results.<sup>7</sup>