



# DIAGNOSIS OF MYCOBACTERIA

Fluorescence microscopy observation with Fluo-RAL *Methylene Blue*



HORIZONTAL SCANNING



VERTICAL SCANNING






NB: use appropriate microscope for LED\*-based Fluorescence Microscopy (WHO recommendations)

(\*) Light Emitting Diode

(\*\*) Laboratory diagnosis of tuberculosis by Sputum Microscopy. The hand book. A publication of the global laboratory initiative, a working group of the stopTB partnership. Global edition. , 2013

(\*\*\*) AFB: Acid-Fast Bacilli

(\*\*\*\*) Confirmation required by another technician, or prepare another smear, stain and read.

Result to report The Handbook/ Stop TB Partnership(**)	Observed 200-250 x magnification 1 length = 30 fields	Observed 400 x magnification 1 length = 40 to 60 fields	Reading
No AFB(***) observed	No AFB / 1 length	No AFB / 1 length	
Confirmation required(****)	1-4 AFB / 1 length	1-2 AFB / 1 length	
Scanty	5-49 AFB / 1 length	3-24 AFB / 1 length	
1+	3-24 AFB / 1 field	1-6 AFB / 1 field	
2+	25-250 AFB / 1 field	7-60 AFB / 1 field	
3+	>250 AFB / 1 field	>60 AFB / 1 field	



# DIAGNOSIS OF MYCOBACTERIA

Fluorescence microscopy observation with Fluo-RAL *Methylene Blue*

## TROUBLESHOOTING

Problem	Possible cause	Remedy
Artefacts	<ul style="list-style-type: none"> <li>• Sample integrity</li> <li>• Glass slide scratched</li> </ul>	<ul style="list-style-type: none"> <li>• Transport &amp; store sample according to WHO &amp;/or local guidelines</li> <li>• Learn to recognize glass artefacts</li> </ul>
Unreadable slide	<ul style="list-style-type: none"> <li>• Poor quality or insufficient amount of specimen</li> </ul>	<ul style="list-style-type: none"> <li>• New sample</li> </ul>
Positive slide reported negative	<ul style="list-style-type: none"> <li>• Incomplete slide reading</li> <li>• Slide exposed to daylight for too long</li> <li>• Smear too thick, washed off during staining/rising</li> </ul>	<ul style="list-style-type: none"> <li>• Read suggested number of fields</li> <li>• Store slides away from light</li> <li>• Prepare new smear less thick and respect the protocol of fixation &gt;20 min 80°C on heating plate</li> </ul>
Low positive or Pale AFB	<ul style="list-style-type: none"> <li>• Too long an interval between staining and reading</li> <li>• Stains kept in the light</li> <li>• Outdated kit lifespan</li> <li>• Auramine staining time less than 10 min</li> <li>• Overdecolourised</li> </ul>	<ul style="list-style-type: none"> <li>• Keep smears covered and examine within 24 hours of staining</li> <li>• Store away from light</li> <li>• Change kit</li> <li>• Stain for a minimum of 10 min</li> <li>• Do not exceed the maximum time (1 min)</li> <li>• Increase discolouration time</li> </ul>
Negative slide reported positive	<ul style="list-style-type: none"> <li>• Insufficient discolouration time</li> <li>• Small artefacts mistaken for AFB</li> <li>• Cross contamination</li> <li>• Smear size too large</li> </ul>	<ul style="list-style-type: none"> <li>• Cf Artefacts</li> <li>• Respect fixation protocol &gt; 20 min 80 °C on heating plate + fixative solution time for a minimum 5 min</li> <li>• WHO recommendation: 1 cm x 2 cm-sized smears</li> </ul>
QC negative spot is positive	<ul style="list-style-type: none"> <li>• Contaminated stain prepared with water containing environmental mycobacteria</li> <li>• Possible artefacts</li> </ul>	<ul style="list-style-type: none"> <li>• Use distilled water</li> <li>• Ready-to use QC slides recommended</li> </ul>
Background too dark	<ul style="list-style-type: none"> <li>• Counterstained too long (or) decolourised too long</li> <li>• Smear too thick</li> </ul>	<ul style="list-style-type: none"> <li>• Do not exceed 1 min for both steps</li> <li>• Prepare new smear</li> </ul>
Too much fluorescence	<ul style="list-style-type: none"> <li>• Insufficient decolourisation</li> <li>• Smear too thick</li> </ul>	<ul style="list-style-type: none"> <li>• Check decolourisation time</li> <li>• Prepare new smear</li> </ul>