

BIOFIRE® SPOTFIRE® Respiratory/Sore Throat (R/ST) Panel



Instructions for Use	https://www.biofiredx.com/e-labeling/ITI0201
Quick Guide	https://www.biofiredx.com/e-labeling/ITI0213
Safety Data Sheet (SDS)	https://www.biofiredx.com/e-labeling/ITI0225
Customer and Technical Support Information	Phone: 1-800-682-2666 (toll free) E-mail: BioFireSupport@bioMerieux.com Website: www.biofiredx.com
*For more information on how to contact Customer and Technical Support, refer to Appendix B.	Or contact the local bioMérieux sales representative or an authorized distributor.

IVD

CE₂₇₉₇

INTENDED PURPOSE

Intended Use

The BIOFIRE® SPOTFIRE® Respiratory/Sore Throat (R/ST) Panel (SPOTFIRE R/ST Panel) test kit is a single product with **two intended uses**, where the intended use for a given test is based upon the sample type selected by the operator based on patient signs and symptoms and sample type collected.

The SPOTFIRE R/ST Panel is a multiplexed polymerase chain reaction (PCR) test intended for use with the BIOFIRE® SPOTFIRE® System for the simultaneous, qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swab (NPS) specimens obtained from **individuals with signs and symptoms of respiratory tract infections, including COVID-19; (Respiratory)** or from **individuals with signs and/or symptoms of pharyngitis** (using a throat swab (TS); **Sore Throat**).

The following organism types and subtypes are identified and differentiated using the SPOTFIRE R/ST Panel:

Viruses (Respiratory and Sore Throat)	Bacteria (Respiratory and Sore Throat)
Adenovirus	<i>Chlamydia pneumoniae</i>
Coronavirus (seasonal)	<i>Mycoplasma pneumoniae</i>
Coronavirus SARS-CoV-2	Bacteria (Respiratory Only)
Human metapneumovirus	<i>Bordetella parapertussis</i>
Human rhinovirus/enterovirus	<i>Bordetella pertussis</i>
Influenza A virus	Bacteria (Sore Throat Only)
Influenza A virus A/ H1-2009	<i>Streptococcus dysgalactiae</i> (group C/G Strep)
Influenza A virus A/H3	<i>Streptococcus pyogenes</i> (group A Strep)
Influenza B virus	
Parainfluenza virus	
Respiratory syncytial virus	

Nucleic acids from the viral and bacterial organisms identified by this test are generally detectable in NPS/TS specimens during the acute phase of infection. The detection and identification of specific viral and bacterial nucleic acids from individuals exhibiting signs and/or symptoms of respiratory infection and/or pharyngitis are indicative of the presence of the identified microorganism and aids in diagnosis if used in conjunction with other clinical and epidemiological information, and laboratory findings. The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

Negative results in the setting of a respiratory illness and/or pharyngitis may be due to infection with pathogens that are not detected by this test, or lower respiratory tract infection that may not be detected by an NPS or TS specimen. Positive results do not rule out coinfection with other organisms. The agent(s) detected by the SPOTFIRE R/ST Panel may not be the definite cause of disease.

Additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence, and radiography) may be necessary when evaluating a patient with possible respiratory tract infection and/or pharyngitis.

Intended User and Use Environment

The SPOTFIRE R/ST Panel is intended for use by trained medical professionals proficient in using the SPOTFIRE System at the point of care (POC) or for laboratory professionals in a laboratory setting.

Refer to the SPOTFIRE System Operator's Manual for the appropriate physical environmental specifications and/or conditions for performing a SPOTFIRE R/ST Panel test.

SUMMARY AND EXPLANATION OF THE TEST

The SPOTFIRE R/ST Panel, designed for use with the SPOTFIRE System, is a polymerase chain reaction (PCR) based sample-to-answer diagnostic test that simultaneously identifies nucleic acids from 15 different bacterial and viral organisms from nasopharyngeal swab (NPS) specimens, or 15 different bacterial and viral organisms from throat swab (TS) specimens, in transport media collected from individuals with signs and symptoms of respiratory infection or pharyngitis, respectively. The SPOTFIRE R/ST Panel uses a single instrument protocol with different reporting of analytes for the two sample types. Sample type is selected at the time of testing and the system's software controls the analyte reporting based on the selected sample type.

Summary of Detected Organisms

Adenoviruses (AdV) are a diverse group of non-enveloped DNA viruses in the family *Adenoviridae* with seven species (A to G).¹ Adenovirus species B, C, and E cause acute respiratory disease, but all types have been associated with human disease.² Other Adenovirus species (A, D, F and G) can cause a variety of illnesses, including cystitis, gastroenteritis, and conjunctivitis³, and may also be found in respiratory specimens. Outbreaks often occur in institutional settings such as military training, long-term care facilities, and pediatric tertiary-care hospitals, due to high rates of transmission in closed populations.⁴⁻⁶ Adenoviruses are shed for long periods of time and persist on surfaces in an infective state.⁶

Coronaviruses (CoV) are related RNA viruses in the family *Coronaviridae*. Human coronaviruses were established as respiratory pathogens in the 1960s and seven serological variants associated with human disease have been characterized to date: four types (**coronaviruses 229E, OC43, HKU1, NL63**) that regularly circulate in human populations, and are a cause of common colds^{7,8}, and three strains (Middle East Respiratory Syndrome Coronavirus (MERS-CoV), Severe Acute Respiratory Syndrome Coronavirus (SARS-CoV), and **Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2)**) that have spread from animal to human populations since 2002⁸ and represent the ongoing public health threat posed by emerging zoonotic pathogens. SARS-CoV was declared contained by the World Health Organization (WHO) in 2003, less than 12 months after its emergence, and no new cases have been reported since 2004. MERS-CoV, which was first described in 2012,⁹ continues to cause occasional outbreaks, which are characterized by animal to human transmissions followed by person-to-person transmission.¹⁰ SARS-CoV-2 is the novel coronavirus that causes COVID-19, an illness that reached a level of pandemic spread in the short time since its emergence in late 2019.^{11,12} Coronaviruses have been linked to croup and exacerbation of asthma.^{13,14} Infections with coronavirus 229E, OC43, HKU1, and NL63 occur more often in the winter, and there appears to be a periodicity of circulation.¹⁵ Illnesses caused by these coronaviruses are generally self-limiting.¹⁶ Though coronaviruses as a group are most commonly associated with upper respiratory tract infections, all human coronaviruses are also associated with lower respiratory tract infection and MERS-CoV, SARS-CoV, and SARS-CoV-2 can cause Acute Respiratory Distress Syndrome (ARDS), as well as significant rates of hospitalization, complications, and death, especially in patients with underlying health conditions.¹⁷

Note: MERS-CoV and SARS-CoV are not detected by the SPOTFIRE R/ST Panel.

Human metapneumovirus (hMPV) is an RNA virus in the family *Paramyxoviridae*.¹⁸ HMPV was discovered in 2001 as a respiratory pathogen in children.¹⁹ Further studies confirmed hMPV infections in persons of all ages.²⁰ The two genotypes, A and B, can circulate at the same time and do not appear to differ in the severity of illness.¹⁸ HMPV is the second leading cause of bronchiolitis in young children.¹⁸ Additionally, infection can result in a broad range of upper and lower respiratory symptoms: cough, rhinorrhea, wheeze, dyspnea, and fever.²¹ HMPV is estimated to be responsible for 5-7% of respiratory tract infections in children and 3% among individuals of all ages.²¹ The seasonal peak of hMPV is winter and early spring and often co-occurs with the seasonal peak of Respiratory syncytial virus (RSV).²²

Influenza A and B are RNA viruses in the *Orthomyxoviridae* family. The dominant type of influenza virus varies often due to antigenic drift.²³ Influenza A can be subtyped by the hemagglutinin (H) and neuraminidase (N) genes; influenza A subtypes **H1N1** and **H3N2** are the strains that most commonly infect humans. More severe disease and increased mortality

are associated with H3N2 subtype.²⁴ During the 2009-10 Influenza season, influenza A (H1N1) pdm09 (H1-2009, also known as “swine flu”) became the dominant circulating influenza virus, accounting for approximately 99% of reported influenza infections and has since replaced pre-2009 H1N1 strains (Table 1).²⁵ Currently, at least four antiviral medications are available for influenza treatment – amantadine, rimantadine, zanamivir and oseltamivir – with type-specific efficacy and drug resistance arising with the spread of new strains of the virus.²⁶ Complications with viral or bacterial pneumonia increase mortality from influenza infections.²⁷

Table 1. Proportions of Influenza Subtype Infections in the United States (as reported by the US Centers for Disease Control)

Flu Season	Influenza A	% of Subtyped Influenza A			Influenza B
		H1	H1-2009	H3 ²	
2023-2024 ³	82.7%	0.0	79.2	20.8	17.3%
2022-2023 ⁴	95.2%	0.0	30.3	69.7	4.8%
2021-2022	94.7%	0.0	12.7	99.9	5.3%
2020-2021 ⁴	63.6%	0.0	40.3	59.7 ⁴	36.4%
2019-2020	58.7%	0.0	92.8	7.2	41.3%

¹ CDC FluView data accessed on February 5, 2024.

² Includes H3N2 and H3N2v subtypes

³ Cumulative results through January 27, 2024.

⁴ Season during which SPOTFIRE R/ST Panel prospective clinical data described in this submission were accumulated.

Parainfluenza viruses (PIVs) are RNA viruses in the *Paramyxoviridae* family. In the 1950s, parainfluenza viruses were determined to be respiratory pathogens different from influenza viruses.²⁸ Parainfluenza viruses are divided into four types (**parainfluenza viruses 1, 2, 3, and 4**). Parainfluenza virus 1 causes biennial epidemics in the fall, with 50% of croup cases attributed to this virus.²⁸ Parainfluenza virus 2 causes epidemics every one to two years, which may alternate with parainfluenza virus 1 circulation.²⁸ Children less than six months old are particularly susceptible to parainfluenza virus 3 infection, with outbreaks occurring in neonatal intensive care units. Parainfluenza virus 3 is associated with the highest mortality and morbidity of all strains²⁹ and epidemics are most common in the spring and summer.²⁸ Parainfluenza virus 4 infection affects all age groups but because of infrequent detection periodicity of infection has not been established.^{30,31}

Respiratory syncytial virus (RSV) is an RNA virus in the *Paramyxoviridae* family and is related to human metapneumoviruses and parainfluenza viruses.³² RSV has two major subtypes (A and B), which vary annually in their prevalence.³³ RSV is the most common cause of severe respiratory disease in infants, with acute bronchiolitis as the major cause of hospitalization.³² RSV is now also recognized as an important pathogen in adults, although adult infections are in generally less severe and limited to the upper respiratory tract.³⁴ Peak activity of RSV is typically in January and February.³⁵

Rhinoviruses and **enteroviruses** are related RNA viruses in the *Picornaviridae* family.³⁶ There are more than 100 serotypes of human rhinovirus based on the serology of the capsid protein.³⁶ Rhinovirus is noted as causing the “common cold”, but may also be involved in precipitating asthma attacks and severe complications.³⁶ Enteroviruses are divided into four species that include a total of at least 89 distinct types. Individual types can be associated with different clinical manifestations, including nonspecific respiratory illnesses in infants or adults.³⁷ Both rhinoviruses and enterovirus are prevalent year round.^{38,39}

Bordetella pertussis is a gram-negative bacterium that is the predominant causative agent of whooping cough or pertussis, a vaccine-preventable, highly infectious disease that is reportable to public health organizations.^{40–42} Pertussis occurs most commonly in children but also occurs in adolescents and adults and outbreaks have been documented in fully vaccinated populations due to waning immunity (immunity has been shown to decrease 5-10 years after vaccination).^{42,43} Early (catarrhal) pertussis disease is non-specific, and classic signs of pertussis (paroxysmal coughing, inspiratory ‘whoop’, post-tussive emesis, as well as apnea or cyanosis in infants) do not arise until approximately two weeks after the initial onset of symptoms. No peak season has been defined for *B. pertussis*.

Bordetella parapertussis is a gram-negative bacterium that causes a pertussis-like disease that is similar to that caused by *B. pertussis* yet milder⁴⁴ and is characterized by early non-specific symptoms followed by classic signs of pertussis.⁴²

Differential diagnosis of *B. parapertussis* with *B. pertussis* is difficult due to the overlap in symptoms and the occurrence of concomitant infections.⁴² Generally less prevalent than *B. pertussis*, *B. parapertussis* is the causative agent in approximately 14% of whooping cough cases, predominantly in children.⁴⁵ Vaccination against pertussis does not protect against *B. parapertussis*. No peak season has been defined for *B. parapertussis*, but outbreaks of whooping cough (*B. parapertussis* and *B. pertussis*) appear to occur on a 3-5 year cycle.⁴⁶

Chlamydia pneumoniae (previously known as *Chlamydophila pneumoniae*) is a gram-negative obligate intracellular bacterium that causes acute respiratory infections and is a common cause of community-acquired atypical (walking) pneumonia and bronchitis.⁴⁷⁻⁴⁹ *C. pneumoniae* has an incubation period of approximately three weeks and can be transmitted from asymptomatic carriers.⁴⁹ Outbreaks occur in schools, military barracks, and nursing homes.⁵⁰ No peak season has been identified for *C. pneumoniae* infections.

Mycoplasma pneumoniae is a gram-negative bacterium that is another bacterial agent of community-acquired atypical pneumonia, occurring frequently in outbreak situations.^{51,52} Incubation time for *M. pneumoniae* infection is approximately 1 to 4 weeks.⁵³ *M. pneumoniae* respiratory disease does not have a defined season of highest incidence but epidemics have a periodicity of 3-7 years.⁵²

Streptococcus pyogenes or Group A Streptococcus (GAS) is a gram-positive beta-hemolytic streptococcus that causes acute pharyngitis, commonly known as strep throat. Patients experience sudden onset of sore throat, painful swallowing and fever. GAS causes pharyngitis in people of all ages but is most common in children and adolescents⁵⁴. GAS is the most common bacterial cause of pharyngitis in both children⁵⁵ and adults⁵⁶. GAS is spread person-to-person through saliva or nasal secretions and is most common in schools, daycare centers, and military training facilities. *S. pyogenes* infections occur most commonly in winter and spring⁵⁴.


Streptococcus dysgalactiae is a gram-positive beta-hemolytic streptococcus belonging to Lancefield groups C and G. *S. dysgalactiae* causes a variety of abscess forming infections in humans including pharyngitis as well as skin and soft tissue infections.^{57,58} For patients presenting with pharyngitis, the clinical signs and symptoms are insufficient in distinguishing between *S. dysgalactiae* and *S. pyogenes*.^{59,60} *S. dysgalactiae* is transmitted person-to-person and occurs in both community acquired outbreaks⁶¹ and clusters of hospital-acquired infections⁶². In contrast to *S. pyogenes*, there is no seasonality of pharyngitis caused by *S. dysgalactiae*.⁶³

PRINCIPLE OF THE PROCEDURE

The SPOTFIRE R/ST Panel pouch is a closed system disposable that stores all the necessary reagents for sample preparation, reverse transcription, polymerase chain reaction (PCR), and detection to isolate, amplify, and detect nucleic acid from multiple respiratory pathogens within a single NPS or TS specimen. After sample collection, the user injects hydration solution and sample combined with Sample Buffer into the pouch, places the pouch into the SPOTFIRE System, and starts a run. The entire run process takes about 15 minutes. Additional detail can be found in the SPOTFIRE System Operator's Manual.

During a run, the SPOTFIRE System:

- Lyses the sample by agitation (bead beating) in addition to chemical lysis mediated by the Sample Buffer.
- Extracts and purifies all nucleic acids from the sample using magnetic bead technology.
- Performs nested multiplex PCR by:
 - First performing reverse transcription, followed by a multiplexed first stage PCR reaction (PCR1).
 - Then performing multiple simultaneous second-stage PCR reactions (PCR2) in the array to amplify sequences within the PCR1 products.
- Uses endpoint melting curve data to detect target-specific amplicons and analyzes the data to generate a result for each analyte.



MATERIALS PROVIDED

The SPOTFIRE R/ST Panel contains materials consisting of primers, buffers, dNTPs, polymerase, molecular grade water, guanidinium chloride (50 - < 60%), Triton-X 100 (10 - < 20%), and LCGreen® Plus.

Each kit contains sufficient reagents to test 30 samples

- Individually packaged SPOTFIRE R/ST Panel pouches

- Sample Preparation Reagent Kits (SPRKs)
 - Single-use Sample Buffer ampoule
 - Single-use pre-filled Hydration Injection Vial (blue)
 - Single-use Sample Injection Vial (red)
 - Individually packaged Fixed-Volume Transfer Pipette
- BIOFIRE® SPOTFIRE® Respiratory/Sore Throat (R/ST) Panel Software
This software is required to run the SPOTFIRE R/ST Panel on the SPOTFIRE System.

MATERIALS REQUIRED BUT NOT PROVIDED

- SPOTFIRE System
 - BIOFIRE® SPOTFIRE® Control Station
 - BIOFIRE® SPOTFIRE® Module(s)
 - BIOFIRE® Pouch Loading Station
- 10% bleach solution or a similar disinfectant

The following materials are examples of those that are compatible with the requirements listed in the Sample Requirements section:

Sample Type	Collection Material	Part Number
Nasopharyngeal Swab (NPS)	Flexible Minitip Size Nylon® flocked swab with 100mm breakpoint	BD 220252 OR Copan 553C/503CS01
	BD™ 3 mL Universal Viral Transport Media	BD 220220 OR Copan 3C047N
	BD™ 3 mL Universal Viral Transport Media with Flexible Minitip Size Nylon® flocked swab with 100mm breakpoint	BD 220531 OR Copan 3C057N
Throat Swab (TS)	Flexible Regular Size Nylon® flocked swab with 80mm breakpoint with 1 mL Amies Medium	Copan Eswab™ 480C

- Additional acceptable media types for NPS collection:
 - Remel MicroTest™ M4RT® Multi-Microbe Media (ThermoFisher part number R12700)
 - Remel MicroTest™ M4RT® Multi-Microbe Media with MicroTip Flocked Swab (ThermoFisher part number R12566)



Note: Compatibility of the SPOTFIRE R/ST Panel with the above-named alternative transport media has been demonstrated analytically, however, clinical performance has not been established.

WARNINGS AND PRECAUTIONS

General Precautions

1. A trained healthcare professional should carefully interpret the results from the SPOTFIRE R/ST Panel in conjunction with a patient's signs and symptoms, results from other diagnostic tests, and relevant epidemiological information.
2. The BIOFIRE System Software displays step-by-step on-screen instructions for the test procedure, including selection of the applicable sample type for the test. Failure to select the correct sample type leads to incomplete reporting of results for all the analytes applicable to the sample type and to reporting of results for some analytes that are not appropriate to the sample type. If the incorrect sample type is selected, do not report the results; repeat the test from the same sample using a new reagent pouch and by selecting the correct sample type from the SPOTFIRE R/ST Panel menu.
3. SPOTFIRE R/ST Panel pouches are only for use with the SPOTFIRE System.
4. Always check the expiration date on the kit. Do not use a kit components or a pouch after its expiration date.
5. A desiccant packet is included in each pouch canister to preserve the stability of the SPOTFIRE R/ST Panel. Do not use a pouch if a desiccant packet is not present in the can.
6. Performance characteristics of the SPOTFIRE R/ST Panel have only been determined with nasopharyngeal swab (NPS) and throat swab (TS) specimens in transport medium.
7. SPOTFIRE R/ST Panel pouches are stored under vacuum in individually wrapped canisters. To preserve the integrity of the pouch vacuum for proper operation, be sure that a Module will be available and operational before unwrapping any pouches for loading.
8. If infection with SARS-CoV-2 is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions. Refer to the CDC Interim Guidelines for Collecting and Handling of Clinical Specimens for COVID-19 Testing for more information. <https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>.

Safety Precautions

1. Wear appropriate Personal Protective Equipment (PPE), including (but not limited to) disposable clean powder-free gloves and lab coats (if available). Protect skin, eyes, and mucous membranes. Change gloves often when handling reagents or samples.
2. Handle all samples and waste materials as if they were capable of transmitting infectious agents. Observe safety guidelines such as those outlined in:
 - CDC/NIH *Biosafety in Microbiological and Biomedical Laboratories*⁶⁴
 - CLSI Document M29 *Protection of Laboratory Workers from Occupationally Acquired Infections*⁶⁵
 - Refer to Interim Laboratory Biosafety Guidelines for Handling and Processing Specimens Associated with Coronavirus Disease 2019 (COVID-19) www.cdc.gov/coronavirus/2019-nCoV/lab-biosafety-guidelines.html or more current guidelines specific for SARS-CoV-2.
3. Follow your institution's safety procedures for handling biological samples.
4. If infection with SARS-CoV-2 is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control

precautions. Viral culture should not be attempted in cases of positive results for SARS-CoV-2 and/or any similar microbial agents unless a facility with an appropriate level of laboratory biosafety (e.g., BSL 3 and BSL 3+, etc.) is available to receive and culture specimens.

5. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to a state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
6. Dispose of materials used in this assay, including reagents, samples, and used buffer vials, according to federal, state, and local regulations.
7. Sample Buffer contains Guanidinium chloride and Triton X100. The following statements apply:
 - Health Hazards
 - Acute Toxicity, oral (Category 4)
 - H302 – Harmful if swallowed.
 - Skin corrosion/irritation (Category 2)
 - H315 - Causes skin irritation.
 - Serious eye damage/eye irritation (Category 1)
 - H318 - Causes serious eye damage.
 - Environment Hazards
 - Hazardous to the aquatic environment, acute aquatic hazard (Category 1)
 - H400 - Very toxic to aquatic life.
 - Hazardous to the aquatic environment, long-term aquatic hazard (Category 1)
 - H410 - Very toxic to aquatic life with long lasting effects.
 - Precautionary Statements
 - P260 – Do not breathe vapor.
 - P273 – Avoid release to the environment.
 - P280 – Wear protective gloves/protective clothing/eye protections/face protection.
 - Response
 - P391 - Collect spillage.
 - P332 + P313 - If skin irritation occurs: Get medical advice/attention.
 - P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 - P301 + P312 - IF SWALLOWED: Call a POISON CENTRE/doctor if you feel unwell.
 - P337 + P313 - If eye irritation persists: Get medical advice/attention.
 - Disposal
 - Dispose of waste and residues in accordance with local authority requirements.

Please refer to the SPOTFIRE R/ST Panel Safety Data Sheet (SDS) for more information:
<https://www.biofiredx.com/e-labeling/ITI0225>

8. Sample Buffer will form hazardous compounds and fumes when mixed with bleach or other disinfectants.



CAUTION: Never add Bleach to Sample Buffer or sample waste.

9. Bleach, a recommended disinfectant, is corrosive and may cause severe irritation or damage to eyes and skin. Vapor or mist may irritate the respiratory tract. Bleach is harmful if swallowed or inhaled.
- Eye contact: Hold eye open and rinse with water for 15-20 minutes. Remove contact lenses after the first 5 minutes and continue rinsing eye. Seek medical attention.
 - Skin contact: Immediately flush skin with plenty of water for at least 15 minutes. If irritation develops, seek medical attention.
 - Ingestion: Do not induce vomiting. Drink a glassful of water. If irritation develops, seek medical attention.
 - Please refer to the appropriate Safety Data Sheet (SDS) for more information.

Laboratory Precautions

1. Preventing organism contamination


Due to the sensitive nature of the SPOTFIRE R/ST Panel, it is important to guard against contamination of the sample and work area by carefully following the testing process outlined in this instruction document, including these guidelines:

- Personnel collecting and/or testing specimens may carry or shed common respiratory pathogens asymptomatically and can inadvertently contaminate the specimen while it is being processed. Careful adherence to the sample processing steps described in this document is recommended to avoid possible contamination. Samples may be processed in a clean biosafety cabinet (if available) or according to local/laboratory guidelines. If a biosafety cabinet is not used, a dead air box (e.g., AirClean PCR workstation), a splash shield (e.g., Bel-Art Scienceware Splash Shields), or a face shield may be used when preparing samples.
- Personnel with active respiratory symptoms (runny nose, cough) should wear a standard surgical mask (or equivalent) and should avoid touching the mask while handling specimens.
- It is recommended to avoid handling specimens or pouches in an area used to routinely process respiratory pathogen culture, and/or immunofluorescence testing, unless the area is thoroughly cleaned first.
- Prior to processing specimens, thoroughly clean both the work area and the Pouch Loading Station using a suitable cleaner such as freshly prepared 10% bleach or a similar disinfectant. To avoid residue build-up and potential damage to the specimen or interference from disinfectants, wipe disinfected surfaces with water.
- Specimens and pouches should be handled and/or tested one-at-a-time. Always change gloves and clean the work area between each pouch and specimen.
- Use clean gloves when removing Sample Buffer ampoules and Sample/Hydration Injection Vials from the SPRK package.
- Avoid collecting or handling specimens in areas that are exposed to vaccine material for pathogens detected by the SPOTFIRE R/ST Panel (e.g., influenza, SARS-CoV-2, *Bordetella pertussis*, and poliovirus (Human rhinovirus/enterovirus)). Vaccines may contain PCR-detectable DNA or RNA. If possible, particular care should be taken to avoid contamination of the specimen or testing areas (especially with nasal spray vaccines such as FluMist® and *B. pertussis* acellular vaccines such as Pentacel®, Daptacel®, and Adacel®; <http://www.cdc.gov/pertussis/clinical/diagnostic-testing/diagnosis-pcr-bestpractices.html>). Contamination of specimens or testing materials with vaccine can cause false-positive results.

2. Preventing amplicon contamination

A common concern with PCR-based assays is false positive results caused by contamination of the work area with PCR amplicon. Because the SPOTFIRE R/ST Panel pouch is a closed system, the risk of amplicon contamination is low provided that pouches remain intact after the test is completed. Adhere to the following guidelines, in addition to those above, to prevent amplicon contamination:

- Discard used pouches in a biohazard container immediately after the run has completed.
- Avoid excessive handling of pouches after test runs.
- Change gloves after handling a used pouch.
- Avoid exposing pouches to sharp edges or anything that might cause a puncture.

 **CAUTION:** If liquid is observed on the exterior of a pouch, the liquid and pouch should be immediately contained and discarded in a biohazard container. The instrument and workspace must be decontaminated as described in the SPOTFIRE System Operator's Manual.

DO NOT PERFORM ADDITIONAL TESTING UNTIL THE AREA HAS BEEN DECONTAMINATED

3. Transport media may contain non-viable organisms and/or nucleic acids at levels that can be detected by the SPOTFIRE R/ST Panel.

The presence of non-viable organisms and/or nucleic acids in transport media may lead to false positive test results.

Precautions Related to Public Health Reporting

Local, state, and federal regulations for notification of reportable disease are continually updated and include a number of organisms for surveillance and outbreak investigations. Additionally, the Centers for Disease Control and Prevention (CDC) recommends that when pathogens from reportable diseases are detected by a culture independent diagnostic test (CIDT), the laboratory should facilitate obtaining the isolate or clinical materials for submission to the appropriate public health laboratory to aid in outbreak detection and epidemiological investigations. Laboratories are responsible for following their state and/or local regulations and should consult their local and/or state public health laboratories for isolate and/or clinical sample submission guidelines.

Pertussis is a nationally notifiable infectious condition in some regions. If *Bordetella pertussis* is detected, notify the local health departments, if applicable.

Positive results for Coronavirus SARS-CoV-2 or suspected novel influenza should be reported to state, local, or federal health departments according to local reporting requirements. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

Precaution Related to REACH Regulation (EC 1907/2006)

This statement only applies to countries within the European Union (EU) with regard to the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) Regulation (EC 1907/2006):

It is recommended that all material associated with the test, including the material used to clean up spills, contaminated packaging, and/or unused and expired IVD tests, is incinerated. Please ensure that you follow local regulations regarding disposal.

REAGENT STORAGE, HANDLING, AND STABILITY

1. Store the test kit, including reagent pouches and buffers, at room temperature (15–25 °C).
2. Avoid storage of any materials near heating or cooling vents or in direct sunlight.
3. Always check the expiration date on the kit. Do not use reagents beyond the expiration date printed on the pouch or kit.
4. A desiccant packet is included in each pouch canister to preserve the stability of the SPOTFIRE R/ST Panel. Do not use a pouch if a desiccant packet is not present in the can.
5. All kit components should be stored and used together. Do not use components from one kit with those of another kit. Discard any extra components from the kit after all pouches have been consumed.
6. Do not remove pouches from their packaging until a sample is ready to be tested. Once the pouch packaging has been opened, the pouch should be loaded as soon as possible (within approximately 30 minutes).
7. Once a pouch has been loaded, the test run should be started as soon as possible (within approximately 60 minutes). Do not expose a loaded pouch to temperatures above 40°C (104°F) prior to testing.

SAMPLE REQUIREMENTS

The following table describes the requirements for specimen collection, preparation, and handling that will help ensure accurate test results.

Specimen Type	<p>Nasopharyngeal swab (NPS) or throat swab (TS) collected according to standard technique and immediately placed in 1 to 3 mL of transport media.</p> <p>Detailed NPS and TS specimen collection instructions can be found in the SPOTFIRE R/ST Panel Quick Guide</p>
Minimum Sample Volume	0.3 mL (300 µL)
Transport and Storage	<p>Specimens should be tested with the SPOTFIRE R/ST Panel as soon as possible.</p> <p>If storage is required, specimens can be held:</p> <ul style="list-style-type: none"> • At room temperature for up to 4 hours (15-25 °C) • Refrigerated for up to 3 days (2-8 °C) • Frozen (≤ -15 °C) (for up to 30 days)^a

^a Frozen storage for up to 30 days was evaluated for NPS and TS sample types. However, longer frozen storage at -70°C or lower may be acceptable. Please follow your institution's rules and protocols regarding sample storage validation.



Note: Specimens should not be centrifuged before testing.



Note: Bleach can damage organisms/nucleic acids within the specimen, potentially causing false negative results. Contact between bleach and specimens during collection, disinfection, and testing procedures should be avoided.

QUALITY CONTROL

Internal Process Controls

Two process controls are included in each pouch:

1. RNA Process Control

The RNA Process Control assay targets an RNA transcript from the yeast *Schizosaccharomyces pombe*. The yeast is present in the pouch in a freeze-dried form and becomes rehydrated when sample is loaded. The control material is carried through all stages of the test process, including lysis, nucleic acid purification, reverse transcription, PCR1, dilution, PCR2, and DNA melting. A positive control result indicates that all steps carried out in the SPOTFIRE R/ST Panel pouch were successful.

2. PCR2 Control

The PCR2 Control assay detects a DNA target that is dried into wells of the array along with the corresponding primers. A positive result indicates that PCR2 was successful.

Both control assays must be positive for the test run to pass. If the controls fail, the sample should be retested using a new pouch.

External Controls

External controls should be used in accordance with laboratory protocols and the appropriate accrediting organization requirements, as applicable. Transport media can be used as an external negative control. Previously characterized positive samples or negative samples spiked with well-characterized organisms can be used as external positive controls. Commercial external control materials may be available from other manufacturers; these should be used in accordance with the manufacturers' instructions and appropriate accrediting organization requirements, as applicable.

PROCEDURE

Refer to the SPOTFIRE R/ST Panel Quick Guide or the SPOTFIRE System Operator's Manual for more detail and pictorial representations of these instructions.

Use clean gloves and other Personal Protective Equipment (PPE) when handling pouches and samples. Only prepare one SPOTFIRE R/ST Panel pouch at a time and change gloves between samples and pouches. Once sample is added to the pouch, promptly transfer to the SPOTFIRE System to start the run. After the run is complete, discard the pouch in a biohazard container.


Step 1: Prepare Pouch

1. Thoroughly clean the work area and the Pouch Loading Station with freshly prepared 10% bleach (or suitable disinfectant) followed by a water rinse.
2. Wearing clean gloves, remove the pouch from its vacuum-sealed package by tearing or cutting the notched outer packaging and opening the protective canister.

Note: The pouch may still be used even if the vacuum seal of the pouch is not intact. Attempt to hydrate the pouch using the steps in the Hydrate Pouch section. If hydration is successful, continue with the run. If hydration fails, discard the pouch and use a new pouch to test the sample.

3. Check the expiration date on the kit. Do not use expired kit components or pouches.
4. Label the SPOTFIRE R/ST Panel pouch with the Sample ID and confirm the pouch and Sample IDs match.


Note: Do not cover the pouch barcode located on the left side of the label with the Sample ID.



5. Insert the pouch into the Pouch Loading Station, aligning the red and blue labels on the pouch with the red and blue arrows on the Pouch Loading Station.
6. Remove injection vials and Sample Buffer ampoule from the SPRK.
7. Place the injection vials into the Pouch Loading Station as follows:
 - Place a **red-capped Sample Injection Vial** into the **red well** of the Pouch Loading Station.
 - Place a **blue-capped Hydration Injection Vial** into the **blue well** of the Pouch Loading Station.
8. Reserve Transfer Pipette and Sample Buffer ampoule for use in Step 3.

Step 2: Hydrate Pouch

1. Unscrew **Hydration Injection Vial**, leaving the blue plastic cover in well of the Pouch Loading Station.
2. Insert the tip of the **Hydration Injection Vial's** blunt needle into the **pouch hydration port** located directly below the blue arrow of the Pouch Loading Station.
3. Forcefully push down in a firm and quick motion to puncture seal until a faint "pop" is heard and there is an ease in resistance. The correct volume of liquid will be pulled into the pouch automatically.
4. Leave **Hydration Injection Vial** in pouch.



Note: Some hydration solution may remain in the **Hydration Injection Vial**

- If the hydration solution is not automatically drawn into the pouch, repeat Step 2 to verify that the seal of the **pouch hydration port** was broken. If hydration solution is again not drawn into the pouch, discard the current pouch, retrieve a new pouch, and repeat from *Step 1: Prepare Pouch*.
5. Verify that the pouch has been hydrated.
 - Flip the barcode label down and check to see that fluid has entered the reagent wells (located at the base of the rigid plastic part of the pouch). Small air bubbles may be seen in some wells.
 - If the pouch fails to hydrate (dry reagents appear as white pellets), repeat Step 2 to verify that the seal of the **pouch hydration port** was broken. If hydration solution is still not drawn into the pouch, discard the current pouch, retrieve a new pouch, and repeat from *Step 1: Prepare Pouch*.
 6. Discard the **Hydration Injection Vial** in an appropriate biohazard container.

Step 3: Prepare Patient or QC Sample


1. Remove the lid from the sample or QC material tube.


 **Note: Positive QC materials are a PCR contamination risk.**

2. Transfer sample to the **Sample Injection Vial** using the following steps:

- Remove the fixed-volume Transfer Pipette from its individual packaging within the SPRK, being careful not to touch the tip.
- Fully squeeze the bulb of the fixed-volume Transfer Pipette and then lower the tip into the sample tube below the surface of the liquid.
- Slowly release the bulb completely, drawing up liquid into the fixed-volume Transfer Pipette.
- Visually confirm that liquid was drawn approximately one-third up the length of the stem of the fixed-volume Transfer Pipette, as shown in the image to the right.



 **Note: In rare cases, the fixed-volume Transfer Pipette may not draw the correct amount of liquid. If the liquid volume is not drawn approximately one-third up the length of the stem of the fixed-volume Transfer Pipette, discard it and retrieve a new fixed-volume Transfer Pipette from another SPRK package.**




 **Note: When performing QC testing, some liquid may be leftover in the QC material tube.**

- Raise the fixed-volume Transfer Pipette out of the sample liquid.
- Dispense the sample into the **Sample Injection Vial** by again fully squeezing the bulb of the fixed-volume Transfer Pipette.
- After dispensing the sample, discard the fixed-volume Transfer Pipette in an appropriate biohazard container.



3. Add Sample Buffer to the **Sample Injection Vial**:

- Open the Sample Buffer ampoule by twisting the tab off the tip.
- Dispense full volume of Sample Buffer by squeezing it into the **Sample Injection Vial**.
- Discard Sample Buffer ampoule in an appropriate biohazard container.
- Tightly close the lid of the **Sample Injection Vial**.



 **Note: To avoid contamination, do not touch the open tip of the Sample Buffer ampoule.**


 **CAUTION—The Sample Buffer is harmful if swallowed and can cause serious eye damage and/or skin irritation. Use appropriate PPE.**

4. Lift entire **Sample Injection Vial** out of the Pouch Loading Station and mix sample by gently inverting 3 times.
5. Return **Sample Injection Vial** to red well of the Pouch Loading Station.



Step 4: Load Patient or QC Sample

1. Slowly twist to unscrew the **Sample Injection Vial** and pause for 5 seconds with the blunt needle remaining in cap to avoid dripping.

 **Note:** *Pausing 5 seconds decreases the risk of dripping and contamination from the sample.*


2. Remove the **Sample Injection Vial**, leaving **red** plastic cover in the well of the Pouch Loading Station.

 **Note:** *To avoid contamination, do not touch the tip of blunt needle*

3. Insert the tip of the blunt needle into the pouch sample port located directly below the red arrow of the Pouch Loading Station.
4. Forcefully push down in a firm and quick motion to puncture seal (a faint “pop” is heard) and sample is pulled into the pouch by vacuum.
5. Wait 5 seconds as the sample mix is pulled into the pouch.

 **Note:** *Some sample mix may remain in the **Sample Injection Vial***

6. Verify that the sample has been loaded.
 - Flip the barcode label down and check to see that fluid has entered the reagent well next to the sample loading port.
 - If the pouch fails to pull sample from the **Sample Injection Vial**, the pouch should be discarded. Retrieve a new pouch and repeat from Step 1: Prepare Pouch.
7. Remove **Hydration Injection Vial** and **Sample Injection Vial** from SPOTFIRE R/ST Panel pouch. Screw each vial back into its plastic cover in the Pouch Loading Station before disposing of the vials according to institutional guidelines.
8. Remove the pouch from the Pouch Loading Station.



Step 5: Start Run

The BIOFIRE® SPOTFIRE® Software includes step-by-step, on-screen instructions that guide the operator through performing a run. Brief instructions for the SPOTFIRE System are given below. Refer to the SPOTFIRE System Operator's Manual for more detailed instructions.

1. Ensure that the SPOTFIRE System is powered on and the software is launched.
2. Select an available Module on the Home Screen and follow on-screen instructions to run test.

 **Note:** *If running QC samples, select the QC icon at the top of the screen and follow on-screen instructions.*

 **Caution:** *QC Testing should only be performed after selecting the QC icon. Failure to do so may lead to erroneous QC results.*

3. Scan the barcode on the pouch using the barcode scanner.



- Pouch identification (Lot Number and Serial Number), Pouch Type, and Protocol information will be automatically entered when the barcode is scanned. If it is not possible to scan the barcode, the pouch Lot Number, Serial Number, Expiration Date, and Pouch Type can be manually entered from the information provided on the pouch label into the appropriate fields.
4. Scan Sample ID.
 - The Sample ID can be entered manually or scanned in by using the barcode scanner when a barcoded Sample ID is used.
 5. Select either “Nasopharyngeal Swab” or “Throat Swab” as the sample type.

 **CAUTION:** Failure to select the appropriate sample type will lead to incorrect test results. Refer to the General Precautions section for more information.

 **Note:** Selecting the appropriate sample type is not required when testing QC samples.

6. Review the entered run information on the screen before inserting the pouch.

 **Note:** The selected Module's front panel LED will blink blue, indicating it is ready to accept a pouch.

7. Insert the pouch into the Module that is blinking blue. The Module will grab onto the pouch and pull it into the chamber and automatically start the run.
 - Once the run has started, the screen displays the panel name and Sample ID, and the minutes remaining on the run.

 **Note:** The selected Module's front panel LED will turn solid green to indicate that the run is in progress.

- When the run is finished, the pouch will automatically eject from the SPOTFIRE System.
8. Use gloves to remove pouch from System and discard the pouch in an appropriate biohazard container.
 9. Results are automatically created upon completion of a run. The test report can be viewed by clicking the Complete tile on the Home Screen or by clicking the appropriate icon for Patient Test Results or QC. The run file is automatically saved in the system database, and the test report can be viewed, printed, and/or saved as a PDF file.

INTERPRETATION OF RESULTS

Assay Interpretation

When PCR2 is complete, the instrument performs a high-resolution DNA melting analysis on the PCR products and records the change in fluorescence signal generated in each well (for more information refer to the SPOTFIRE System Operator's Manual). The software then performs several analyses and assigns a final assay result. The steps in the analyses are described below.

Analysis of melt curves. The software evaluates the DNA melt curve for each well of the PCR2 array to determine if a PCR product was present in that well. If the melt profile indicates the presence of a PCR product, then the analysis software calculates the melting temperature (T_m) of the curve and compares it against the expected T_m range for the assay. If the software determines that the T_m falls inside the assay-specific T_m range, the melt curve is called positive. If the software determines that the melt curve is not in the appropriate T_m range, the melt curve is called negative.

Analysis of replicates. Once melt curves have been identified, the software evaluates the three replicates for each assay to determine the assay result. For an assay to be called positive, at least two of the three associated melt curves must be

called positive, and the T_m for at least two of the three positive melt curves must be similar. Assays that do not meet these criteria are called negative.

Organism Interpretation

For most organisms detected by the SPOTFIRE R/ST Panel, the organism is reported as Positive if a single corresponding assay is positive. For example, Human metapneumovirus will have a test report result of Human metapneumovirus Positive if the hMPV assay is positive (at least two of the three hMPV assay wells on the array have similar positive melt peaks with T_m values that are within the assay-specific T_m range). The test results for Adenovirus, Coronavirus (seasonal), Coronavirus SARS-CoV-2, Influenza A virus, and Parainfluenza depend on the interpretation of results from more than one assay. Interpretation and actions for the multi-assay results are provided below.

Adenovirus

Four assays are included for the detection of Adenovirus. If one assay or any combination of assays is positive, the test report result will be Adenovirus Positive. If all assays are negative, the test report result will be Adenovirus Negative.

Coronavirus (seasonal)

The SPOTFIRE R/ST Panel uses four assays for the detection of seasonal Coronaviruses. The software interprets each of these assays and the results are combined as a final test result. If one assay or any combination of assays is positive, the test report result will be Coronavirus (seasonal) Positive. If all assays are negative, the test report will be Coronavirus (seasonal) Negative.

Coronavirus SARS-CoV-2

Two assays are included for the detection of SARS-CoV-2. One assay targets the spike protein (S) gene and the other targets the membrane protein (M) gene. If either one or both assays are positive, the test report result will be Coronavirus SARS-CoV-2 Positive. If both assays are negative, the test report result will be Coronavirus SARS-CoV-2 Negative.

Influenza A

The Influenza A virus result is determined by four assays: two pan-influenza A assays and two subtype assays that differentiate between hemagglutinin (HA) types (H1-2009 and H3). Results from all four assays are evaluated to determine the reported Influenza A virus result (See Table 2). A positive Influenza A virus result requires at least two of the four assays to be positive, including at least one pan-FluA assay. The subtype (H1-2009 or H3) is reported as Positive if the Influenza A virus result is Positive and the corresponding HA assay is positive. If all four assays are negative, the test report result will be Influenza A virus Negative. If only one assay is positive or if only the two HA assays are positive, the test report result will be Influenza A virus Uncertain. An Influenza A virus Uncertain result could occur when the amount of the virus in the specimen is low and not detected by one or more required assays (Table 2). An Influenza A virus Uncertain result could also indicate the presence of an atypical influenza A subtype (e.g. avian H7N9 or H5N1 types), or a novel influenza A strain. Specimens with an Uncertain result should be retested once.

Table 2. Possible Influenza A Virus Assay Results, Interpretations, and Actions Required

Result Summary	Interpretation Results		Assay Result Combinations			Action Required
	Interpretation	Result	Pan FluA- Assays (N=2)	H1-2009 Subtype Assay	H3 Subtype Assay	
POSITIVE: Influenza A virus (subtype H1-2009)	Influenza A virus	Positive	≥1 Positive	Positive	Negative	None
	Influenza A virus A/H1-2009	Positive				
	Influenza A virus A/ H3	Negative				
POSITIVE: Influenza A virus (subtype H3)	Influenza A virus	Positive	≥1 Positive	Negative	Positive	
	Influenza A virus A/H1-2009	Negative				
	Influenza A virus A/H3	Positive				
POSITIVE: Influenza A virus (multiple subtypes)	Influenza A virus	Positive	≥1 Positive	Positive	Positive	
	Influenza A virus A/H1-2009	Positive				
	Influenza A virus A/H3	Positive				
POSITIVE: Influenza A virus (no subtype identified)	Influenza A virus	Positive	2 Positive	Negative	Negative	Retest
	Influenza A virus A/H1-2009	Negative				
	Influenza A virus A/H3	Negative				
UNCERTAIN: Influenza A virus ^a	Influenza A virus	Uncertain	1 Positive Negative Negative Negative	Negative Positive Negative Positive	Negative Negative Positive Positive	
	Influenza A virus A/H1-2009	Negative				
	Influenza A virus A/H3	Negative				
NEGATIVE	Influenza A virus	Negative	Negative	Negative	Negative	None
	Influenza A virus A/H1-2009	Negative				
	Influenza A virus A/H3	Negative				

^aAny one of the four listed combinations of assay results will generate an Influenza A virus Uncertain interpretation with both the Influenza A virus A/H1-2009 and Influenza A virus A/H3 interpretations reported as Negative.

Influenza A virus (no subtype identified)

If both FluA-pan assays are positive, but none of the hemagglutinin subtyping assays are positive, then the result is Influenza A virus (no subtype identified). This result could occur when the amount of the virus in the specimen is low and not detected by the subtyping assays. This result could also indicate the presence of a novel influenza A strain. In both cases, the sample in question should be retested. If the retest provides the same Influenza A virus (no subtype identified) result, contact the appropriate public health authorities for confirmatory testing.

Parainfluenza virus


The SPOTFIRE R/ST Panel uses four assays for the detection of Parainfluenza virus. The software interprets each of these assays and the results are combined as a final test result. If one assay or any combination of assays is positive, the test report result will be Parainfluenza virus Positive. If all assays are negative, the test report will be Parainfluenza virus Negative.

SPOTFIRE R/ST Panel Test Report

The SPOTFIRE R/ST Panel test report is automatically displayed upon completion of a run and can be printed or saved as a PDF file. Each report contains a Run Summary, a Result Summary, and a Run Details section.

Run Summary

The Run Summary section of the test report provides the test type (Respiratory Menu or Sore Throat Menu), time and date of the run, Sample ID, and the identity of the operator that performed the test.

 **Caution:** Confirm that the test report displays the appropriate test menu (i.e., “Respiratory” or “Sore Throat”) and that this matches the ordered test prior to reporting results.

Results Summary

The Result Summary section of the test report lists the overall results of the test. Possible test results include Negative, Positive, Uncertain (Influenza A virus only), and Invalid. An Action Bar will appear underneath the test results only when further action is necessary.

The result for each organism tested by the panel is also shown. A check mark next to an organism indicates a positive result, while a question mark next to an organism indicates an uncertain result (Influenza A virus only). When no symbol is present, the result was negative. If the run result is Invalid, this section is not displayed.

Table 3 provides an explanation for each interpretation of patient test results and any follow-up necessary to obtain a final result. Table 4 provides an explanation for each interpretation of a QC test result and any follow-up necessary to obtain a final result. Table 5 provides an explanation for Invalid results.

Table 3. Interpretation of Results

Result	Explanation	Action
NEGATIVE	Test controls Pass AND Test is negative for all organisms	Report the results
POSITIVE: [Organism name(s)]	Test controls Pass AND Test is positive for the organism(s) listed	Report the results
POSITIVE: Multiple Organisms (4+)*	Test controls Pass AND Test is positive for four or more organisms	Detection of four or more organisms is possible but rare. If contamination is suspected, clean the area and retest the sample, then report the results of the retest. If additional guidance is needed, contact your local bioMérieux subsidiary or distributor.
UNCERTAIN: Influenza A virus	Test controls Pass AND Results for Influenza A virus are inconclusive	Retest the original sample ONCE and report the results of the retest.
POSITIVE: Influenza A virus (no subtype identified)	Test controls Pass AND Results for Influenza A virus do not identify a specific subtype	Uncommon result, retest ONCE. If the retest provides the same results, contact the appropriate public health authorities for confirmatory testing.
POSITIVE: Influenza A virus (multiple subtypes) Influenza B virus	Test controls Pass AND Results indicate the presence of multiple Influenza A virus subtypes and Influenza B virus	Detection of multiple Influenza infections is possible but rare. This result could be caused by recent administration of a nasal vaccine (e.g. FluMist®) or environmental contamination with an Influenza vaccine. If contamination is suspected, clean the area (refer to Laboratory Precautions section) and retest the sample, then report the results of the retest. If additional guidance is needed, contact your local bioMérieux subsidiary or distributor.

* In rare cases, POSITIVE: Multiple Organisms (4+) and UNCERTAIN: Influenza A virus may both appear in the Result Summary section when exactly three organisms (each indicated by a checkmark) in addition to an Influenza A virus uncertain result (indicated by a question mark) are shown in the results section.

Table 4. Interpretation of QC Results

QC Result	Explanation	Action
PASS	Test controls Pass AND A Positive QC Test or Negative QC Test has Pass results	Report the results
FAIL	Test controls Pass AND A Positive QC Test has negative results	Retest the Positive QC material ONCE. If the failure persists, contact Customer Support for further instruction.
	Test controls Pass AND A Negative QC Test has positive results	If contamination is suspected, clean the area according to instructions on p.1 of the SPOTFIRE R/ST Panel Quick Guide and retest using a new Negative QC vial. If the failure persists, contact your local bioMérieux subsidiary or distributor for further instruction.

Table 5. Interpretation of Invalid Results

QC Result	Explanation	Action
INVALID: [Failure Reason]	Run does not complete due to any of the following failures: <ul style="list-style-type: none"> • Instrument error • Aborted run • Run incomplete • Software error • Operational conditions out of range • Internal control failure 	Follow on-screen instructions. If the failure persists, contact your local bioMérieux subsidiary or distributor.

Run Details Summary

The Run Details section provides additional information about the run including pouch information (pouch type, lot number, and serial number), control results, sample type, and the Module used to perform the test.

The Run Details will automatically be shown when printed and can also be accessed from the on-screen report by selecting the “Show Details” button.

The “Controls” field displays Pass only if the run completed successfully (no instrument or software errors) and both of the pouch control assays (RNA Process Control and PCR2 Control) were successful. The Control field will display Fail if the run was completed successfully (no instrument or software errors) but one or both pouch control assays failed. If the control result is Fail, then the results for all the tests on the panel are not displayed and the specimen will need to be retested with a new pouch. Table 6 provides additional information for each of the possible control field results.

Table 6. Interpretation of Controls Field on the SPOTFIRE R/ST Panel Test Report

Control Result	Explanation
Pass	The run was successfully completed AND Both pouch controls were successful.
Fail	The run was successfully completed BUT At least one of the pouch controls (RNA Process Control and/or PCR2 Control) failed.
Invalid	The controls are invalid because the run did not complete. (Typically, this indicates a software or hardware error).

LIMITATIONS

1. For prescription use only.
2. This product can be used only with the SPOTFIRE System.
3. The SPOTFIRE R/ST Panel is a qualitative test and does not provide a quantitative value for the virus(es) and/or bacteria detected in the specimen.
4. Results from this test must be correlated with the clinical history, epidemiological data, and other data available to the clinician evaluating the patient.
5. The SPOTFIRE R/ST Panel has been evaluated for use with human specimen material only.
6. The SPOTFIRE R/ST Panel has not been validated for testing of specimens other than nasopharyngeal swab (NPS) or throat swab (TS) specimens in transport medium.
7. The SPOTFIRE R/ST Panel has not been validated for the testing of pooled specimens or the screening of specimens from asymptomatic individuals that are not suspected of COVID-19 infection.
8. The performance of SPOTFIRE R/ST Panel has not been established for specimens collected from individuals without signs or symptoms of respiratory and/or pharyngitis infection.
9. The performance of the SPOTFIRE R/ST Panel has not been specifically evaluated for specimens from immunocompromised individuals.
10. The performance of the SPOTFIRE R/ST Panel has not been specifically evaluated in a population known to be vaccinated against illnesses caused by any of the SPOTFIRE R/ST Panel analytes (e.g. SARS-CoV-2 (COVID-19), influenza, pertussis (whooping cough), or RSV).
11. The effect of antibiotic or antiviral treatment on test performance has not been evaluated.
12. The performance of the SPOTFIRE R/ST Panel has not been specifically evaluated for specimens in patients with Multisystem Inflammatory Syndrome in Children (MIS-C) or similar syndromes.
13. The performance of the SPOTFIRE R/ST Panel has not been established with potentially interfering medications for the treatment of influenza or cold viruses. The effect of interfering substances has only been evaluated for those listed in the Interference section. Interference from substances that were not evaluated could lead to erroneous results.
14. The performance of the SPOTFIRE R/ST Panel has not been established for monitoring treatment of infection with any of the panel organisms.
15. The performance of SPOTFIRE R/ST Panel has not been established for screening of blood or blood products.
16. False positive and false negative results can be the result of a variety of sources and causes. It is important that these results be used in conjunction with other clinical, epidemiological, or laboratory information.
17. The detection of viral and bacterial nucleic acid is dependent upon proper specimen collection, handling, transportation, storage, and preparation. Failure to observe proper procedures in any one of these steps can lead to incorrect results. There is a risk of false positive or false negative values resulting from improperly collected, transported, or handled specimens.
18. A negative SPOTFIRE R/ST Panel result does not exclude the possibility of viral or bacterial infection. Negative test results may occur due to the presence of sequence variants (or mutation) in the region targeted by the assay, the presence of inhibitors, technical error, sample mix-up, an infection caused by an organism not detected by the panel, or lower respiratory tract infection that is not detected by a nasopharyngeal or throat swab specimen. Test results may also be affected by concurrent antiviral/antibacterial therapy or levels of organism in the specimen that

are below the limit of detection for the test. Negative results should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

19. If four or more organisms are reported as positive in a specimen, retesting is recommended to confirm the polymicrobial result.
20. Viral and bacterial nucleic acids may persist in vivo independent of organism viability. Detection of organism target(s) does not imply that the corresponding organisms are infectious or are the causative agents for clinical symptoms.
21. Positive and negative predictive values are highly dependent on prevalence. The likelihood of a negative result being false is higher during peak activity when prevalence of disease is high. The likelihood of a positive result being false is higher during periods when prevalence is moderate to low.
22. Performance characteristics for influenza A were established when influenza A/H1-2009 and A/H3 were predominant. When other influenza A viruses are emerging, performance characteristics may differ.
23. Due to the small number of positive results observed for certain organisms during the prospective clinical study, performance characteristics for *Bordetella parapertussis*, *Bordetella pertussis*, *Chlamydia pneumoniae*, Coronavirus 229E, Human metapneumovirus, Influenza A H3, Influenza A H1-2009, Influenza B, *Mycoplasma pneumoniae*, and parainfluenza virus (serotypes 1, 2, and 4) established primarily with retrospective clinical specimens.
24. The SPOTFIRE R/ST Panel influenza A subtyping assays target the influenza A hemagglutinin (H) gene only. The SPOTFIRE R/ST Panel does not detect or differentiate the influenza A neuraminidase (N) subtypes.
25. The SPOTFIRE R/ST Panel may not be able to distinguish between existing viral strains and new variants as they emerge. For example, the SPOTFIRE R/ST Panel can detect influenza A H3N2v (first recognized in August 2011) but will not be able to distinguish this variant from influenza A H3N2 seasonal. If variant virus infection is suspected, clinicians should contact their state or local health department to arrange specimen transport and request a timely diagnosis at a state public health laboratory.
26. Recent administration of nasal vaccines (e.g. FluMist®) prior to NPS or TS specimen collection could lead to accurate virus detection by the SPOTFIRE R/ST Panel of the viruses contained in the vaccine but would not represent infection by those agents.
27. Due to the genetic similarity between human rhinovirus and enterovirus, the SPOTFIRE R/ST Panel cannot reliably differentiate them. A SPOTFIRE R/ST Panel rhinovirus/enterovirus Positive result should be followed-up using an alternate method (e.g. cell culture or sequence analysis) if differentiation between the viruses is required.
28. The SPOTFIRE R/ST Panel detects a single-copy target (*fim2*, present at one copy per cell) in *B. pertussis*.
 - Other PCR tests for *B. pertussis* target the multi-copy IS481 insertion sequence and are therefore capable of detecting lower levels of *B. pertussis* (i.e. more sensitive).
 - The IS481 sequence is also present in *B. holmesii* (and to a lesser extent in *B. bronchiseptica*), whereas the SPOTFIRE R/ST Panel assay (Fim2) was designed to be specific for *B. pertussis*.
 - Due to lower sensitivity, the SPOTFIRE R/ST Panel *B. pertussis* assay is less susceptible than IS481 assays to the detection of very low levels of contaminating *B. pertussis* vaccine material. However, care must always be taken to avoid contamination of specimens with vaccine material as higher levels may still lead to false positive results with the SPOTFIRE R/ST Panel test (see contamination prevention guidelines).
29. Some strains of *B. bronchiseptica* (rarely isolated from humans) do carry IS1001 insertion sequences identical to those carried by most strains of *B. parapertussis*. These sequences will be amplified by the IS1001 assay and reported by SPOTFIRE R/ST Panel as *Bordetella parapertussis*.

30. The SPOTFIRE R/ST Panel Human rhinovirus/enterovirus assay may amplify off-target sequences found in strains of *B. pertussis*, *B. bronchiseptica* and *B. parapertussis*. Bovine and canine picornaviruses may also be detected and reported as Human rhinovirus/enterovirus when present at high concentration. Refer to Analytical Specificity (Cross-Reactivity & Exclusivity).
31. *Chlamydia gallinacea* may be detected and reported as *Chlamydia pneumoniae* when present at high concentration. Refer to Analytical Specificity (Cross-Reactivity & Exclusivity).
32. There is a risk of false positive results due to contamination with organisms, nucleic acids, vaccine material, amplified products, or from non-specific signals in the assay. Particular attention should be given to the Laboratory Precautions noted under the Warnings and Precautions section.
33. Transport media may contain non-viable organisms and/or nucleic acid at levels that can be detected by the SPOTFIRE R/ST Panel.
34. There is a risk of false positive results due to non-specific amplification and cross-reactivity with organisms found in the respiratory tract. Observed and predicted cross-reactivity for the SPOTFIRE R/ST Panel is described in the Analytical Specificity (Cross-Reactivity & Exclusivity) section. Erroneous results due to cross-reactivity with organisms that were not evaluated or new variant sequences that emerge is also possible.
35. Primers for both SPOTFIRE R/ST Panel SARS-CoV-2 assays share substantial sequence homology with the Bat coronavirus RaTG13 (accession: MN996532) and cross-reactivity with this closely related viral sequence is predicted. In addition, the SARS-CoV2-2 assay may cross-react with Pangolin coronavirus (accession: MT084071) and two other bat SARS-like coronavirus sequences (accession MG772933 and MG772934). It is unlikely that these viruses would be found in a human clinical specimen; but if present, the cross-reactive product(s) produced by the SPOTFIRE R/ST Panel will be reported as Coronavirus SARS-CoV-2.
36. The clinical performance has not been established for all circulating variants of SARS-CoV-2 but is anticipated to be reflective of the prevalent variants in circulation at the time and location of the clinical evaluation. Performance at the time of testing may vary depending on the variants circulating, including newly emerging strains of SARS-CoV-2 and their prevalence, which change over time. Refer to *In Silico* Reactivity Predictions for Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) Assays for more information.
37. The SPOTFIRE R/ST Panel pan-influenza A, A/H1-2009, A/H3 and influenza B assays are based on those included in other previously FDA-cleared BIOFIRE panels. Positive agreement with the SPOTFIRE R/ST Panel for detection of influenza A, A/H1-2009, A/H3 and influenza B was established using retrospective, archived specimens. The performance observed with prospectively collected specimens may differ.
38. Negative results for Influenza A virus A/H1-2009 or Influenza A virus A/H3 do not preclude infection with other Influenza A viruses.
39. The FluA-H1-2009 assay may react with the H1 hemagglutinin gene sequences from viruses of swine origin. The SPOTFIRE R/ST Panel will report either Influenza A virus (Subtype H1-2009) or Influenza A virus (No Subtype Identified), depending on the strain and concentration in the sample. In addition, certain strains of avian H1N2 may also be detected and reported as H1-2009.
40. Additional follow-up testing by culture is required if the SPOTFIRE R/ST Panel *Streptococcus pyogenes* assay result is negative and clinical symptoms persist, or in the event of an outbreak of acute rheumatic fever (ARF).

EXPECTED VALUES

In the prospective clinical evaluation of the SPOTFIRE R/ST Panel, 1120 nasopharyngeal swab (NPS) specimens and 877 throat swab (TS) specimens were collected from consented volunteers or obtained as residual, leftover specimens from subjects of all ages and tested at five study sites across the United States and one site in the United Kingdom over approximately 19 months (December 2020 to September 2021 for NPS and TS specimens, and September 2022 to May 2023 for TS specimens only). Expected value summaries (as determined by the SPOTFIRE R/ST Panel) are stratified by specimen enrollment site in Table 7 and Table 8 for NPS and TS specimens, respectively. Expected values are stratified by subject age group in Table 9 and Table 10 for NPS and TS specimens, respectively.

Table 7. Expected Value (EV; as determined by SPOTFIRE R/ST Panel positives / specimens analyzed) Summary by Site for NPS Specimens Collected During the SPOTFIRE R/ST Panel Prospective Clinical Evaluation (December 2020 to September 2021)

SPOTFIRE R/ST Panel R Menu Result	Overall		Site 1		Site 2		Site 3		Site 4		Site 5	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Viruses												
Adenovirus	56/1115	5.0%	3/361	0.8%	0/36	0%	34/508	6.7%	19/147	12.9%	0/63	0%
Coronavirus SARS-CoV-2	77/1110	6.9%	27/361	7.5%	7/36	19.4%	8/502	1.6%	8/147	5.4%	27/64	42.2%
Coronavirus (seasonal)	114/1115	10.2%	24/361	6.6%	4/36	11.1%	73/508	14.4%	13/147	8.8%	0/63	0%
Human metapneumovirus	1/1115	0.1%	1/361	0.3%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
Human rhinovirus/enterovirus	417/1115	37.4%	59/361	16.3%	7/36	19.4%	265/508	52.2%	85/147	57.8%	1/63	1.6%
Influenza A virus	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
Influenza A virus A/H1-2009	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
Influenza A virus A/H3	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
Influenza B virus	0/1110	0%	0/361	0%	0/36	0%	0/502	0%	0/147	0%	0/64	0%
Parainfluenza virus	107/1115	9.6%	20/361	5.5%	3/36	8.3%	72/508	14.2%	11/147	7.5%	1/63	1.6%
Respiratory syncytial virus	28/1115	2.5%	7/361	1.9%	0/36	0%	16/508	3.1%	5/147	3.4%	0/63	0%
Bacteria												
<i>Bordetella parapertussis</i>	0/1110	0%	0/361	0%	0/36	0%	0/502	0%	0/147	0%	0/64	0%
<i>Bordetella pertussis</i>	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
<i>Chlamydia pneumoniae</i>	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%
<i>Mycoplasma pneumoniae</i>	0/1115	0%	0/361	0%	0/36	0%	0/508	0%	0/147	0%	0/63	0%

^a #/SA = number of SPOTFIRE R/ST Panel positives divided by the number of specimens analyzed (i.e., with valid comparator results) for that analyte

^b EV = Expected Value

Table 8. Expected Value (EV; as determined by SPOTFIRE R/ST Panel positives / specimens analyzed) Summary by Site for TS Specimens Collected During the SPOTFIRE R/ST Panel Prospective Clinical Evaluation (December 2020 to September 2021 and September 2022 to May 2023)

SPOTFIRE R/ST Panel ST Menu Result	Overall		Site 1		Site 2		Site 3		Site 4		Site 5		Site 6	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Viruses														
Adenovirus	72/875	8.2%	0/152	0%	4/97	4.1%	30/324	9.3%	26/253	10.3%	1/22	4.5%	11/27	40.7%
Coronavirus SARS-CoV-2	46/876	5.3%	13/152	8.6%	5/97	5.2%	8/324	2.5%	18/253	7.1%	2/22	9.1%	0/28	0%
Coronavirus (seasonal)	48/875	5.5%	4/152	2.6%	3/97	3.1%	13/324	4.0%	23/253	9.1%	0/22	0%	5/27	18.5%
Human metapneumovirus	15/875	1.7%	3/152	2.0%	1/97	1.0%	5/324	1.5%	6/253	2.4%	0/22	0%	0/27	0%
Human rhinovirus/enterovirus	245/875	28.0%	17/152	11.2%	15/97	15.5%	125/324	38.6%	77/253	30.4%	1/22	4.5%	10/27	37.0%
Influenza A virus	35/875	4.0%	18/152	11.8%	6/97	6.2%	8/324	2.5%	3/253	1.2%	0/22	0%	0/27	0%
Influenza A virus A/H1-2009	12/875	1.4%	9/152	5.9%	1/97	1.0%	0/324	0%	2/253	0.8%	0/22	0%	0/27	0%
Influenza A virus A/H3	23/875	2.6%	9/152	5.9%	5/97	5.2%	8/324	2.5%	1/253	0.4%	0/22	0%	0/27	0%
Influenza B virus	4/876	0.5%	1/152	0.7%	0/97	0%	3/324	0.9%	0/253	0%	0/22	0%	0/28	0%

SPOTFIRE R/ST Panel ST Menu Result	Overall		Site 1		Site 2		Site 3		Site 4		Site 5		Site 6	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Parainfluenza virus	52/875	5.9%	2/152	1.3%	5/97	5.2%	33/324	10.2%	12/253	4.7%	0/22	0%	0/27	0%
Respiratory syncytial virus	23/875	2.6%	4/152	2.6%	0/97	0%	11/324	3.4%	8/253	3.2%	0/22	0%	0/27	0%
Bacteria														
<i>Chlamydia pneumoniae</i>	0/875	0%	0/152	0%	0/97	0%	0/324	0%	0/253	0%	0/22	0%	0/27	0%
<i>Mycoplasma pneumoniae</i>	0/875	0%	0/152	0%	0/97	0%	0/324	0%	0/253	0%	0/22	0%	0/27	0%
<i>Streptococcus dysgalactiae</i> (group C/G Strep)	32/869	3.7%	4/152	2.6%	6/97	6.2%	5/320	1.6%	16/251	6.4%	1/21	4.8%	0/28	0%
<i>Streptococcus pyogenes</i> (group A Strep)	212/869	24.4%	15/152	9.9%	27/97	27.8%	92/320	28.8%	68/251	27.1%	6/21	28.6%	4/28	14.3%

^a #/SA = number of SPOTFIRE R/ST Panel positives divided by the number of specimens analyzed (i.e., with valid comparator results) for that analyte

^b EV = Expected Value

Table 9 Expected Value (EV; as determined by SPOTFIRE R/ST Panel positives / specimens analyzed) Summary by Age Group for NPS Specimens Collected During the SPOTFIRE R/ST Panel Prospective Clinical Evaluation (December 2020 to September 2021)

SPOTFIRE R/ST Panel R Menu Result	Overall		≤5 years		6-18 years		19-40 years		41-60 years		61+ years	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Viruses												
Adenovirus	56/1115	5.0%	53/455	11.6%	3/257	1.2%	0/159	0%	0/146	0%	0/98	0%
Coronavirus SARS-CoV-2	77/1110	6.9%	9/450	2.0%	7/257	2.7%	27/158	17.1%	21/147	14.3%	13/98	13.3%
Coronavirus (seasonal)	114/1115	10.2%	71/455	15.6%	22/257	8.6%	10/159	6.3%	9/146	6.2%	2/98	2.0%
Human metapneumovirus	1/1115	0.1%	1/455	0.2%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
Human rhinovirus/enterovirus	417/1115	37.4%	236/455	51.9%	139/257	54.1%	23/159	14.5%	14/146	9.6%	5/98	5.1%
Influenza A virus	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
Influenza A virus A/H1-2009	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
Influenza A virus A/H3	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
Influenza B virus	0/1110	0%	0/450	0%	0/257	0%	0/158	0%	0/147	0%	0/98	0%
Parainfluenza virus	107/1115	9.6%	89/455	19.6%	4/257	1.6%	4/159	2.5%	8/146	5.5%	2/98	2.0%
Respiratory syncytial virus	28/1115	2.5%	21/455	4.6%	6/257	2.3%	0/159	0%	1/146	0.7%	0/98	0%
Bacteria												
<i>Bordetella parapertussis</i>	0/1110	0%	0/450	0%	0/257	0%	0/158	0%	0/147	0%	0/98	0%
<i>Bordetella pertussis</i>	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
<i>Chlamydia pneumoniae</i>	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%
<i>Mycoplasma pneumoniae</i>	0/1115	0%	0/455	0%	0/257	0%	0/159	0%	0/146	0%	0/98	0%

^a #/SA = number of SPOTFIRE R/ST Panel positives divided by the number of specimens analyzed (i.e., with valid comparator results) for that analyte

^b EV = Expected Value

Table 10 Expected Value (EV; as determined by SPOTFIRE R/ST Panel positives / specimens analyzed) Summary by Age Group for TS Specimens Collected During the SPOTFIRE R/ST Panel Prospective Clinical Evaluation (December 2020 to September 2021 and September 2022 to May 2023)

SPOTFIRE R/ST Panel ST Menu Result	Overall		≤5 years		6-18 years		19-40 years		41-60 years		61+ years	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Viruses												
Adenovirus	72/875	8.2%	25/137	18.2%	41/455	9.0%	5/186	2.7%	1/72	1.4%	0/25	0%
Coronavirus SARS-CoV-2	46/876	5.3%	2/138	1.4%	22/455	4.8%	13/186	7.0%	4/72	5.6%	5/25	20.0%
Coronavirus (seasonal)	48/875	5.5%	3/137	2.2%	35/455	7.7%	7/186	3.8%	3/72	4.2%	0/25	0%
Human metapneumovirus	15/875	1.7%	4/137	2.9%	7/455	1.5%	2/186	1.1%	2/72	2.8%	0/25	0%
Human rhinovirus/enterovirus	245/875	28.0%	58/137	42.3%	152/455	33.4%	23/186	12.4%	9/72	12.5%	3/25	12.0%
Influenza A virus	35/875	4.0%	2/137	1.5%	10/455	2.2%	15/186	8.1%	5/72	6.9%	3/25	12.0%
Influenza A virus A/H1-2009	12/875	1.4%	0/137	0%	2/455	0.4%	5/186	2.7%	3/72	4.2%	2/25	8.0%

SPOTFIRE R/ST Panel ST Menu Result	Overall		≤5 years		6-18 years		19-40 years		41-60 years		61+ years	
	#/SA ^a	EV ^b	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV	#/SA	EV
Influenza A virus A/H3	23/875	2.6%	2/137	1.5%	8/455	1.8%	10/186	5.4%	2/72	2.8%	1/25	4.0%
Influenza B virus	4/876	0.5%	1/138	0.7%	2/455	0.4%	1/186	0.5%	0/72	0%	0/25	0%
Parainfluenza virus	52/875	5.9%	25/137	18.2%	18/455	4.0%	6/186	3.2%	3/72	4.2%	0/25	0%
Respiratory syncytial virus	23/875	2.6%	6/137	4.4%	13/455	2.9%	0/186	0%	3/72	4.2%	1/25	4.0%
Bacteria												
<i>Chlamydia pneumoniae</i>	0/875	0%	0/137	0%	0/455	0%	0/186	0%	0/72	0%	0/25	0%
<i>Mycoplasma pneumoniae</i>	0/875	0%	0/137	0%	0/455	0%	0/186	0%	0/72	0%	0/25	0%
<i>Streptococcus dysgalactiae</i> (group C/G Strep)	32/869	3.7%	4/136	2.9%	15/451	3.3%	11/185	5.9%	2/72	2.8%	0/25	0%
<i>Streptococcus pyogenes</i> (group A Strep)	212/869	24.4%	35/136	25.7%	130/451	28.8%	39/185	21.1%	8/72	11.1%	0/25	0%

^a #/SA = number of SPOTFIRE R/ST Panel positives divided by the number of specimens analyzed (i.e., with valid comparator results) for that analyte

^b EV = Expected Value

Observed multiple detections (as determined by the SPOTFIRE R/ST Panel) during the prospective clinical evaluation are presented in Table 11. At least one analyte was detected in a total of 677 NPS specimens (60.4% positivity rate; 677/1120); at least one analyte was detected in a total of 604 TS specimens (68.9% positivity rate; 604/877). Polymicrobial detections of up to four organisms were observed in NPS specimens and polymicrobial detections of up to three organisms were observed in TS specimens.

Table 11 Expected Values (Multiple Detections as Determined by the SPOTFIRE R/ST Panel) for the SPOTFIRE R/ST Panel Prospective Clinical Evaluation (December 2020 to September 2021 for NPS and TS specimens, and September 2022 to May 2023 for TS specimens only)

SPOTFIRE R/ST Panel Result	Expected Value (as Determined by Testing of 1120 Prospective NPS Specimens)		Expected Value (as Determined by Testing of 877 Prospective TS Specimens)	
	Number Detected and Reported	% of Total (% of Positives)	Number Detected and Reported	% of Total (% of Positives)
Detected (at least one organism result)	677	60.4% (100%)	604	68.9% (100%)
One organism result	571	51.0% (84.3%)	441	50.3% (73.0%)
Two organism results	90	8.0% (13.3%)	146	16.6% (24.2%)
Three organism results	15	1.3% (2.2%)	17	1.9% (2.8%)
Four organism results	1	0.1% (0.1%)	0	0% (0%)

The SPOTFIRE R/ST Panel reported a total of 106 NPS specimens with discernable multiple organism detections (9.5% of all NPS specimens, 106/1120; 15.7% of positive NPS specimens, 106/677) and a total of 163 TS specimens with discernable multiple organism detections (18.6% of all TS specimens, 163/877; 27.0% of positive TS specimens, 163/604). The resulting co-detection combinations as reported by the SPOTFIRE R/ST Panel are presented in Table 12 and Table 13 for NPS and TS specimens, respectively. These tables also indicate the number of specimens with FP results for each co-detection combination, as well as the specific analyte(s) that were discrepant.

Table 12. Multiple Detection Combinations in NPS Specimens as Reported by the SPOTFIRE R/ST Panel

Distinct Co-Detection Combinations				Total Specimens with Co-Detections	Number of Specimens with False Positive Co-Detections	False Positive Analyte(s)
Analyte 1	Analyte 2	Analyte 3	Analyte 4			
Adenovirus	Human metapneumovirus	Human rhinovirus/enterovirus	Parainfluenza virus	1	1	Adenovirus, Parainfluenza virus
Adenovirus	Coronavirus (seasonal)	Human rhinovirus/enterovirus		6	3	Adenovirus (2), Human rhinovirus/enterovirus (1)
Adenovirus	Coronavirus (seasonal)	Respiratory syncytial virus		1	1	Coronavirus (seasonal)
Adenovirus	Human rhinovirus/enterovirus	Parainfluenza virus		1	0	-
Coronavirus SARS-CoV-2	Human rhinovirus/enterovirus	Parainfluenza virus		1	0	-
Coronavirus (seasonal)	Human rhinovirus/enterovirus	Parainfluenza virus		4	4	Coronavirus (seasonal) (1), Human rhinovirus/enterovirus (2), Parainfluenza virus (2)
Human rhinovirus/enterovirus	Parainfluenza virus	Respiratory syncytial virus		2	2	Human rhinovirus/enterovirus (2), Parainfluenza virus (1), Respiratory syncytial virus (1)
Adenovirus	Coronavirus (seasonal)			5	1	Adenovirus
Adenovirus	Human rhinovirus/enterovirus			19	13	Adenovirus (11), Human rhinovirus/enterovirus (5)
Adenovirus	Parainfluenza virus			4	2	Adenovirus (2)
Adenovirus	Respiratory syncytial virus			2	2	Adenovirus (2)
Coronavirus SARS-CoV-2	Human rhinovirus/enterovirus			9	5	Coronavirus SARS-CoV-2 (1), Human rhinovirus/enterovirus (5)
Coronavirus SARS-CoV-2	Respiratory syncytial virus			1	0	-
Coronavirus (seasonal)	Human rhinovirus/enterovirus			11	5	Coronavirus (seasonal) (2), Human rhinovirus/enterovirus (3)
Coronavirus (seasonal)	Parainfluenza virus			9	1	Coronavirus (seasonal)
Coronavirus (seasonal)	Respiratory syncytial virus			1	0	-
Human rhinovirus/enterovirus	Parainfluenza virus			20	8	Human rhinovirus/enterovirus (6), Parainfluenza virus (3)
Human rhinovirus/enterovirus	Respiratory syncytial virus			8	5	Human rhinovirus/enterovirus (4), Respiratory syncytial virus (1)
Parainfluenza virus	Respiratory syncytial virus			1	0	-
Total Co-Detections				106	53	62/229 ^a
Total Double Detections				90	42	47/180
Total Triple Detections				15	10	13/45
Total Quadruple Detections				1	1	2/4

^a Of the 62 discrepant analytes (out of 229 total analytes), 43 (69.4%) were confirmed as being present in the specimen during discrepancy investigation.

Table 13 Multiple Detection Combinations in TS Specimens as Reported by the SPOTFIRE R/ST Panel

Distinct Co-Detection Combinations			Total Specimens with Co-Detections	Number of Specimens with False Positive Co-Detections	False Positive Analyte(s)
Analyte 1	Analyte 2	Analyte 3			
Adenovirus	Coronavirus (seasonal)	<i>S. dysgalactiae</i>	1	0	-
Adenovirus	Coronavirus (seasonal)	<i>S. pyogenes</i>	1	0	-
Adenovirus	Human rhinovirus/enterovirus	Influenza A virus A/H3	1	0	-
Adenovirus	Human rhinovirus/enterovirus	Respiratory syncytial virus	1	0	-
Adenovirus	Human rhinovirus/enterovirus	<i>S. pyogenes</i>	5	2	Adenovirus (1), Human rhinovirus/enterovirus (1), <i>S. pyogenes</i> (1)
Adenovirus	Parainfluenza virus	<i>S. pyogenes</i>	1	1	Parainfluenza Virus
Adenovirus	Respiratory syncytial virus	<i>S. pyogenes</i>	1	1	<i>S. pyogenes</i>
Coronavirus SARS-CoV-2	Human rhinovirus/enterovirus	<i>S. pyogenes</i>	1	0	-
Coronavirus (seasonal)	Parainfluenza virus	<i>S. pyogenes</i>	1	0	-
Human metapneumovirus	Human rhinovirus/enterovirus	<i>S. pyogenes</i>	1	1	Human rhinovirus/enterovirus
Human rhinovirus/enterovirus	Parainfluenza virus	<i>S. pyogenes</i>	2	1	Human rhinovirus/enterovirus, <i>S. pyogenes</i>
Human rhinovirus/enterovirus	Respiratory syncytial virus	<i>S. pyogenes</i>	1	0	-
Adenovirus	Coronavirus (seasonal)		1	0	-
Adenovirus	Human metapneumovirus		1	0	-
Adenovirus	Human rhinovirus/enterovirus		7	3	Adenovirus (2), Human rhinovirus/enterovirus (1)
Adenovirus	Influenza A virus A/H1-2009		1	0	-
Adenovirus	Parainfluenza Virus		3	2	Adenovirus (2)
Adenovirus	<i>S. pyogenes</i>		17	4	Adenovirus (1), <i>S. pyogenes</i> (3)
Coronavirus SARS-CoV-2	Human metapneumovirus		1	0	-
Coronavirus SARS-CoV-2	Human rhinovirus/enterovirus		3	0	-
Coronavirus SARS-CoV-2	Parainfluenza Virus		1	1	Coronavirus SARS-CoV-2
Coronavirus SARS-CoV-2	Respiratory syncytial virus		1	0	-
Coronavirus SARS-CoV-2	<i>S. dysgalactiae</i>		1	1	<i>S. dysgalactiae</i>
Coronavirus SARS-CoV-2	<i>S. pyogenes</i>		8	2	<i>S. pyogenes</i> (2)
Coronavirus (seasonal)	Human rhinovirus/enterovirus		4	2	Coronavirus (seasonal) (1), Human rhinovirus/enterovirus (1)

Distinct Co-Detection Combinations			Total Specimens with Co-Detections	Number of Specimens with False Positive Co-Detections	False Positive Analyte(s)
Analyte 1	Analyte 2	Analyte 3			
Coronavirus (seasonal)	Influenza A virus A/H3		2	2	Coronavirus (seasonal) (2)
Coronavirus (seasonal)	<i>S. pyogenes</i>		7	2	Coronavirus (seasonal) (1), <i>S. pyogenes</i> (1)
Human metapneumovirus	Human rhinovirus/enterovirus		2	1	Human rhinovirus/enterovirus
Human metapneumovirus	<i>S. pyogenes</i>		2	0	-
Human rhinovirus/enterovirus	Influenza A virus A/H3		1	1	Human rhinovirus/enterovirus
Human rhinovirus/enterovirus	Parainfluenza Virus		7	1	Parainfluenza virus
Human rhinovirus/enterovirus	Respiratory syncytial virus		5	1	Respiratory syncytial virus
Human rhinovirus/enterovirus	<i>S. dysgalactiae</i>		5	1	<i>S. dysgalactiae</i>
Human rhinovirus/enterovirus	<i>S. pyogenes</i>		48	16	Human rhinovirus/enterovirus (12), <i>S. pyogenes</i> (4)
Influenza A virus A/H3	Respiratory syncytial virus		1	0	-
Influenza A virus A/H3	<i>S. pyogenes</i>		2	1	<i>S. pyogenes</i>
Influenza B virus	<i>S. pyogenes</i>		2	0	-
Parainfluenza Virus	<i>S. dysgalactiae</i>		1	0	-
Parainfluenza Virus	<i>S. pyogenes</i>		7	2	<i>S. pyogenes</i> (2)
Respiratory syncytial virus	<i>S. dysgalactiae</i>		1	0	-
<i>S. dysgalactiae</i>	<i>S. pyogenes</i>		4	4	<i>S. dysgalactiae</i> (2), <i>S. pyogenes</i> (2)
Total Co-Detections			163	53	55/343 ^b
Total Double Detections			146	47	47/292
Total Triple Detections			17	6	8/51

^a FP *Streptococcus* results based on culture comparator method.

^b Of the 55 discrepant analytes (out of 343 total analytes), 29 (52.7%) were confirmed as being present in the specimen during discrepancy investigation

PERFORMANCE CHARACTERISTICS

Clinical Performance

Prospective Clinical Evaluation

The clinical performance (encompassing both accuracy and ease of use) of the SPOTFIRE R/ST Panel was established during a prospective multi-center study that was further supplemented with archived specimens. Six geographically distinct urgent care or emergency department study sites representative of the intended use setting (five in the US and one in the UK) participated in these studies from December 2020 to September 2021 (NPS and TS specimens) and from September 2022 to May 2023 (TS specimens only). All SPOTFIRE R/ST Panel testing was performed according to the manufacturer's instructions by minimally trained operators. No hands-on training was provided to the SPOTFIRE R/ST Panel test operators; rather, training was limited to written materials (i.e. quick reference guides) that were intended to be included with the SPOTFIRE System.

A total of 1215 NPS specimens and 1165 TS were enrolled from consented volunteers or obtained as residual leftover specimens from subjects of all ages for the prospective clinical study; 95 of these NPS specimens and 288 of these TS specimens were excluded. The most common reasons for specimen exclusion were the SPOTFIRE R/ST Panel pouch was later determined to be expired or the specimen was unable to be tested within the designated timeframe. The final data set consisted of 1120 NPS specimens and 877 TS specimens. Across the six study sites, 352 specimens (259 NPS specimens and 93 TS specimens) were initially collected and immediately frozen for later testing at the source study site. The remaining 1645 specimens (861 NPS specimens and 784 TS specimens) were collected and tested fresh (without freezing). No difference in performance was observed when fresh and frozen specimen results were compared. Therefore, the data collected from 352 valid frozen specimens are combined with data from the valid 1645 fresh specimens for all analyses.

Table 14 provides a summary of demographic information for the specimens included in the prospective specimen studies.

Table 14. Demographic Summary for Prospective SPOTFIRE R/ST Panel Clinical Evaluation

Category		Prospective NPS Specimens	Prospective TS Specimens
Sex	Male	587 (52.4%)	361 (41.2%)
	Female	533 (47.6%)	516 (58.8%)
Age	≤5 years	457 (40.8%)	138 (15.7%)
	6-18 years	258 (23.0%)	456 (52.0%)
	19-40 years	160 (14.3%)	186 (21.2%)
	41-60 years	147 (13.1%)	72 (8.2%)
	61+ years	98 (8.8%)	25 (2.9%)
	Total	1120	877

For most analytes, the performance of the SPOTFIRE R/ST Panel was evaluated by comparing the test results with those from FDA-cleared multiplexed respiratory pathogen panels. The comparator methods for *Streptococcus dysgalactiae* (group C/G Strep), and *Streptococcus pyogenes* (group A Strep) were two analytically validated PCR assays followed by bidirectional sequencing (these analytes were not detected by the FDA-cleared panels), and culture isolation of beta-hemolytic streptococci followed by molecular species identification using PCR confirmed by bidirectional sequencing (performed on cultured isolates).

The performance for prospective and archived studies is summarized in Table 15 and Table 16 for NPS and TS specimens, respectively. Sensitivity or positive percent agreement (PPA) for each analyte was calculated as $100\% \times (TP / (TP + FN))$. True positive (TP) indicates that both the SPOTFIRE R/ST Panel and the comparator method had a positive result for the

specific analyte, and false negative (FN) indicates that the SPOTFIRE R/ST Panel was negative while the comparator result was positive. Specificity or negative percent agreement (NPA) was calculated as $100\% \times (TN / (TN + FP))$. True negative (TN) indicates that both the SPOTFIRE R/ST Panel and the comparator method had negative results, and false positive (FP) indicates that the SPOTFIRE R/ST Panel was positive while the comparator method was negative. The exact binomial two-sided 95% confidence interval (95%CI) was calculated. Investigations of discrepant results are summarized in the footnotes.

Table 15. SPOTFIRE R/ST Panel Clinical Performance Summary for NPS Specimens

SPOTFIRE R/ST Panel R Menu Result	Positive Percent Agreement			Negative Percent Agreement		
	TP/(TP + FN)	%	95%CI	TN/(TN + FP)	%	95%CI
Viruses						
Adenovirus ^a	32/33	97.0	84.7-99.5%	1058/1082	97.8	96.7-98.5%
Coronavirus SARS-CoV-2 ^b	71/73	97.3	90.5-99.2%	1031/1037	99.4	98.7-99.7%
Coronavirus (seasonal) ^c	101/102	99.0	94.7-99.8%	1000/1013	98.7	97.8-99.2%
Human metapneumovirus	1/1	100	-	1114/1114	100	99.7-100%
Human rhinovirus/enterovirus ^d	345/348	99.1	97.5-99.7%	695/767	90.6	88.3-92.5%
Influenza A virus	0/0	-	-	1115/1115	100	99.7-100%
Influenza A virus A/H1-2009	0/0	-	-	1115/1115	100	99.7-100%
Influenza A virus A/H3	0/0	-	-	1115/1115	100	99.7-100%
Influenza B virus	0/0	-	-	1110/1110	100	99.7-100%
Parainfluenza virus ^e	96/98	98.0	92.9-99.4%	1006/1017	98.9	98.1-99.4%
Respiratory syncytial virus ^f	26/27	96.3	81.7-99.3%	1086/1088	99.8	99.3-99.9%
Bacteria						
<i>Bordetella parapertussis</i>	0/0	-	-	1110/1110	100	99.7-100%
<i>Bordetella pertussis</i>	0/0	-	-	1115/1115	100	99.7-100%
<i>Chlamydia pneumoniae</i>	0/0	-	-	1115/1115	100	99.7-100%
<i>Mycoplasma pneumoniae</i>	0/0	-	-	1115/1115	100	99.7-100%

^a Adenovirus was not detected in the single FN specimens upon SPOTFIRE R/ST Panel retest. Adenovirus was detected in 21/24 FP specimens using an additional molecular method.

^b SARS-CoV-2 was detected in 1/2 FN specimens upon SPOTFIRE R/ST Panel retest. SARS-CoV-2 was detected in 2/6 FP specimens using an additional molecular method.

^c Coronavirus (seasonal) was detected in the single FN specimen upon SPOTFIRE R/ST Panel retest. Coronavirus (seasonal) was detected in 8/13 FP specimens using an additional molecular method.

^d Human rhinovirus/enterovirus was detected in 1/3 FN specimens upon SPOTFIRE R/ST Panel retest. Human rhinovirus/enterovirus was detected in 48/72 FP specimens using an additional molecular method.

^e Parainfluenza virus was detected in both FN specimens during discrepancy investigation upon SPOTFIRE R/ST Panel retest. Parainfluenza virus was detected in 9/11 FP specimens using an additional molecular method.

^f Respiratory syncytial virus was detected in the single FN specimen upon SPOTFIRE R/ST Panel retest. Respiratory syncytial virus was detected in 1/2 FP specimens using an additional molecular method.

Table 16. SPOTFIRE R/ST Panel Clinical Performance Summary for TS Specimens^a

SPOTFIRE R/ST Panel ST Menu Result	Positive Percent Agreement			Negative Percent Agreement		
	TP/(TP + FN)	%	95%CI	TN/(TN + FP)	%	95%CI
Viruses						
Adenovirus ^b	60/65	92.3	83.2-96.7%	798/810	98.5	97.4-99.2%
Coronavirus SARS-CoV-2 ^c	43/46	93.5	82.5-97.8%	827/830	99.6	98.9-99.9%
Coronavirus (seasonal) ^d	41/42	97.6	87.7-99.6%	826/833	99.2	98.3-99.6%
Human metapneumovirus ^e	15/17	88.2	65.7-96.7%	858/858	100	99.6-100%
Human rhinovirus/enterovirus ^f	202/213	94.8	91.0-97.1%	619/662	93.5	91.4-95.1%
Influenza A virus	35/35	100	90.1-100%	840/840	100	99.5-100%
Influenza A virus A/H1-2009	12/12	100	75.8-100%	863/863	100	99.6-100%
Influenza A virus A/H3	23/23	100	85.7-100%	852/852	100	99.6-100%
Influenza B virus	4/4	100	51.0-100%	872/872	100	99.6-100%
Parainfluenza virus ^g	50/52	96.2	87.0-98.9%	821/823	99.8	99.1-99.9%
Respiratory syncytial virus ^h	21/24	87.5	69.0-95.7%	849/851	99.8	99.1-99.9%

SPOTFIRE R/ST Panel ST Menu Result		Positive Percent Agreement			Negative Percent Agreement		
		TP/(TP + FN)	%	95%CI	TN/(TN + FP)	%	95%CI
Bacteria							
<i>Chlamydia pneumoniae</i>		0/0	-	-	875/875	100	99.6-100%
<i>Mycoplasma pneumoniae</i>		0/0	-	-	875/875	100	99.6-100%
<i>Streptococcus dysgalactiae</i> (group C/G Strep)	PCR ⁱ	30/30	100	88.6-100%	843/847	99.5	98.8-99.8%
	Culture ^j	26/26	100	87.1-100%	837/843	99.3	98.5-99.7%
<i>Streptococcus pyogenes</i> (group A Strep)	PCR ^k	209/217	96.3	92.9-98.1%	654/660	99.1	98.0-99.6%
	Culture ^l	174/177	98.3	95.1-99.4%	654/692	94.5	92.6-96.0%

^a The performance measures of sensitivity and specificity only refer to the *Streptococcus* analytes for which culture was used as the reference method. Performance measures of PPA and NPA refer to all other analytes, for which molecular assays were used as comparator methods.

^b Adenovirus was detected in 3/5 FN specimens using an additional molecular method. Adenovirus was detected in 4/12 FP specimens using an additional molecular method.

^c SARS-CoV-2 was detected in 1/3 FN specimens upon SPOTFIRE R/ST Panel retest. The three FP specimens were negative for SARS-CoV-2 when tested with additional molecular methods.

^d Coronavirus (seasonal) was detected in the single FN specimen using an additional molecular method. Coronavirus (seasonal) was detected in the 4/7 FP specimens using an additional molecular method.

^e Human metapneumovirus was detected in 1/2 FN specimens upon SPOTFIRE R/ST Panel retest.

^f Human rhinovirus/enterovirus was detected in 7/11 FN specimens during discrepancy investigation: four using an additional molecular method and three upon SPOTFIRE R/ST Panel retest. Human rhinovirus/enterovirus was detected in 14/43 FP specimens using an additional molecular method.

^g Parainfluenza virus was detected in 1/2 FN specimens upon SPOTFIRE R/ST Panel retest. Parainfluenza virus was detected in 1/2 FP specimens using an additional molecular method.

^h Respiratory syncytial virus was detected in all three FN specimen upon SPOTFIRE R/ST Panel retest. Respiratory syncytial virus was detected in 1/2 FP specimens using an additional molecular method.

ⁱ *S. dysgalactiae* was detected in 3/4 FP specimens using an additional molecular method.

^j *S. dysgalactiae* was detected in 5/6 FP specimens using an additional molecular method.

^k *S. pyogenes* was detected in 7/8 FN specimens during discrepancy investigation: four using an additional molecular method and three upon SPOTFIRE R/ST Panel retest. *S. pyogenes* was detected in 2/6 FP specimens using an additional molecular method.

^l *S. pyogenes* was detected in all three FN specimens during discrepancy investigation: one using an additional molecular method and two upon SPOTFIRE R/ST Panel retest. *S. pyogenes* was detected in 34/38 FP specimens using an additional molecular method.

The overall success rate for initial specimen tests was 96.4% (2044/2120). Fourteen (14) tests (14/2120; 0.7%) did not complete on the initial test attempt, resulting in an instrument success rate of 99.3% (2106/2120) for initial specimen tests. Retests were not possible due to insufficient specimen volume. Of the 2106 tests that successfully produced a completed run on the initial test, 2044 had valid internal pouch controls. This represents a 97.1% (2044/2106) success rate for internal pouch controls in completed runs in the initial specimen tests.

Testing of Preselected Archived Specimens

A number of analytes on the SPOTFIRE R/ST Panel were of low prevalence during the prospective study and were not encountered in large enough numbers to adequately demonstrate system performance. To supplement the results of the prospective clinical study, an evaluation of preselected archived retrospective NPS and TS specimens was performed.

A total of 562 frozen archived NPS specimens and 136 frozen archived TS specimens were obtained from 15 external laboratories world-wide and retrospectively tested at the four US clinical sites. Of these, 542 NPS specimens and 128 TS specimens had valid results that were included in performance analysis. The analyte composition of the archived specimens was confirmed using the same comparator methods as the prospective study (described above) for the analyte result to be included in the performance analysis.

The specimens were randomized such that the users performing both the confirmation and the SPOTFIRE R/ST Panel testing were blinded to the expected test result. A summary of the available demographic information of the tested specimens is provided in Table 17, and the results of the SPOTFIRE R/ST Panel performance for these archived NPS and TS specimens are shown in Table 18 and Table 19.

Table 17. Demographic Summary for Valid Archived NPS Specimens

Category		Archived NPS Specimens	Archived TS Specimens
Sex	Male	254 (46.9%)	59 (46.1%)
	Female	185 (34.1%)	69 (53.9%)
	Unknown	103 (19.0%)	0 (0%)
Age	≤5 years	234 (43.2%)	22 (17.2%)
	6-18 years	98 (18.1%)	85 (66.4%)
	19-40 years	36 (6.6%)	16 (12.5%)
	41-60 years	35 (6.5%)	5 (3.9%)
	61+ years	39 (7.2%)	0 (0%)
	Unknown	100 (18.5%)	0 (0%)
	Total	542	128

Table 18. SPOTFIRE R/ST Panel Archived Performance Summary for NPS Specimens

Analyte	Positive Percent Agreement			Negative Percent Agreement		
	TP/(TP + FN)	%	95% CI	TN/(TN + FP)	%	95% CI
Viruses						
Adenovirus ^a	31/31	100	89.0-100%	439/453	96.9	94.9-98.2%
Coronavirus SARS-CoV-2	0/0	-	-	0/0	-	-
Coronavirus (seasonal) ^b	95/96	99.0	94.3-99.8%	381/388	98.2	96.3-99.1%
Human metapneumovirus ^c	32/33	97.0	84.7-99.5%	451/451	100	99.2-100%
Human rhinovirus/enterovirus ^d	29/30	96.7	83.3-99.4%	439/454	96.7	94.6-98.0%
Influenza A virus ^e	58/59	98.3	91.0-99.7%	423/423	100	99.1-100%
Influenza A virus A/H1-2009 ^e	31/32	96.9	84.3-99.4%	450/450	100	99.2-100%
Influenza A virus A/H3	27/27	100	87.5-100%	455/455	100	99.2-100%
Influenza B virus	30/30	100	88.6-100%	28/28	100	87.9-100%
Parainfluenza virus ^f	116/118	98.3	94.0-99.5%	359/366	98.1	96.1-99.1%
Respiratory syncytial virus ^g	37/37	100	90.6-100%	440/447	98.4	96.8-99.2%
Bacteria						
<i>Bordetella parapertussis</i> ^h	24/25	96.0	80.5-99.3%	33/33	100	89.6-100%
<i>Bordetella pertussis</i> ⁱ	27/28	96.4	82.3-99.4%	452/456	99.1	97.8-99.7%
<i>Chlamydia pneumoniae</i> ^j	30/30	100	88.6-100%	452/454	99.6	98.4-99.9%
<i>Mycoplasma pneumoniae</i> ^k	33/33	100	89.6-100%	446/451	98.9	97.4-99.5%

^a Adenovirus was detected in 6/12 FP specimens during discrepancy investigation: one was detected by standard of care and five were detected using an additional molecular method; two additional FP specimens were unable to be investigated.

^b The single FN specimen was unable to be investigated. Seasonal coronavirus (229E/HKU1/NL63/OC43) was detected in 3/6 FP specimens during discrepancy investigation using an additional molecular method; one additional FP specimen was unable to be investigated.

^c Human metapneumovirus was detected in the single FN specimen by standard of care.

^d The single FN specimen was unable to be investigated. Human rhinovirus/enterovirus was detected in 4/14 FP specimens during discrepancy investigation using an additional molecular method; one additional FP specimen was unable to be investigated.

^e Influenza A virus A/H1-2009 was detected in the single FN specimen by standard of care.

^f Parainfluenza virus was detected in the single FN specimen by standard of care; one additional FN specimen was unable to be investigated. Parainfluenza virus was detected in all seven FP specimens during discrepancy investigation: six were detected by standard of care and one was detected using an additional molecular method.

^g Respiratory syncytial virus was detected in 4/6 FP specimens during discrepancy investigation using an additional molecular method; one additional FP specimen was unable to be investigated.

^h *Bordetella parapertussis* was detected in the single FN specimen by standard of care.

ⁱ *Bordetella pertussis* was detected in the single FN specimen by standard of care. *Bordetella pertussis* was detected in 3/4 FP specimens by standard of care.

^j *Chlamydia pneumoniae* was detected in both FP specimens by standard of care.

^k *Mycoplasma pneumoniae* was detected in all five FP specimens during discrepancy investigation: four were detected by standard of care and one was detected using an additional molecular method.

Table 19. SPOTFIRE R/ST Panel Archived Performance Summary for TS Specimens

Analyte	Positive Percent Agreement			Negative Percent Agreement		
	TP/(TP + FN)	%	95% CI	TN/(TN + FP)	%	95% CI
Viruses						
Adenovirus ^a	9/11	81.8	52.3-94.9%	44/48	91.7	80.4-96.7%
Coronavirus SARS-CoV-2	0/0	-	-	0/0	-	-
Coronavirus (seasonal)	12/12	100	75.8-100%	47/47	100	92.4-100%
Human metapneumovirus	4/4	100	51.0-100%	55/55	100	93.5-100%
Human rhinovirus/enterovirus ^b	2/2	100	34.2-100%	55/57	96.5	88.1-99.0%
Influenza A virus	11/11	100	74.1-100%	44/44	100	92.0-100%
Influenza A virus A/H1-2009	7/7	100	64.6-100%	48/48	100	92.6-100%
Influenza A virus A/H3	4/4	100	51.0-100%	51/51	100	93.0-100%
Influenza B virus	20/20	100	83.9-100%	0/0	-	-
Parainfluenza virus	1/1	100	-	58/58	100	93.8-100%
Respiratory syncytial virus ^c	2/2	100	34.2-100%	56/57	98.2	90.7-99.7%
Bacteria						
<i>Chlamydia pneumoniae</i>	2/2	100	34.2-100%	57/57	100	93.7-100%
<i>Mycoplasma pneumoniae</i>	4/4	100	51.0-100%	55/55	100	93.5-100%
<i>Streptococcus dysgalactiae</i> (group C/G Strep)	0/0	-	-	0/0	-	-
<i>Streptococcus pyogenes</i> (group A Strep) ^d	38/39	97.4	86.8-99.5%	10/10	100	72.2-100%

^a Adenovirus was detected in 2/2 FN specimens during discrepancy investigation by standard of care. Adenovirus was detected in 1/2 FP specimens during discrepancy investigation using standard of care; two additional FP specimens were unable to be investigated.

^b Human rhinovirus/enterovirus virus was detected in the single FP specimen by standard of care; one additional FP specimen was unable to be investigated.

^c Respiratory syncytial virus was detected in the single FP specimen by standard of care.

^d *Streptococcus pyogenes* was detected in the single FN by standard of care.

Contrived Testing

Several analytes were not observed in TS specimens in sufficient numbers to demonstrate SPOTFIRE R/ST Panel performance in the prospective and archived specimen studies. Therefore, contrived specimens (N=431; at least 50 for each analyte) were made from unique, analyte-negative clinical TS specimens that were spiked with a variety of different isolates/strains for each organism at concentrations that spanned the detection range of each assay. At least half (50%) of the contrived specimens had analyte concentrations at 2 × the limit of detection (LoD). Contrived specimens were randomized and coded along with 29 negative (unspiked) specimens such that the analyte status of each specimen was unknown to the users performing the testing. The coded contrived specimens were distributed to prospective clinical study sites for testing.

The results of the SPOTFIRE R/ST Panel performance for contrived TS specimens are shown in Table 20.

Table 20. SPOTFIRE R/ST Panel Contrived TS Specimen Performance Summary

Analyte	PPA			NPA		
	TP/(TP + FN)	%	95% CI	TN/(TN + FP)	%	95% CI
Viruses						
Adenovirus	50/50	100	92.9-100%	380/381	99.7	98.5-100%
Coronavirus (seasonal)	129/130	99.2	95.8-99.9%	301/301	100	98.7-100%
Human metapneumovirus	45/49	91.8	80.8-96.8%	382/382	100	99.0-100%
Influenza A virus	93/93	100	96.0-100%	332/332	100	98.9-100%
Influenza A virus A/H1-2009	46/46	100	92.3-100%	379/379	100	99.0-100%
Influenza A virus A/H3	47/47	100	92.4-100%	378/378	100	99.0-100%
Influenza B virus	47/49	95.9	86.3-98.9%	333/333	100	98.9-100%
Parainfluenza virus	122/128	95.3	90.2-97.8%	303/303	100	98.7-100%
Respiratory syncytial virus	49/50	98.0	89.5-99.6%	381/381	100	99.0-100%
Bacteria						
<i>Chlamydia pneumoniae</i>	50/50	100	92.9-100%	381/381	100	99.0-100%

Analyte	PPA			NPA		
	TP/(TP + FN)	%	95% CI	TN/(TN + FP)	%	95% CI
<i>Mycoplasma pneumoniae</i>	48/50	96.0	86.5-98.9%	381/381	100	99.0-100%

ANALYTICAL PERFORMANCE CHARACTERISTICS

Limit of Detection

The limit of detection (LoD) for SPOTFIRE R/ST Panel analytes was estimated by testing dilutions of contrived samples containing known concentrations of organism in both VTM and Amies media. The LoD concentrations were confirmed by testing at least 20 replicates at the estimated LoD. Confirmation of LoD required detection in at least 95% of replicates tested. The confirmed LoD concentrations for the SPOTFIRE R/ST Panel are listed in Table 21.

Table 21. Limit of Detection (LoD) for the SPOTFIRE R/ST Panel Analytes

Analyte	Isolate Source ID	LoD Concentration ^a
Viruses		
Adenovirus	Species A Serotype 31 ZeptoMetrix 0810073CF	8.2E-03 TCID ₅₀ /mL ^b (2.0E+02 copies/mL)
	Species B Serotype 3 ZeptoMetrix 0810062CF	8.0E-01 TCID ₅₀ /mL (8.4E+02 copies/mL)
	Species C Serotype 2 WHO I.S. NIBSC 16-324	8.2E+02 IU/mL ^c (8.2E+02 copies/mL ^c)
	Species D Serotype 37 ZeptoMetrix 0810119CF	1.1E-02 TCID ₅₀ /mL (4.5E+02 copies/mL)
	Species E Serotype 4 ZeptoMetrix 0810070CF	1.8E-02 TCID ₅₀ /mL (1.0E+04 copies/mL)
	Species F Serotype 41 ZeptoMetrix 0810085CF	1.4E-02 TCID ₅₀ /mL (1.0E+02 copies/mL)
Coronavirus (seasonal)	229E ATCC VR-740	6.5E-01 TCID ₅₀ /mL (1.1E+01 copies/mL)
	HKU1 Clinical Specimen ^d	1.8E+04 copies/mL
	OC43 ZeptoMetrix 0810024CF	1.6E-02 TCID ₅₀ /mL (6.3E+01 copies/mL)
	NL63 ZeptoMetrix 0810228CF	2.5E-03 TCID ₅₀ /mL (4.7E+00 copies/mL)
Coronavirus SARS-CoV-2	USA-WA1/2020 (heat inactivated) ATCC VR-1986HK	1.1E-01 TCID ₅₀ /mL (2.5E+02 copies/mL)
	1 st WHO International Standard England/02/2020 NIBSC 20/146	2.5E+02 IU/mL ^b
Human metapneumovirus	A1-16 Iowa 10/2003 ZeptoMetrix 0810161CF	3.2E+00 TCID ₅₀ /mL (2.4E+02 copies/mL)
	B1-3 Peru2-2002 ZeptoMetrix 0810156CF	2.5E-01 TCID ₅₀ /mL (5.4E+02 copies/mL)
	A2-27 Iowa A/2004 ZeptoMetrix 0810164CF	5.8E-01 TCID ₅₀ /mL (1.8E+03 copies/mL)
	B2-18 IA18-2003 ZeptoMetrix 0810162CF	2.0E+00 TCID ₅₀ /mL (7.7E+02 copies/mL)
Human rhinovirus/enterovirus	Human Rhinovirus 1A ZeptoMetrix 0810012CFN	2.1E-01 TCID ₅₀ /mL (1.1E+00 copies/mL)
	Enterovirus D68 US/MO/14-18947 ATCC VR-1823	1.1E+01 TCID ₅₀ /mL (5.4E+01 copies/mL)
Influenza A virus Subtype H1-2009	Influenza A H1N1 pdm A/Michigan/45/15 ZeptoMetrix 0810538CF	8.2E-01 TCID ₅₀ /mL (3.4E+02 copies/mL)

Analyte	Isolate Source ID	LoD Concentration ^a
Influenza A virus Subtype H3	Influenza A H3N2 A/Hong Kong/4801/14 ZeptoMetrix 0810526CF	8.6E-01 TCID ₅₀ /mL (3.4E+02 copies/mL)
Influenza B virus	B/Florida/02/06 (Victoria Lineage) ZeptoMetrix 0810037CF	3.3E-02 TCID ₅₀ /mL (1.6E+02 copies/mL)
	B/Nevada/03/2011 (Victoria Lineage) BEI NR-44023	1.6E+00 CEID ₅₀ /mL (4.3E+00 copies/mL)
	B/Florida/04/06 (Yamagata Lineage) ZeptoMetrix 0810255CF	4.0E-01 TCID ₅₀ /mL (3.2E+01 copies/mL)
Parainfluenza virus	Serotype 1 ZeptoMetrix 0810014CF	4.6E+00 TCID ₅₀ /mL ^e (1.4E+03 copies/mL)
	Serotype 2 ZeptoMetrix 0810015CF	1.4E+01 TCID ₅₀ /mL (1.6E+02 copies/mL)
	Serotype 3 ZeptoMetrix 0810016CF	2.6E+01 TCID ₅₀ /mL ^e (6.1E+01 copies/mL)
	Serotype 4 ZeptoMetrix 0810060CF	2.0E+02 TCID ₅₀ /mL ^e (1.1E+03 copies/mL)
Respiratory syncytial virus	Type A 2006 ZeptoMetrix 0810040ACF	6.2E-02 TCID ₅₀ /mL (2.2E+01 copies/mL)
	Type B 3/2015 Isolate #1 ZeptoMetrix 0810479CF	2.8E-02 TCID ₅₀ /mL (2.4E+01 copies/mL)
Bacteria		
<i>Chlamydia pneumoniae</i>	AR-39 ATCC 53592	2.0E+01 IFU/mL (1.4E+02 copies/mL)
<i>Mycoplasma pneumoniae</i>	M129 ZeptoMetrix 0801579	1.0E+01 CCU/mL (2.1E+03 copies/mL)
Bacteria (Respiratory only^f)		
<i>Bordetella parapertussis</i>	E595 ZeptoMetrix 0801462	4.0E+01 CFU/mL
<i>Bordetella pertussis</i>	A639 ZeptoMetrix 0801459	3.3E+02 CFU/mL (3.8E+02 copies/mL)
Bacteria (Sore Throat only^g)		
<i>Streptococcus dysgalactiae</i> (Group C/G Strep)	ssp. <i>equisimilis</i> Z068 ZeptoMetrix 0801516	3.3E+02 CFU/mL (3.4E+03 copies/mL)
<i>Streptococcus pyogenes</i> (Group A Strep)	SF-130 T1 ATCC 12344	4.5E+02 cells/mL (2.9E+03 copies/mL)

^a LoD concentration may vary from what is listed based on the accuracy and precision of the quantification method.

^b LoD confirmation (≥95% detection) was achieved at a two-fold lower level in VTM.

^c IU = International Units; BioFire Diagnostics quantified the WHO International Standard by quantitative real-time PCR to demonstrate that 8.2E+02 IU/mL=8.2E+02 copies/mL

^d Testing for Coronavirus HKU1 utilized a clinical specimen due to the lack of availability of a cultured isolate. Viral concentration was determined in RNA copies/mL by quantitative real-time RT-PCR.

^e LoD confirmation (≥95% detection) was achieved at a 2 to 3-fold lower level in Amies media.

^f Results only reported when the Nasopharyngeal Swab sample type is selected.

^g Results only reported when the Throat Swab sample type is selected.

NOTE: LoD concentrations in copies/mL in Table 21 above are based on extraction of nucleic acids from isolate cultures followed by quantitative real-time PCR (qPCR) or digital PCR (dPCR). The accuracy of concentrations may be affected by extraction efficiency, standard curve accuracy (qPCR only), assay conditions, inhibitors, and/or sequence variance. The quantification has not been compared to a reference material or other quantification methods.

NOTE: LoD concentrations of cultured viruses provided in units of TCID₅₀ (50% Tissue Culture Infectious Dose) or CEID₅₀ (50% Chicken Embryo Infectious Dose) are not a direct count of viral particles or nucleic acid, but an indirect measure of viral concentration based on infectivity and cytotoxicity. TCID₅₀/mL and CEID₅₀/mL will therefore vary depending on technique and methodology (including cell type, culture media and conditions, cytotoxicity of the virus, etc.). It is not appropriate to make determinations on relative sensitivity of detection for different cultures and/or different molecular assays based on LoD values measured in TCID₅₀/mL or CEID₅₀/mL.

Analytical Reactivity (Inclusivity)

The analytical reactivity (inclusivity) of the SPOTFIRE R/ST Panel assays was assessed by testing over 190 viral and bacterial isolates that represented clinical and genetic diversity and included the available phylogenetic, geographic, and temporal diversity of each SPOTFIRE R/ST Panel analyte. Isolates were tested in triplicate at concentrations near the LoD.

A summary of analytical reactivity is included in Table 22 to Table 35.

NOTE: Influenza A assays will react variably with non-human influenza A viruses and rarely encountered human influenza A viruses that are not H1-2009 or H3; generally producing Uncertain: Influenza A virus or Influenza A virus (No Subtype Found) results.

NOTE: The SPOTFIRE R/ST Panel assays may react with vaccines that contain specific segments of the pathogen genome or full genome or vaccines containing attenuated/inactivated pathogen, including vaccines for SARS-CoV- 2, influenza A (various subtypes), influenza B, poliovirus (Human rhinovirus/enterovirus), and Bordetella pertussis. Care should be taken to minimize contamination of samples with vaccines, and clinical history of vaccine administration should be considered in the interpretation of results, particularly for vaccines administered by nasal spray.

Table 22. Summary of Reactivity to Adenovirus Isolates

Species	Type	Source/Isolate ID	[Strain/Location/Year]	Result
A	12	ATCC VR-863	[Huie/Massachusetts]	Adenovirus Positive
	18	ATCC VR-19	[Washington DC/1954]	
	31	ZeptoMetrix 0810073CF	-	
B	3	ZeptoMetrix 0810062CF	-	
	7	ATCC VR-7	[Gomen/California/1954]	
	7a	ZeptoMetrix 0810021CF	-	
	7d/d2	UIRF	[Iowa/2001]	
	11	ATCC VR-12	[Slobitski/Massachusetts]	
	14	ATCC VR-15	[De Wit/Netherlands/1955]	
	16	ATCC VR-17	[CH.79/Saudi Arabia/1955]	
	21	ATCC VR-1833	[128/Saudi Arabia/1956]	
	34	ATCC VR-716	[Compton/1972]	
	35	ATCC VR-718	[Holden]	
C	50	ATCC VR-1602	[Wan, RIVM no. 88-1773]	
	1	ZeptoMetrix 0810050CF	-	
	2	W.H.O NIBSC 16/324	-	
	5	ZeptoMetrix 0810020CF	-	
D	6	ATCC VR-6	[Tonsil 99]	
	8	ZeptoMetrix 0810069CF	-	
	20	ZeptoMetrix 0810115CF	[KB]	
E	37	ZeptoMetrix 0810119CF	-	
	4	ZeptoMetrix 0810070CF	-	
	4	ATCC VR-1572	[RI-67/Missouri/1952-1953]	
F	4a	UIRF	[S.Carolina/2004]	
	40	NCPV 0101141v	-	
	40	ZeptoMetrix 0810084CF	-	
	41	ATCC VR-930	[Tak (73-3544)]	
	41	ZeptoMetrix 0810085CF	[Tak]	

Table 23. Summary of Reactivity to Bordetella parapertussis Isolates

Source/Isolate ID	Strain/Location/Year	Result
ZeptoMetrix 0801462	[E595]	Bordetella parapertussis Positive
ATCC 9305	[517]	
ATCC 53892	[PT28G]	
ATCC 53893	[PT 26/28G]	
ATCC 15237	[NCTC 10853]	
ATCC 15311	[NCTC 5952]	
ATCC BAA-587	[12822/Germany/1993]	
ZeptoMetrix 0801461	[A747]	
ZeptoMetrix 0801643	[C510]	

Source/Isolate ID	Strain/Location/Year	Result
ZeptoMetrix 0801644	[E838]	

Table 24. Summary of Reactivity to *Bordetella pertussis* Isolates

Source/Isolate ID	Strain/Location/Year	Result
ZeptoMetrix 0801459	[A639]	<i>Bordetella pertussis</i> Positive
ATCC 10380	[10-536]	
ATCC 51445	[CNCTC Hp 12/63 [623]]	
ATCC 8467	[F]	
ATCC 9340	[5 [17921]	
ATCC 9797	[18323 [NCTC 10739]	
ATCC BAA-1335	[MN2531]	
ATCC BAA-589	[Tohama]	
ZeptoMetrix 0801460	[E431]	

Table 25. Summary of Reactivity to *Chlamydia pneumoniae* Isolates

Source/Isolate ID	Strain/Location/Year	Result
ATCC 53592	[AR-39]	<i>Chlamydia pneumoniae</i> Positive
ATCC VR-1310	[CWL-029]	
ATCC VR-1356	[TWAR strain 2023]	
ATCC VR-1360	[CM-1]	
ATCC VR-1435	[J-21]	
ATCC VR-1452	[A03]	
ATCC VR-2282	[TWAR strain, TW-183/Taiwan/1965]	

Table 26. Summary of Reactivity to Coronavirus Isolates

Type	Source/Isolate ID	Strain/Location/Year	Result
NL63	ZeptoMetrix 0810228CF	-	Coronavirus (seasonal) Positive
	BEI NR-470	[Amsterdam/2003]	
229E	ATCC VR-740	-	
	ZeptoMetrix 0810229CF	-	
HKU1	Clinical Specimen	[Columbus OH, 2016]	
	Clinical Specimen	[South Carolina/2010]	Coronavirus SARS-CoV-2 Positive
	Clinical Specimen	[France/2016]	
	Clinical Specimen	[France/2016]	
OC43	ZeptoMetrix 0810024CF	-	
	ATCC VR-759	-	
SARS-CoV-2^a	ATCC VR-1986HK	[USA-WA1/2020]	
	ATCC VR-1991D	[Hong Kong/VM20001061/2020]	
	ATCC VR-1992D	[2019-nCoV/Italy-INMI1]	
	ATCC VR-1994D	[Germany/BavPat1/2020]	
	ATCC VR-3326D	[USA/CA_CDC_5574/2020]	
	BEI NR-52499 ^b	[England/02/2020]	
	BEI NR-52501 ^c	[Singapore/2/2020]	
	BEI NR-52503 ^d	[USA-IL1/2020]	
	BEI NR-52505 ^e	[USA-AZ1/2020]	
	BEI NR-52507 ^f	[USA-CA3/2020]	
	BEI NR-52510 ^g	[Chile/Santiago_op4d1/2020]	
	BEI NR-53518 ^h	[New York-PV08410/2020]	
	LGC SeraCare AccuPlex™ 0505-0298 ⁱ	[Omicron B.1.1.529 Variant]	

^a See **Table 36** for additional SARS-CoV-2 reactivity predictions based on *in silico* analysis.

^b The following reagent was deposited by Professor Maria Zambon and obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate England/02/2020, NR-52499.

^c The following reagent was contributed by Duke-National University of Singapore, Programme in Emerging Infectious Diseases for distribution through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate Singapore/2/2020, NR-52501.

^d The following reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate USA-IL1/2020, NR-52503.

^e The following reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate USA-AZ1/2020, NR-52505.

^f The following reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate USA-CA3/2020, NR-52507.

^g The following reagent was obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate Chile/Santiago_op4d1/2020, NR-52510.

^h The following reagent was obtained through BEI Resources, NIAID, NIH: Genomic RNA from SARS-Related Coronavirus 2, Isolate New York-PV08410/2020, NR-53518.

ⁱ Recombinant alphavirus that contains the full SARS-CoV-2 genome with mutations identified in the S and N gene of the SARS-CoV-2 Omicron variant B.1.1.529.

Table 27. Summary of Reactivity to Human Rhinovirus and Enterovirus Isolates

Species	Serotype	Source/Isolate ID	[Strain/Location/Year]	Result
Human Rhinovirus				
A	1	ZeptoMetrix 0810012CFN	[1A]	Human rhinovirus/enterovirus Positive
	77	ATCC VR-1187	[130-63]	
	85	ATCC VR-1195	[50-525-CV54]	
	34	ATCC VR-1365	[137-3]	
	57	ATCC VR-1600	[Ch47]	
	7	ATCC VR-1601	[68-CV11]	
	16	ATCC VR-283	[11757]	
	2	ATCC VR-482	[HGP]	
B	17	ATCC VR-1663	[33342]	
	14	ATCC VR-284	[1059]	
	42	ATCC VR-1950	[56822]	
	3	ATCC VR-483	[FEB]	
	27	ATCC VR-1137	[5870]	
	83	ATCC VR-1193	[Baylor 7]	
Enterovirus				
A	Enterovirus 71	ATCC VR-1432	[71 H]	Human rhinovirus/enterovirus Positive
	Coxsackievirus 10	ATCC VR-168	[NY/1950]	
B	Coxsackievirus 9	ZeptoMetrix 0810017CF	-	
	Echovirus 11	ZeptoMetrix 0810023CF	-	
	Coxsackievirus B3	ZeptoMetrix 0810074CF	-	
	Coxsackievirus B4	ZeptoMetrix 0810075CF	-	
	Echovirus 6	ZeptoMetrix 0810076CF	-	
	Echovirus 9	ZeptoMetrix 0810077CF	-	
C	Coxsackievirus A24	ATCC VR-583	[DN-19/TX/1963]	
	Coxsackievirus A21	ATCC VR-850	[Kuykendall/CA/1952]	
D	Enterovirus D68	ATCC VR-1823	[US/MO/14-18947]	

Table 28. Summary of Reactivity to Human Metapneumovirus Isolates

Genotype	Serotype	Source/Isolate ID	[Location/Year]	Result
A1	9	ZeptoMetrix 0810160CF	[Iowa 3/2002]	Human metapneumovirus Positive
	16	ZeptoMetrix 0810161CF	[Iowa 10/2003]	
A2	20	ZeptoMetrix 0810163CF	[Iowa 14/2003]	
	27	ZeptoMetrix 0810164CF	[Iowa 27/2004]	
B1	3	ZeptoMetrix 0810156CF	[Peru2-2002]	
	5	ZeptoMetrix 0810158CF	[Peru 3/2003]	
B2	4	ZeptoMetrix 0810157CF	[Peru 1/2002]	
	8	ZeptoMetrix 0810159CF	[Peru 6/2003]	
	18	ZeptoMetrix 0810162CF	[IA18-2003]	
	Unknown	BEI NR-22232	[TN/91-316]	

Table 29. Summary of Reactivity to Influenza A Isolates

Type	Host	Source/Isolate ID	Strain/Location/Year	Result
H1N1pdm09	Human	ZeptoMetrix 0810538CF	[Michigan/45/15]	Influenza A virus Subtype H1-2009 Positive
		BEI NR-19823	[Netherlands/2629/2009]	
		BEI NR-42938	[Georgia/F32551/2012]	
		BEI NR-44345	[Hong Kong/H090-761-V1(0)/2009]	
		ZeptoMetrix 0810109CFJ	[Canada/6294/2009]	
		ZeptoMetrix 0810165CF	[California/07/2009]	
		ZeptoMetrix 0810166CF	[Mexico/4108/2009]	
H3N2	Human	ZeptoMetrix 0810249CF	[SwineNY/03/2009]	Influenza A virus Subtype H3 Positive
		ATCC VR-544	[Hong Kong/8/1968]	
		ATCC VR-547	[Aichi/2/1968]	
		ATCC VR-776	[Alice]	
		ATCC VR-810	[Port Chalmers/1/1973]	
		ATCC VR-822	[Victoria/3/1975]	
		ZeptoMetrix 0810138CF	[Brisbane/10/2007]	
		ZeptoMetrix 0810238CF	[Texas/50/2012]	
		ZeptoMetrix 0810252CF	[Wisconsin/67/2005]	
H1N1	Human	ZeptoMetrix 0810526CF	[Hong Kong/4801/14]	Influenza A virus Positive (No Subtype Found)
		ZeptoMetrix 0810036CFN	[New Caledonia/20/1999]	
		ZeptoMetrix 0810244CF	[Solomon Islands/3/2006]	
		ZeptoMetrix 0810244CF	[Brisbane/59/2007]	

Type	Host	Source/Isolate ID	Strain/Location/Year	Result
	Swine	ATCC VR-333	[Swine/Iowa/15/1930]	
		ATCC VR-897	[Swine/A/New Jersey/8/76]	
		ATCC VR-99	[Swine/1976/1931]	
H2N2	Human	BEI NR-2775 ^a	[A/Japan/305/1957]	Uncertain: Influenza A virus
H5N3	Avian	BEI NR-9682 ^b	[A/Duck/Singapore/645/97]	
H1N2	Recombinant	BEI NR-3478 ^c	[Kilbourne F63 A/NWS/1934 (HA) x A/Rockefeller Institute/5/1957 (NA)]	
H10N7	Avian	BEI NR-2765 ^d	[A/Chicken/Germany/N/49]	

^a The following reagent was obtained through BEI Resources, NIAID, NIH: Genomic RNA from Influenza A virus, A/Japan/305/1957 (H2N2), NR-2775.

^b The following reagent was obtained through BEI Resources, NIAID, NIH: Genomic RNA from Influenza A virus, A/duck/Singapore/645/1997 (H5N3), Wild Type, NR-9682.

^c The following reagent was obtained through BEI Resources, NIAID, NIH: Kilbourne F63: A/NWS/1934 (HA) x A/Rockefeller Institute/5/1957 (NA) (H1N2), Reassortant NWS-F, NR-3478.

^d The following reagent was obtained through BEI Resources, NIAID, NIH: Genomic RNA from Influenza A virus, A/chicken/Germany/N/1949 (H10N7), NR-2765.

Table 30. Summary of Reactivity to Influenza B Isolates

Lineage	Source/Isolate ID	[Strain/Location/Year]	Result
Yamagata	ZeptoMetrix 0810255CF	[Florida/04/06]	Influenza B virus Positive
	ZeptoMetrix 0810239CF	[2/Massachusetts/2012]	
	ZeptoMetrix 0810241CF	[1/Wisconsin/2010]	
	ZeptoMetrix 0810256CF	[07/Florida/2004]	
Victoria	ZeptoMetrix 0810037CF	[B/Florida/02/06]	
	BEI NR-44023	[B/Nevada/03/2011]	
	ATCC VR-823	[5/Hong Kong/1972]	
	CDC 2005743348	[1/Ohio/2005]	
Unknown	ZeptoMetrix 0810258CF	[2506/Malaysia/2004]	
	ATCC VR-101	[Lee/1940]	
	ATCC VR-102	[Allen/1945]	
	ATCC VR-103	[GL/1739/1954]	
	ATCC VR-295	[2/Taiwan/1962]	
	ATCC VR-296	[1/Maryland/1959]	
	ATCC VR-786	[Brigit/Russia/1969]	

Table 31. Summary of Reactivity to *Mycoplasma pneumoniae* Isolates

Source/Isolate ID	Strain/Location/Year	Result
ZeptoMetrix 0801579	[M129]	<i>Mycoplasma pneumoniae</i> Positive
ATCC 29085	[PI 1428]	
ATCC 29342	[M129-B7]	
ATCC 15492	[Mac]	
ATCC 15531-TTR	[FH strain of Eaton Agent [NCTC 10119]	
ATCC 15293	[M52]	
ATCC 15377	[Bru]	
ATCC 39505	[Mutant 22]	
ATCC 49894	[UTMP-10P]	

Table 32. Summary of Reactivity to Parainfluenza Virus Isolates

Serotype	Subtype	Source/Isolate ID	Strain/Location/Year	Result
1		ZeptoMetrix 0810014CF	-	Parainfluenza virus Positive
		ATCC VR-94	[C-35/1957]	
		BEI NR-48680	[FRA/29221106/2009]	
2		ZeptoMetrix 0810015CF	-	
		ATCC VR-92	[Greer/1955]	
3		ZeptoMetrix 0810016CF	-	
		ATCC VR-93	[C-243/1957]	
		BEI NR-3233	[NIH 47885 Wash/47885/57]	
4	A	ZeptoMetrix 0810060CF	-	
		ATCC VR-1378	[M-25/1958]	
	B	ZeptoMetrix 0810060CF	-	
		ATCC VR-1377	[CH-19503/1962]	
		ZeptoMetrix 0810060BCF	-	

Table 33. Summary of Reactivity to Respiratory Syncytial Virus Isolates

Type	Source/Isolate ID	Strain/Location/Year	Result
A	ZeptoMetrix 0810040ACF	[2006]	Respiratory syncytial virus Positive
	ATCC VR-26	[Long/Maryland/1956]	
	ATCC VR-1540	[A2/Melbourne/1961]	

Type	Source/Isolate ID	Strain/Location/Year	Result
B	ZeptoMetrix 0810474CF	[2/2015 Isolate #2]	
	ZeptoMetrix 0810452CF	[12/2014 Isolate #2]	
	ZeptoMetrix 0810479CF	[3/2015 Isolate #1]	
	ZeptoMetrix 0810040CF	[Ch-93 (18)-18]	
	ATCC VR-1400	[WV/14617/1985]	
	ATCC VR-955	[9320/Massachusetts/1977]	
	ATCC VR-1580	[18537/WashingtonDC/1962]	
	ZeptoMetrix 0810451CF	[11/2014 Isolate #2]	

Table 34. Summary of Reactivity to *Streptococcus dysgalactiae* ssp. *equisimilis* Isolates

Source/Isolate ID	Strain/Location/Year	Result
ZeptoMetrix 0801516	[Z068]	<i>Streptococcus dysgalactiae</i> (Group C/G Strep) Positive
ATCC 12388	[C74]	
NCTC 8543	[655 Chestle]	
CCUG 45841	-	
CCUG 45898	-	
ATCC 12394	[D166B]	
ATCC 9542	[1180]	
ATCC BAA-337	[AC-2074]	
ATCC 10009 ^a	[13-166]	<i>Streptococcus dysgalactiae</i> (Group C/G Strep) Negative

^a A deletion in the gene target was identified that prevents amplification/detection of this isolate and other *S. dysgalactiae* isolates of animal origin.

Table 35. Summary of Reactivity to *Streptococcus pyogenes* Isolates

Source ID	[Strain/Location/Year]	Result
ATCC 12344	[SF-130 T-type 1]	<i>Streptococcus pyogenes</i> (Group A Strep) Positive
ATCC BAA-947	[MGAS 5005 M-type 1]	
ATCC 700294	[SF370 M1 GAS M-type 1]	
ATCC 12384	[C203 T-type 3]	
ATCC BAA-595	[MGAS 315, M type 3]	
ATCC 51500	[DLS 88002, Weller], M type 3]	
NCTC 8193	[T5 B, M type 5]	
ATCC 12348	[S43 T-type 6]	
ATCC BAA-1065	[MGAS 2096 [A374], M type 12]	
ATCC 12356	[Typing strain J17E [A. Coburn R9], M type 17]	
ATCC BAA-572	[MGAS 8232, M type 18]	
ATCC 8133	[Typing strain T23 [F. Griffith strain Barts 102], M type 23]	
ATCC 12360	[Typing strain J17F [A. Coburn R17], M type 26]	
NCTC 8310	[Coggins, M type 29]	
NCTC 8195	[Quinn (D24/46), M type 30]	
NCTC 8229	[C 95/12B, M type 39]	
ATCC 12372	[Typing strain C143 [C143], M type 40]	
ATCC 12377	[C105 [20RS14], M type 46]	
NCTC 10880	[R66/3489 (378), M type 62]	
ATCC 19615	[Bruno]	
ATCC 49399	[AC A62]	
ATCC 700466	CDC-SS-872 [R67/3884]	
ATCC 25663	[P20080]	
ATCC 49117	[397]	
ATCC 21060 ^a	[Su]	<i>Streptococcus pyogenes</i> (Group A Strep) Negative

^a A deletion in the gene target was identified that prevents amplification/detection of this isolate.

In Silico Reactivity Predictions for Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) Assays

Evaluation of analytical reactivity for the SPOTFIRE R/ST Panel SARS-CoV-2 assays (SARS-CoV2-1 and SARS-CoV2-2) was based on *in silico* analysis of all available sequences in the GISAID database as of December 31 2023.

This analysis determined that the >99.99% of 13,682,523 sequences will be detected by one or both SPOTFIRE R/ST Panel SARS-CoV-2 assays based on homology and mismatch location with one or both sets of primers. A limitation on detection (both assays impaired) is predicted for less than 0.006% of the sequences evaluated (824 / 13,682,523 (Table 36).

The sequences evaluated include lineages and variants of concern (VOC) or variants under investigation (VUI) that may have important epidemiological, immunological, or pathogenic properties from a public health perspective, such as Delta and Omicron variants. Variants evaluated are listed in the *BioFire® Respiratory Panels SARS-CoV-2 Reactivity Tech Note* technical note at www.biofire.com/support/documents.

All lineages and variants of public health interest identified as of December 2023 are predicted to be detected; new sequences and variants will continue to be monitored for impacts on detection by the SPOTFIRE R/ST Panel assays.

Table 36. In silico Prediction of SARS-CoV-2 Detection by SPOTFIRE R Panel Assays

+/+ indicates detected by both assays with no impairment, +/- indicates detection by one assay with no impairment and potential for impaired detection by the other assay, -/- indicates potential for impaired detection by both assays

In silico prediction				
Predicted Assay Result		SARS-CoV2-1		Total Sequences
Number of Sequences		+	-	
SARS-CoV2-2	+	13,423,817	203,572	13,681,699/13,682,523
	-	54,310	824 ^a	99.994%

^a 196 unique sequences out of 824 total sequences

Periodic updates to the *in silico* analysis are performed based on the most currently available GISAID sequences and can also be accessed in the *BioFire® Respiratory Panels SARS-CoV-2 Reactivity Tech Note* at www.biofire.com/support/documents.

Analytical Specificity (Cross-Reactivity & Exclusivity)

The potential for non-specific amplification and detection by the SPOTFIRE R/ST Panel assays was evaluated by *in silico* analysis of available sequences and by testing of high concentrations of on-panel and off-panel organisms. The organisms evaluated included relevant bacteria, fungi, and viruses that are either phylogenetically related to organisms detected by the SPOTFIRE R/ST Panel or pathogenic/commensal organisms that may be present in NPS or TS specimens. Each organism was tested in triplicate at the highest possible concentration (generally $\geq 1.0 \times 10^7$ units/mL for bacteria and $\geq 1.0 \times 10^5$ units/mL for viruses).

In silico analysis and testing identified a risk of SARS-CoV-2 assay cross-reactivity with a few sequences of SARS-like viruses isolated from bats and pangolin as well as intra-panel cross-reactivity with *Bordetella* species and Influenza A subtypes of swine origin. A summary of potential cross-reactivity is provided in Table 37. The on-panel and off-panel isolates and concentrations tested are listed in Table 38 and Table 39, respectively.

Table 37. Predicted and Observed Cross-Reactivity of the SPOTFIRE R/ST Panel

Cross-Reactive Organism/Sequence	SPOTFIRE R/ST Panel Analyte Result	Description
Bat coronavirus Pangolin coronavirus Bat SARS-like coronavirus	Coronavirus SARS-CoV-2 ^a	The SARS-CoV-2 assays can amplify select sequences from closely related sarbecoviruses isolated from bats and pangolin.
<i>Bordetella bronchiseptica</i> (isolates containing IS 1001 gene)	<i>Bordetella parapertussis</i>	Some strains of <i>B. bronchiseptica</i> contain IS 1001 insertion sequences identical to those present in <i>B. parapertussis</i> . These sequences will be amplified by

		the IS1001 assay and reported as <i>Bordetella parapertussis</i> .
<i>Bordetella bronchiseptica</i>^a <i>Bordetella parapertussis</i>^b <i>Bordetella pertussis</i>	Human rhinovirus/enterovirus^c	The HRV/EV assay designed to target Human Rhinovirus/Enterovirus may amplify the oxidoreductase gene from <i>Bordetella</i> species (<i>B. pertussis</i> , <i>B. parapertussis</i> , and <i>B. bronchiseptica</i>) when organisms are present at a high concentration. Cross-reactivity with <i>B. pertussis</i> was observed at a concentration of $\geq 1.3 \times 10^9$ CFU/mL.
Bovine picornavirus Canine picornavirus	Human rhinovirus/enterovirus^a	Bovine and canine picornaviruses may be detected and reported as Human rhinovirus/enterovirus when present at high concentration.
Influenza A H1N1 (swine origin)	Influenza A Subtype H1-2009^a	The FluA-H1-2009 assay may react with the H1 hemagglutinin gene sequences from viruses of swine origin. The SPOTFIRE R/ST Panel will report either Influenza A virus (Subtype H1-2009) or Influenza A virus (no subtype identified), depending on the strain and concentration in the sample.
<i>Chlamydia gallinacea</i>	<i>Chlamydia pneumoniae</i>^a	<i>Chlamydia gallinacea</i> may be detected and reported as <i>Chlamydia pneumoniae</i> when present at high concentration.

^a Indicated cross-reactivity is predicted based on *in silico* analysis.

^b Cross-reactivity between the HRV/EV assay and *B. bronchiseptica* and *B. parapertussis* is predicted based on *in silico* analysis but was not observed when testing organisms at the highest possible concentrations (8.3×10^9 cells/mL for *B. bronchiseptica* and 4.6×10^9 CFU/mL for *B. parapertussis*).

^c Cross-reactivity between the HRV/EV assay and *B. parapertussis* or *B. pertussis* will have positive results reported for Human rhinovirus/enterovirus and *B. parapertussis* or *B. pertussis*, respectively. Cross-reactivity between the HRV/EV assay and *B. bronchiseptica* isolates that do not carry the IS1001 sequence will only have positive results reported for Human rhinovirus/enterovirus. Cross-reactivity between the HRV/EV assay and *B. bronchiseptica* isolates that contain the IS1001 sequence will have positive results reported for Human rhinovirus/enterovirus and *B. parapertussis*.

Table 38. Summary of Results for On-Panel Organisms Tested During Evaluation of Analytical Specificity of the SPOTFIRE R/ST Panel

Organism	Isolate ID	Concentration Tested	Observed Cross-Reactivity
Bacteria			
<i>Bordetella parapertussis</i>	ZeptoMetrix 0801462	4.6E+09 CFU/mL	None
<i>Bordetella pertussis</i>	ZeptoMetrix 0801459	1.3E+10 CFU/mL	Human rhinovirus/enterovirus^a
		1.3E+09 CFU/mL	None
<i>Chlamydia pneumoniae</i>	ATCC 53592	2.9E+07 IFU/mL	None
<i>Mycoplasma pneumoniae</i>	ZeptoMetrix 0801579	2.5E+07 CCU/mL	None
<i>Streptococcus dysgalactiae</i> ssp. <i>equisimilis</i> (isolated from human) ^b	ZeptoMetrix 0801516	7.8E+08 CFU/mL	None
<i>Streptococcus pyogenes</i>	ATCC 12344	3.4E+08 cells/mL	None
Viruses			
Adenovirus A	ZeptoMetrix 0810073CF	1.4E+05 TCID ₅₀ /mL	None
Adenovirus B	ZeptoMetrix 0810062CF	1.2E+07 TCID ₅₀ /mL	None
Adenovirus C	ZeptoMetrix 0810110CF	2.2E+06 TCID ₅₀ /mL	None
Adenovirus D	ZeptoMetrix 0810119CF	1.7E+05 TCID ₅₀ /mL	None
Adenovirus E	ZeptoMetrix 0810070CF	1.4E+05 TCID ₅₀ /mL	None
Adenovirus F	ZeptoMetrix 0810085CF	1.1E+06 TCID ₅₀ /mL	None
Coronavirus 229E	ATCC VR-740	8.9E+06 TCID ₅₀ /mL	None
Coronavirus HKU1	Clinical Specimens	4.5E+07 copies/mL	None
Coronavirus NL63	ZeptoMetrix 0810228CF	5.0E+05 TCID ₅₀ /mL	None
Coronavirus OC43	ZeptoMetrix 0810024CF	3.6E+05 TCID ₅₀ /mL	None
Coronavirus SARS-CoV-2 (heat-inactivated)	ATCC VR-1986HK	7.6E+07 copies/mL	None
Enterovirus D68	ATCC VR-1823	1.6E+07 TCID ₅₀ /mL	None
Human metapneumovirus A1	ZeptoMetrix 0810161CF	2.5E+05 TCID ₅₀ /mL	None
Human metapneumovirus A2	ZeptoMetrix 0810164CF	3.6E+05 TCID ₅₀ /mL	None
Human metapneumovirus B1	ZeptoMetrix 0810156CF	1.6E+04 TCID ₅₀ /mL	None
Human metapneumovirus B2	ZeptoMetrix 0810162CF	1.3E+06 TCID ₅₀ /mL	None
Human rhinovirus A1	ZeptoMetrix 0810012CFN	1.3E+06 TCID ₅₀ /mL	None
Influenza A H1N1pdm09	ZeptoMetrix 0810538CF	1.4E+05 TCID ₅₀ /mL	None
Influenza A H3N2	ZeptoMetrix 0810526CF	7.2E+05 TCID ₅₀ /mL	None
Influenza B (Victoria Lineage)	BEI NR-44023	2.8E+08 CEID ₅₀ /mL	None
	ZeptoMetrix 0810037CF	2.5E+05 TCID ₅₀ /mL	None
Influenza B (Yamagata Lineage)	ZeptoMetrix 0810256CF	2.1E+04 TCID ₅₀ /mL	None
Parainfluenza virus 1	ZeptoMetrix 0810014CF	4.2E+05 TCID ₅₀ /mL	None
Parainfluenza virus 2	ZeptoMetrix 0810015CF	1.2E+07 TCID ₅₀ /mL	None
Parainfluenza virus 3	ZeptoMetrix 0810016CF	3.4E+07 TCID ₅₀ /mL	None
Parainfluenza virus 4	ZeptoMetrix 0810060CF	3.4E+07 TCID ₅₀ /mL	None
Respiratory syncytial virus A	ZeptoMetrix 0810040ACF	4.2E+05 TCID ₅₀ /mL	None

Respiratory syncytial virus B	ZeptoMetrix 0810479CF	4.2E+05 TCID ₅₀ /mL	None
-------------------------------	-----------------------	--------------------------------	------

^a The HRV/EV assay may amplify off-target sequences found in strains of *Bordetella* species (*B. pertussis*, *B. parapertussis*, and *B. bronchiseptica*) when present at a concentration $\geq 1.3 \times 10^4$ CFU/mL.

^b The Sdysgalactiae assay was designed to target a genomic region present in *Streptococcus dysgalactiae* ssp. *equisimilis* isolates of human origin.

Table 39. Summary of Results for Off-Panel Organisms Tested During Evaluation of Analytical Specificity of the SPOTFIRE R/ST Panel

Organism	Isolate ID	Concentration Tested	Observed Cross-Reactivity
Bacteria			
<i>Arcanobacterium bernardiae</i>	ATCC BAA-441	1.6E+09 cells/mL	None
<i>Arcanobacterium haemolyticum</i>	ATCC 9345	1.5E+08 cells/mL	None
<i>Arcanobacterium pyogenes</i>	ATCC 49698	6.7E+09 cells/mL	None
<i>Bacillus cereus</i>	ATCC 7064	8.3E+09 cells/mL	None
<i>Bordetella bronchiseptica</i>	ATCC 10580	8.3E+09 cells/mL	None
	ATCC 4617	7.9E+09 cells/mL	None
	ATCC 19395	7.9E+09 cells/mL	None
	NRRL B-59914	7.1E+09 cells/mL	None
	NRRL B-59909 ^a	2.8E+01 cells/mL ^b	<i>Bordetella parapertussis</i>^a
<i>Bordetella holmesii</i>	ATCC 700052	2.8E+00 cells/mL ^b	None
<i>Burkholderia cepacia</i>	ATCC 51671	8.3E+09 cells/mL	None
<i>Campylobacter rectus</i>	ATCC 33238	7.9E+09 cells/mL	None
<i>Chlamydia trachomatis</i>	ZeptoMetrix 0801775	7.6E+07 cells/mL	None
<i>Corynebacterium diphtheriae</i>	ATCC 27010	1.3E+08 IFU/mL	None
<i>Corynebacterium pseudodiphtheriticum</i>	ATCC 10700	8.0E+09 cells/mL	None
<i>Enterococcus casseliflavus</i>	ATCC 49605	8.7E+09 cells/mL	None
<i>Enterococcus faecalis</i>	ZeptoMetrix 0801637	8.0E+09 cells/mL	None
<i>Escherichia coli</i>	ATCC BAA-2196	8.0E+09 CFU/mL	None
<i>Fusobacterium necrophorum</i> ssp. <i>funduliforme</i>	ATCC 51357	7.2E+09 cells/mL	None
<i>Fusobacterium nucleatum</i>	ATCC 25586	4.4E+08 cells/mL	None
<i>Fusobacterium varium</i>	ATCC 27725	4.9E+08 cells/mL	None
<i>Gemella haemolysans</i>	ATCC 10379	1.6E+08 cells/mL	None
<i>Gemella morbillorum</i>	ATCC 27824	4.0E+09 cells/mL	None
<i>Granulicatella adiacens</i>	ATCC 49175	1.0E+08 cells/mL	None
<i>Haemophilus influenzae</i>	ATCC 10211	1.3E+09 cells/mL	None
<i>Haemophilus parahaemolyticus</i>	ATCC 49700	8.3E+09 cells/mL	None
<i>Klebsiella pneumoniae</i>	CDC AR#0115	8.7E+09 cells/mL	None
<i>Lactobacillus rhamnosus</i>	ATCC 7469	7.3E+09 CFU/mL	None
<i>Lactococcus lactis</i>	ATCC 29146	7.9E+09 cells/mL	None
<i>Legionella pneumophila</i>	ATCC 33215	6.2E+09 cells/mL	None
<i>Leptotrichia buccalis</i>	ATCC 14201	7.0E+09 cells/mL	None
<i>Moraxella catarrhalis</i>	ATCC 43627	4.4E+08 cells/mL	None
<i>Mycobacterium tuberculosis</i>	ZeptoMetrix 0801660	7.2E+09 cells/mL	None
<i>Mycoplasma buccale</i>	Mycoplasma Experience NC10136	6.1E+06 CFU/mL	None
<i>Mycoplasma faucium</i>	Mycoplasma Experience NC10174	1.4E+07 CFU/mL	None
<i>Mycoplasma fermentans</i>	Mycoplasma Experience NC10117	1.4E+06 CFU/mL	None
<i>Mycoplasma genitalium</i>	Mycoplasma Experience NC10195	2.8E+07 CFU/mL	None
<i>Mycoplasma hominis</i>	Mycoplasma Experience NC10111	1.8E+06 CFU/mL	None
<i>Mycoplasma lipophilum</i>	Mycoplasma Experience NC10173	1.2E+07 CFU/ml	None
<i>Mycoplasma orale</i>	Mycoplasma Experience NC10112	1.5E+06 CFU/mL	None
<i>Mycoplasma salivarium</i>	Mycoplasma Experience NC10113	2.2E+07 CFU/mL	None
<i>Neisseria elongata</i>	ATCC 25295	4.4E+06 CFU/mL	None
<i>Neisseria gonorrhoeae</i>	ZeptoMetrix 0801482	8.5E+09 cells/mL	None
<i>Neisseria lactamica</i>	ATCC 23971	4.9E+07 CFU/mL	None
<i>Neisseria meningitidis</i>	ATCC 13113	2.7E+09 cells/mL	None
<i>Neisseria sicca</i>	ATCC 9913	7.4E+09 cells/mL	None
<i>Neisseria subflava</i>	ATCC 49275	7.2E+09 cells/mL	None
<i>Parvimonas micra</i> ^c	ATCC 33270	8.0E+09 cells/mL	None
<i>Pneumocystis carinii</i>	ATCC PRA-159	6.0E+07 cells/mL	None
<i>Porphyromonas endodontalis</i>	ATCC 35406	1.0E+07 nuclei/mL	None
<i>Porphyromonas gingivalis</i>	ATCC BAA-308	1.6E+07 cells/mL	None
<i>Prevotella histicola</i>	BEI HM-471	5.0E+08 cells/mL	None
<i>Prevotella melaninogenica</i>	ATCC 25845	9.0E+08 cell/mL	None
<i>Prevotella oralis</i>	ATCC 33322	6.9E+08 cells/mL	None
<i>Pseudomonas aeruginosa</i>	CDC AR#0092	6.2E+08 cells/mL	None
<i>Rhodococcus equi</i>	ATCC 33706	8.3E+09 cells/mL	None
<i>Serratia marcescens</i>	ATCC 27137	7.3E+09 cells/mL	None
<i>Staphylococcus aureus</i>	ATCC BAA-1700	8.9E+09 cells/mL	None

Organism	Isolate ID	Concentration Tested	Observed Cross-Reactivity
<i>Staphylococcus epidermidis</i>	ATCC 12228	8.0E+09 cells/mL	None
<i>Staphylococcus haemolyticus</i>	ATCC 29968	8.0E+09 cells/mL	None
<i>Staphylococcus intermedius</i>	ATCC 29663	8.2E+09 cells/mL	None
<i>Staphylococcus saprophyticus</i>	ATCC 15305	8.1E+09 cells/mL	None
<i>Streptococcus agalactiae</i>	ATCC 13813	6.0E+09 cells/mL	None
<i>Streptococcus anginosus</i>	ATCC 700231	7.1E+09 cells/mL	None
<i>Streptococcus constellatus</i> ssp. <i>pharyngis</i>	NCTC 13122	5.6E+08 cells/mL	None
<i>Streptococcus dysgalactiae</i> ssp. <i>dysgalactiae</i> ^d	ATCC 43078	6.7E+09 cells/mL	None
	NCTC 4669	7.4E+09 cells/mL	None
	NCTC 4335	8.4E+09 cells/mL	None
	NCTC 4670	6.6E+09 cells/mL	None
	CCUG 27665	7.4E+09 cells/mL	None
	CCUG 28112	6.7E+09 cells/mL	None
	CCUG 28114	7.5E+09 cells/mL	None
<i>Streptococcus dysgalactiae</i> ssp. <i>equisimilis</i> (isolated from pig) ^e	CCUG 28117	7.1E+09 cells/mL	None
<i>Streptococcus dysgalactiae</i> ssp. <i>equisimilis</i> (isolated from horse) ^e	CCUG 27664	7.5E+09 cells/mL	None
	ATCC 10009	6.9E+09 cells/mL	None
<i>Streptococcus gallolyticus</i>	ATCC 43143	2.8E+09 cells/mL	None
<i>Streptococcus gordonii</i>	ATCC 10558	4.5E+09 cells/mL	None
<i>Streptococcus intermedius</i>	ATCC 27335	2.9E+09 cells/mL	None
<i>Streptococcus mitis</i>	ATCC 15914	3.2E+09 cells/mL	None
<i>Streptococcus mutans</i>	ATCC 25175	2.3E+09 cells/mL	None
<i>Streptococcus oralis</i>	ATCC 10557	1.1E+09 cells/mL	None
<i>Streptococcus parasanguinis</i>	ATCC 15912	7.8E+09 cells/mL	None
<i>Streptococcus pneumoniae</i>	ATCC 49619	2.5E+08 cells/mL	None
<i>Streptococcus salivarius</i>	ATCC 13419	6.6E+09 cells/mL	None
<i>Streptococcus sanguinis</i>	ATCC 10556	1.1E+09 cells/mL	None
<i>Tannerella forsythia</i>	ATCC BAA-2717	2.6E+08 cells/mL	None
<i>Treponema denticola</i>	ATCC 33520	2.2E+08 cells/mL	None
<i>Ureaplasma urealyticum</i>	ATCC 27618	5.7E+07 cells/mL	None
<i>Veillonella parvula</i>	ATCC 10790	4.7E+08 cells/mL	None
Fungi			
<i>Candida albicans</i>	ATCC MYA-2876	2.8E+08 cells/mL	None
<i>Saccharomyces cerevisiae</i>	ATCC 18824	1.9E+08 cells/mL	None
Viruses			
Cytomegalovirus	ZeptoMetrix 0810003CF	1.9E+05 TCID ₅₀ /mL	None
Epstein-Barr virus	ZeptoMetrix 0810008CF	5.9E+06 copies/mL	None
Human herpes simplex virus 1	ATCC VR-260	8.9E+06 TCID ₅₀ /mL	None
Measles virus	ZeptoMetrix 0810025CF	2.5E+05 TCID ₅₀ /mL	None
Middle East respiratory syndrome coronavirus (heat-inactivated)	ZeptoMetrix 0810575CFHI	1.2E+05 TCID ₅₀ /mL	None
Mumps virus	ZeptoMetrix 0810079CF	2.0E+06 TCID ₅₀ /mL	None
Severe acute respiratory syndrome coronavirus (purified genomic RNA)	BEI NR-52346	5.3E+05 genomes/mL	None

^a *Bordetella bronchiseptica* strain NRRL B-59909 contains genomic insertions of the IS1001 gene and is reactive with the SPOTFIRE R/ST Panel *Bordetella parapertussis* IS1001 assay.

^b Positive results were observed for *Bordetella parapertussis* in all replicates tested at concentrations $\geq 2.8 \times 10^1$ cells/mL and one out of three replicates tested at 2.8×10^0 cells/mL.

^c *Parvimonas micra* was formerly classified as *Micromonas micros* and *Peptostreptococcus micros*.

^d The Sdysgalactiae assay was designed to target a genomic region present in *Streptococcus dysgalactiae* ssp. *equisimilis* isolates of human origin. The *Streptococcus dysgalactiae* ssp. *dysgalactiae* isolates listed in the table do not contain the genomic region targeted by the Sdysgalactiae assay; no cross-reactivity was observed.

^e The Sdysgalactiae assay was designed to target a genomic region present in *Streptococcus dysgalactiae* ssp. *equisimilis* isolates of human origin. The *Streptococcus dysgalactiae* ssp. *equisimilis* isolates of animal origin listed in the table do not contain the region targeted by the Sdysgalactiae assay; no cross-reactivity was observed.

Near-LoD/Reproducibility Evaluation

A near-LoD/reproducibility evaluation was performed to demonstrate that the SPOTFIRE R/ST Panel could reproducibly provide accurate results for weak-positive and negative samples when used by minimally trained operators. Contrived samples were tested at three of the prospective clinical study sites and additionally on three unique SPOTFIRE Systems at BioFire Diagnostics (BioFire) by trained BioFire personnel. The contrived samples contained combinations of SPOTFIRE R/ST Panel analytes prepared at or near (1× to 3×) the LoD. For testing performed at clinical sites, samples were tested over five testing events (non-consecutive days) by two operators during the course of their normal workday routine. Each site was equipped with a single SPOTFIRE System. Testing at all three sites was performed with a single reagent lot. For each testing event, each operator ran two replicate pouches for a total of 20 replicates per site and 60 total replicates across all three sites. For testing performed at BioFire, samples were tested over five consecutive days, by two operators per system, using three different reagent lots. Each day of testing, the two operators each tested three replicates on each system for a total of 30 replicates per system and 90 total replicates across all systems. When combined, each analyte was tested in a total of 150 replicates by at least 12 different operators across six different SPOTFIRE Systems.

A summary of results (percent (%) agreement with the expected positive or negative result) for each analyte (by site and system) is provided in Table 40. The SPOTFIRE R/ST Panel reported the expected positive results for panel analytes in 94% -100% of samples and the expected negative results for all analytes in 100% of samples. Comparison of the positive percent agreement between user groups (98.8% for trained operators at BioFire versus 98.9% for minimally trained operators) demonstrates that the accuracy of the SPOTFIRE R/ST Panel is not dependent upon the specific expertise of the user.

Table 40. Reproducibility of Results on the SPOTFIRE R/ST Panel and SPOTFIRE System

Analyte Isolate (Source ID)		Concentration Tested (test level)	Expected Result	SPOTFIRE System testing								All Sites /Systems [95% Confidence Interval]
				BioFire Dx				Clinical				
				System A	System B	System C	Total	Site 1	Site 2	Site 3	Total	
Adenovirus Species B (ZeptoMetrix 0810062CF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		2.4E+00 TCID ₅₀ /mL (3x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Bordetella parapertussis (ZeptoMetrix 0801462)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		1.2E+02 CFU/mL (3x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Bordetella pertussis (ZeptoMetrix 0801459)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		9.9E+02 CFU/mL (3x LoD)	Positive	29/30 (96.7%)	29/30 (96.7%)	30/30 (100%)	88/90 (97.8%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	148/150 98.7% [95.3-99.8%]
Chlamydia pneumoniae (ATCC 53592)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		2.0E+01 IFU/mL (1x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Coronavirus (seasonal)	No Analyte		Negative	60/60 (100%)	60/60 (100%)	60/60 (100%)	180/180 (100%)	40/40 (100%)	40/40 (100%)	40/40 (100%)	120/120 (100%)	300/300 100% [98.8-100%]
	229E (ATCC VR-740)	2.0E+00 TCID ₅₀ /mL (3x LoD)	Positive	30/30 (100%)	30/30 (100%)	29/30 (96.7%)	89/90 (98.9%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	149/150 99.3% [96.3-100%]
	OC43 (ZeptoMetrix 0810024CF)	1.6E-02 TCID ₅₀ /mL (1x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
	NL63 (ZeptoMetrix 0810228CF)	2.5E-03 TCID ₅₀ /mL (1x LoD)	Positive	27/30 (90.0%)	30/30 (100%)	30/30 (100%)	87/90 (96.7%)	18/20 (90.0%)	20/20 (100%)	18/20 (90.0%)	56/60 (93.3%)	143/150 95.3% [90.6-98.1%]
Coronavirus SARS-CoV-2		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]

Analyte Isolate (Source ID)		Concentration Tested (test level)	Expected Result	SPOTFIRE System testing								All Sites /Systems [95% Confidence Interval]
				BioFire Dx				Clinical				
				System A	System B	System C	Total	Site 1	Site 2	Site 3	Total	
Severe Acute Respiratory Syndrome Coronavirus 2 (ATCC VR-1986HK)		2.5E+02 copies/mL (1× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Human metapneumovirus 3 Type B1 (ZeptoMetrix 0810156CF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		7.5E-01 TCID ₅₀ /mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	19/20 (95.0%)	20/20 (100%)	19/20 (95.0%)	58/60 (96.7%)	148/150 98.7% [95.3-99.8%]
Human rhinovirus/ enterovirus Enterovirus D68 US/MO/14-18947 (ATCC VR-1823)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		1.1E+01 TCID ₅₀ /mL (1× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Influenza A virus	No Analyte		Negative	90/90 (100%)	90/90 (100%)	90/90 (100%)	270/270 (100%)	60/60 (100%)	60/60 (100%)	60/60 (100%)	180/180 (100%)	450/450 100% [99.2-100%]
	Influenza A H1N1pdm (ZeptoMetrix 0810538CF)	2.5E+00 TCID ₅₀ /mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
	Influenza A H3N2 (ZeptoMetrix 0810526CF)	2.6E+00 TCID ₅₀ /mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	19/20 (95.0%)	19/20 (95.0%)	20/20 (100%)	58/60 (96.7%)	148/150 98.7% [95.3-99.8%]
Influenza A virus Subtype H1-2009 Influenza A H1N1pdm (ZeptoMetrix 0810538CF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		2.5E+00 TCID ₅₀ /mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
Influenza A virus Subtype H3 Influenza A H3N2 (ZeptoMetrix 0810526CF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		2.6E+00 TCID ₅₀ /mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	19/20 (95.0%)	19/20 (95.0%)	20/20 (100%)	58/60 (96.7%)	148/150 98.7% [95.3-99.8%]

Analyte Isolate (Source ID)		Concentration Tested (test level)	Expected Result	SPOTFIRE System testing								All Sites /Systems [95% Confidence Interval]
				BioFire Dx				Clinical				
				System A	System B	System C	Total	Site 1	Site 2	Site 3	Total	
Influenza B virus (ZeptoMetrix 0810037CF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		9.9E-02 TCID ₅₀ /mL (3x LoD)	Positive	29/30 (96.7%)	30/30 (100%)	30/30 (100%)	89/90 (98.9%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	149/150 99.3% [96.3-100%]
Mycoplasma pneumoniae (ZeptoMetrix 0801579)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		1.0E+01 CCU/mL (1x LoD)	Positive	29/30 (96.7%)	30/30 (100%)	30/30 (100%)	89/90 (98.9%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	149/150 99.3% [96.3-100%]
Parainfluenza virus	No Analyte		Negative	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
	Parainfluenza virus 1 (ZeptoMetrix 0810014CF)	4.6E+00 TCID ₅₀ /mL (1x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	150/150 100% [97.6-100%]
	Parainfluenza virus 2 (ZeptoMetrix 0810015CF)	4.2E+01 TCID ₅₀ /mL (3x LoD)	Positive	30/30 (100%)	29/30 (96.7%)	28/30 (93.3%)	87/90 (96.7%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	147/150 98.0% [94.3-99.6%]
	Parainfluenza virus 3 (ZeptoMetrix 0810016CF)	8.8E+00 TCID ₅₀ /mL (1x LoD)	Positive	28/30 (93.3%)	30/30 (100%)	30/30 (100%)	88/90 (97.8%)	19/20 (95.0%)	18/20 (90.0%)	20/20 (100%)	57/60 (95.0%)	145/150 96.7% [92.4-98.9%]
	Parainfluenza virus 4 (ZeptoMetrix 0810060CF)	2.0E+02 TCID ₅₀ /mL (1x LoD)	Positive	29/30 (96.7%)	30/30 (100%)	30/30 (100%)	89/90 (98.9%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	149/150 99.3% [96.3-100%]
Respiratory syncytial virus (ZeptoMetrix 0810040ACF)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
		6.2E-02 TCID ₅₀ /mL (1x LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	19/20 (95.0%)	20/20 (100%)	59/60 (98.3%)	149/150 99.3% [96.3-100%]
Streptococcus dysgalactiae (Group C/G Strep)		No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]

Analyte Isolate (Source ID)	Concentration Tested (test level)	Expected Result	SPOTFIRE System testing								All Sites /Systems [95% Confidence Interval]
			BioFire Dx				Clinical				
			System A	System B	System C	Total	Site 1	Site 2	Site 3	Total	
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i> (ZeptoMetrix 0801516)	3.3E+02 CFU/mL (1× LoD)	Positive	30/30 (100%)	29/30 (96.7%)	27/30 (90.0%)	86/90 (95.6%)	20/20 (100%)	20/20 (100%)	20/20 (100%)	60/60 (100%)	146/150 97.3% [93.3-99.3%]
<i>Streptococcus pyogenes</i> (Group A Strep) (ATCC 12344)	No Analyte	Negative	120/120 (100%)	120/120 (100%)	120/120 (100%)	360/360 (100%)	80/80 (100%)	80/80 (100%)	80/80 (100%)	240/240 (100%)	600/600 100% [99.4-100%]
	1.4E+03 cells/mL (3× LoD)	Positive	30/30 (100%)	30/30 (100%)	30/30 (100%)	90/90 (100%)	20/20 (100%)	19/20 (95.0%)	20/20 (100%)	59/60 (98.3%)	149/150 99.3% [96.3-100%]
Total positive agreement (%) by system/user group			681/690 98.7%	687/690 99.6%	684/690 99.1%	2052/2070 99.1%	454/460 98.7%	454/460 98.7%	457/460 99.3%	1365/1380 98.9%	3417/3450 99.0% [98.7-99.3%]
Overall positive agreement (%) [95% Confidence Interval]											

Interference

Potentially interfering substances that could be present in NPS or TS specimens or may be introduced into specimens during collection or subsequent handling and testing were evaluated for their effect on SPOTFIRE R/ST Panel performance. The substances tested included endogenous substances that may be found at normal or elevated levels in clinical specimens (e.g. blood, mucus/mucin, human genomic DNA), various commensal or infectious microorganisms, medications, washes or topical applications for the nasal passage or throat, various swabs and transport media for specimen collection, and substances used to clean, decontaminate, or disinfect work areas.

Each substance was added to contrived samples containing representative organisms at concentrations near (3×) the LoD. The concentration of substance added to the samples was equal to or greater than the highest level expected to be in NPS or TS specimens.

Valid and accurate results were obtained for each sample containing substances and microorganisms at the concentrations listed in Table 41.

Table 41. Substances Tested on the SPOTFIRE R/ST Panel - No Interference Observed


Substance Tested	Concentration Tested
Endogenous Substances	
Human Whole Blood (with Na Citrate)	10% (v/v)
Human Sputum/Mucus	1% (v/v)
Human Genomic DNA	20 ng/μL
Exogenous Substances^a	
Promethazine hydrochloride	1.04 μmol/L (3.34E-04 mg/mL)
Acetaminophen (paracetamol)	1.0E+03 μmol/L (1.5E-01 mg/mL)
Acetylsalicylic acid (Aspirin)	167 μmol/L (3.0E-02 mg/mL)
Ibuprofen	1060 μmol/L (2.2E-01 mg/mL)
Albuterol sulfate (common ingredient in rescue inhalers)	0.188 μmol/L (5.4E-05 mg/mL)
Triple antibiotic ointment (neomycin/polymyxin B/bacitracin)	2% w/v
Mucinex [®] Severe Nasal Congestion Relief Clear & Cool Nasal Spray (Oxymetazoline hydrochloride 0.05%)	1% v/v
Saline nasal spray (sodium chloride 0.65%, disodium phosphate, phenylcarbinol, monosodium phosphate and benzalkonium chloride solution)	1% v/v
Vicks [®] VapoRub [®] Cough Suppressant Topical Analgesic (Camphor 4.8%, eucalyptus oil 1.2%, and menthol 2.6%)	1% w/v
Vaseline [®] Petroleum Jelly (100% white petrolatum)	1% w/v
Orajel [™] (benzalkonium chloride 0.13%, benzocaine 20%, menthol 0.5%, zinc chloride 0.15%)	2% w/v
Chloraseptic [®] Sore Throat Spray (Phenol 1.4%)	1% v/v
Vicks [®] Formula 44 [™] DM (dextromethorphan hydrobromide 0.67 mg/mL, guaifenesin 13 mg/mL) (cough syrup)	1% v/v
Phenylephrine hydrochloride (common ingredient in nasal decongestants)	1% w/v
Nasal spray (fluticasone propionate 50 mcg)	1% v/v
Sucrets [®] Sore Throat (dyclonine hydrochloride 2.0 mg/lozenge)	1% w/v
Benadryl [®] Allergy Liqui-gels [®] (diphenhydramine hydrochloride 25 mg/capsule)	1% v/v
Zicam [®] Cold Remedy (Galphimia Glauca 4x, Luffa Operculata 4x, Sabadilla 4x)	1% v/v
Cold-eeze [®] (zinc gluconate 2.3%)	1% w/v
HALLS lozenge (menthol 5 mg/lozenge)	1% w/v
Listerine [®] Cool Mint [®] (menthol 0.042%, thymol 0.064%, methyl salicylate 0.06%, eucalyptol 0.092%)	6.5% v/v
Copenhagen [®] Snuff (Tobacco)	1% w/v
JUICE HEAD (30% propylene glycol, 70% vegetable glycerin) (e-juice)	1% v/v
Technique-Specific Substances	
Rayon swab (COPAN Diagnostics Inc.)	1 swab
Nylon flocked swab (COPAN Diagnostics Inc.)	1 swab
Polyester swab (COPAN Diagnostics Inc.)	1 swab
Calcium Alginate swab (Puritan [®])	1 swab
Cary-Blair	90% v/v
Dulbecco's Modified Eagles-Medium (DMEM)	90% v/v
MicroTest [™] M4 [™] Viral Transport Media (VTM)	90% v/v
Viral Preservative Media (VPM)	90% v/v
Phosphate Buffered Saline (PBS)	90% v/v
PrimeStore [®] MTM Molecular Transport Media	90% v/v
Stuart Transport Medium	90% v/v
eNAT [™] Molecular Transport Medium	90% v/v
Bleach	1% v/v, 2% v/v ^b
Ethanol	7% v/v
Disinfecting wipes (ammonium chloride)	0.25 – 0.5 inch square/sample
DNAZap [™]	1% v/v
RNaseZap [™]	1% v/v
Competing Microorganisms	
On-Panel	
Coronavirus 229E	1.5E+07 copies/mL
Enterovirus D68	7.8 E+07 copies/mL
Parainfluenza virus 3	8.0E+06 copies/mL
Respiratory syncytial virus A	1.5E+07 copies/mL
Adenovirus A31	1.6E+07 copies/mL
<i>Bordetella pertussis</i>	1.6E+09 copies/mL
<i>Streptococcus pyogenes</i>	2.2E+08 copies/mL
Off-Panel	
Cytomegalovirus (CMV)	4.2E+04 TCID ₅₀ /mL
Herpes Simplex Virus 1	9.0E+06 TCID ₅₀ /mL
<i>Staphylococcus aureus</i>	7.4E+08 CFU/mL
<i>Streptococcus pneumoniae</i>	2.5E+07 CFU/mL
<i>Haemophilus influenzae</i>	8.3E+08 CFU/mL

Substance Tested	Concentration Tested
<i>Candida albicans</i>	2.8E+07 CFU/mL

^a Nasal influenza vaccines (e.g. FluMist®) were not evaluated but are predicted to be reactive with the Influenza A (subtype) and Influenza B assays.

^b Incubation of sample with 1% (v/v) bleach for 15 minutes, 4 hours, or ~18.5-hour (overnight) or 2% (v/v) bleach for 15 minutes did not result in interference.

 **NOTE:** Avoid contact between samples and bleach prior to testing (bleach can damage nucleic acids and prevent amplification and detection by the panel).

 **NOTE:** Compatibility of the SPOTFIRE R/ST Panel with NPS in PrimeStore® MTM has not been evaluated in the intended use setting. PrimeStore® MTM and Sample Buffer contain guanidine salts that will react with bleach to form a toxic gas. Use caution if using bleach for disinfection purposes when collecting or testing NPS or TS specimens.

Carry-Over

When using sensitive molecular tests to analyze patient samples, cross-contamination, or carry-over, of analyte material between samples is a concern. This study evaluated the risk of sample-to-sample carry-over during pouch loading and testing for contrived liquid samples.

To simulate a high-risk scenario, sample loading alternated between samples with high levels of organism (N=5 for each organism; refer to Table 42 for test concentrations) and negative samples. All samples were loaded using the same Pouch Loading Station. Analyte positivity was evaluated for all pouches to determine the risk of false positive results due to sample-to-sample carry-over. No unexpected positive results were observed in this study.






















The data support that sample-to-sample carry-over poses a negligible risk to the accuracy of the SPOTFIRE R/ST Panel test results when the test is used according to the provided instructions.

Table 42: Test Concentration for Evaluation of Carry-over

Organism	Source ID	Concentration Tested (multiple of LoD)
Enterovirus	ATCC VR-1823	1.6E+07 TCID ₅₀ /mL (1,500,000x)
Influenza A H3N2	ZeptoMetrix 0810526CF	7.2E+04 TCID ₅₀ /mL (84,000x)
<i>Chlamydia pneumoniae</i>	ATCC 53592	2.9E+06 IFU/mL (150,000x)
Coronavirus SARS-CoV-2	ATCC VR-1986HK	8.4E+02 TCID ₅₀ /mL (7,600x)

APPENDIX A

Symbols Glossary

ISO 15223-1					
Medical devices - Symbols to be used with medical devices labels, labeling and information to be supplied					
5.1.1 	Manufacturer	5.1.2 	Authorized representative in the European Community	5.1.4 	Use-By date (YYYY-MM-DD)
5.1.5 	Batch Code (Lot Number)	5.1.6 	Catalog Number	5.1.7 	Serial Number
5.2.8 	Do Not Use if Package Is Damaged	5.3.2 	Keep Away from Sunlight	5.3.7 	Temperature Limit
5.4.2 	Do Not Reuse	5.4.3 	Consult Instructions for Use	5.5.1 	In vitro Diagnostic Medical Device
5.5.5 	Contains Sufficient For <n> Tests		5.7.10 	Unique Device Identifier	
Use of Symbols in Labeling – 81 FR 38911, Docket No. (FDA-2013-N-0125)					
Rx Only	Prescription Use Only				
United Nations Globally Harmonized System of Classification and Labeling of chemicals (GHS) (ST/SG/AC.10/30)					
	Serious eye damage, Category. 1		Acute toxicity, oral, Category. 4 & Skin corrosion, irritation, Category 2		Acute aquatic hazard, Category 1 & Long-term aquatic hazard, Category 1
European Union In Vitro Diagnostic Directive (IVDD 98/79/EC) and European In Vitro Diagnostic Regulation (IVDR 2017/746)					
			European Union Conformity		
Manufacturer Symbols (BioFire Diagnostics, LLC)					
	The NOTE symbols explain how to perform the SPOTFIRE R/ST Panel test more efficiently.				
	A BIOFIRE SPOTFIRE Respiratory/Sore Throat (R/ST) Panel			European Union Product Importer	

APPENDIX B

Contact and Legal Information

Customer and Technical Support	
Reach Us on the Web http://www.BioFireDX.com	Reach Us by Phone 1-800-682-2666
Reach Us by E-mail BioFireSupport@bioMerieux.com	Reach Us by Mail 515 Colorow Drive Salt Lake City, UT 84108 USA
Reach Us by Fax (801) 588-0507	



BioFire Diagnostics, LLC
 515 Colorow Drive
 Salt Lake City, UT 84108
 USA



Qarad EC-REP BV
 Pas 257
 B-2440 Geel, Belgium



bioMérieux SA
 376, Chemin de l'Orme
 69280 Marcy l'Etoile-
 France

NOTE FOR CUSTOMERS WITHIN THE EUROPEAN UNION (EU): Any serious incident that has occurred in relation to the device must be reported to BioFire Diagnostics, LLC or local bioMérieux sales representative and the competent authority of the Member State in which the user and/or the patient is established.

© Copyright 2007–2024, BioFire Diagnostics, LLC. All rights reserved.

The information contained in this document is subject to change without notice. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of BioFire Diagnostics, LLC.

BioFire FilmArray Software, Detector, and Metacall software Modules © 2002–2024 BioFire Diagnostics, LLC.

BioFire Diagnostics, BioFire, the BioFire logo, FilmArray and LCGreen are trademarks of BioFire Diagnostics, LLC or BioFire Defense, LLC and are registered trademarks in the United States.

All other names of products and brands appearing in this manual are trademarks or registered trademarks of their respective owners.

The purchase of this product includes a limited, nontransferable license under specific claims of one or more U.S. patents as listed on BioFire Diagnostics' Web site (<http://www.biofiredx.com/legal-notice/>) and owned by BioFire and the University of Utah Research Foundation.

Warranty Information

Product warranty information is available online at:

<http://www.biofiredx.com/support/documents/>

For warranty information for customers outside the United States, contact the local bioMérieux sales representative or an authorized distributor.

APPENDIX C

References

1. Jones, M. S. *et al.* New adenovirus species found in a patient presenting with gastroenteritis. *J. Virol.* **81**, 5978–5984 (2007).
2. Lenaerts, L., De Clercq, E. & Naesens, L. Clinical features and treatment of adenovirus infections. *Rev. Med. Virol.* **18**, 357–374 (2008).
3. Centers for Disease Control and Prevention. Centers for Disease Control and Prevention, National Center for Immunization and Respiratory Diseases (NCIRD), Division of Viral Diseases (DVD) Web site. <http://www.cdc.gov/ncidod/dvrd/revb/respiratory/eadfeat.htm>.
4. Calder, J. A. M. *et al.* Adenovirus type 7 genomic-type variant, New York City, 1999. *Emerg. Infect. Dis.* **10**, 149–152 (2004).
5. Metzgar, D. *et al.* Abrupt emergence of diverse species B adenoviruses at US military recruit training centers. *J. Infect. Dis.* **196**, 1465–1473 (2007).
6. Russell, K. L. *et al.* Transmission dynamics and prospective environmental sampling of adenovirus in a military recruit setting. *J. Infect. Dis.* **194**, 877–885 (2006).
7. Corman, V. M., Muth, D., Niemeyer, D. & Drosten, C. Hosts and Sources of Endemic Human Coronaviruses. in *Advances in Virus Research* vol. Advances in Virus Research 163–188 (Elsevier, 2018).
8. Yang, Y. *et al.* The deadly coronaviruses: The 2003 SARS pandemic and the 2020 novel coronavirus epidemic in China. *J. Autoimmun.* **109**, 102434 (2020).
9. Zaki, A. M., Van Boheemen, S., Bestebroer, T. M., Osterhaus, A. D. & Fouchier, R. A. Isolation of a novel coronavirus from a man with pneumonia in Saudi Arabia. *N. Engl. J. Med.* **367**, 1814–1820 (2012).
10. Chan, J. F. *et al.* Middle East respiratory syndrome coronavirus: another zoonotic betacoronavirus causing SARS-like disease. *Clin. Microbiol. Rev.* **28**, 465–522 (2015).
11. Wu, Y. *et al.* SARS-CoV-2 is an appropriate name for the new coronavirus. *The Lancet* **395**, 949–950 (2020).
12. Cucinotta, D. & Vanelli, M. WHO Declares COVID-19 a Pandemic. *Acta Bio Medica Atenei Parm.* **91**, 157–160 (2020).
13. van der Hoek, L. *et al.* Croup is associated with the novel coronavirus NL63. *PLoS Med.* **2**, e240 (2005).
14. Kahn, J. S. & McIntosh, K. History and recent advances in coronavirus discovery. *Pediatr. Infect. Dis. J.* **24**, S223–227, discussion S226 (2005).
15. Kuypers, J. *et al.* Clinical disease in children associated with newly described coronavirus subtypes. *Pediatrics* **119**, e70–76 (2007).
16. Yuen, K.-S., Ye, Z.-W., Fung, S.-Y., Chan, C.-P. & Jin, D.-Y. SARS-CoV-2 and COVID-19: The most important research questions. *Cell Biosci.* **10**, 40 (2020).
17. Petrosillo, N., Viceconte, G., Ergonul, O., Ippolito, G. & Petersen, E. COVID-19, SARS and MERS: are they closely related? *Clin. Microbiol. Infect.* S1198743X20301713 (2020) doi:10.1016/j.cmi.2020.03.026.
18. Kahn, J. S. Epidemiology of human metapneumovirus. *Clin. Microbiol. Rev.* **19**, 546–557 (2006).
19. van den Hoogen, B. G. *et al.* A newly discovered human pneumovirus isolated from young children with respiratory tract disease. *Nat Med* **7**, 719–724 (2001).
20. Falsey, A. R., Erdman, D., Anderson, L. J. & Walsh, E. E. Human metapneumovirus infections in young and elderly adults. *J. Infect. Dis.* **187**, 785–790 (2003).
21. VAN DEN HOOGEN, B. G., OSTERHAUS, D. M. E. & FOUCHIER, R. A. M. Clinical impact and diagnosis of human metapneumovirus infection. *Pediatr. Infect. Dis. J.* **23**, S25–S32 (2004).
22. Esper, F. *et al.* A 1-year experience with human metapneumovirus in children aged < 5 years. *J. Infect. Dis.* **189**, 1388–1396 (2004).
23. Budd, A. *et al.* Chapter 6: Influenza. in *Manual for the Surveillance of Vaccine-Preventable Diseases* (2017).
24. Bammer, L., Fukuda, A., Klimov, and N. Cox. Influenza. in *VPD Surveillance Manual* (2002).
25. Update: influenza activity - United States, August 30, 2009–March 27, 2010, and composition of the 2010–11 influenza vaccine. *MMWR Morb. Mortal. Wkly. Rep.* **59**, 423–430 (2010).
26. Bammer, L., Fukuda, A., Klimov, and N. Cox. Influenza. in *VPD Surveillance Manual* (2002).
27. Morens, D. M., Taubenberger, J. K. & Fauci, A. S. Predominant role of bacterial pneumonia as a cause of death in pandemic influenza: implications for pandemic influenza preparedness. *J Infect Dis* **198**, 962–70 (2008).
28. Henrickson, K. J. Parainfluenza viruses. *Clin. Microbiol. Rev.* **16**, 242–264 (2003).
29. Senchi, K., Matsunaga, S., Hasegawa, H., Kimura, H. & Ryo, A. Development of oligomannose-coated liposome-based nasal vaccine against human parainfluenza virus type 3. *Front Microbiol* **4**, 346 (2013).
30. Lau, S. K. P. *et al.* Human parainfluenza virus 4 outbreak and the role of diagnostic tests. *J. Clin. Microbiol.* **43**, 4515–4521 (2005).
31. Fry, A. M. *et al.* Seasonal trends of human parainfluenza viral infections: United States, 1990–2004. *Clin. Infect. Dis.* **43**, 1016–1022 (2006).
32. Mohapatra, S. S. & Boyapalle, S. Epidemiologic, experimental, and clinical links between respiratory syncytial virus infection and asthma. *Clin. Microbiol. Rev.* **21**, 495–504 (2008).
33. Anderson, L. J., Hendry, R. M., Pierik, L. T., Tsou, C. & McIntosh, K. Multicenter study of strains of respiratory syncytial virus. *J. Infect. Dis.* **163**, 687–692 (1991).
34. Falsey, A. R. & Walsh, E. E. Respiratory syncytial virus infection in adults. *Clin. Microbiol. Rev.* **13**, 371–384 (2000).
35. Hall, C. B. Respiratory syncytial virus and parainfluenza virus. *N. Engl. J. Med.* **344**, 1917–1928 (2001).
36. Anzueto, A. & Niederman, M. S. Diagnosis and treatment of rhinovirus respiratory infections. *Chest* **123**, 1664–1672 (2003).

37. Jacques, J. *et al.* Epidemiological, molecular, and clinical features of enterovirus respiratory infections in French children between 1999 and 2005. *J. Clin. Microbiol.* **46**, 206–213 (2008).
38. Jacobs, S. E., Lamson, D. M., George, K. S. & Walsh, T. J. Human rhinoviruses. *Clin. Microbiol. Rev.* **26**, 135–162 (2013).
39. Sawyer, M. H. Enterovirus infections: diagnosis and treatment. *Curr. Opin. Pediatr.* **13**, 65–69 (2001).
40. Centers for Disease Control and Prevention. Centers for Disease Control and Prevention, Pertussis (Whooping Cough) Web Site.
41. World Health Organization. WHO Immunization, Vaccines, and Biologics; Pertussis Web Site.
42. Mattoo, S. & Cherry, J. D. Molecular pathogenesis, epidemiology, and clinical manifestations of respiratory infections due to *Bordetella pertussis* and other *Bordetella* subspecies. *Clin. Microbiol. Rev.* **18**, 326–382 (2005).
43. Srugo, I. *et al.* Pertussis infection in fully vaccinated children in day-care centers, Israel. *Emerg. Infect. Dis.* **6**, 526–529 (2000).
44. Heininger, U. *et al.* Clinical characteristics of illness caused by *Bordetella parapertussis* compared with illness caused by *Bordetella pertussis*. *Pediatr. Infect. Dis. J.* **13**, 306–309 (1994).
45. Cherry, J. D. & Seaton, B. L. Patterns of *Bordetella parapertussis* respiratory illnesses: 2008–2010. *Clin. Infect. Dis. Off. Publ. Infect. Dis. Soc. Am.* **54**, 534–537 (2012).
46. Bouchez, V. & Guiso, N. *Bordetella pertussis*, *B. parapertussis*, vaccines and cycles of whooping cough. *Pathog. Dis.* **73**, ftv055 (2015).
47. Hahn, D. L., Azenabor, A. A., Beatty, W. L. & Byrne, G. I. Chlamydia pneumoniae as a respiratory pathogen. *Front. Biosci. J. Virtual Libr.* **7**, e66–76 (2002).
48. Grayston, J. T. Chlamydia pneumoniae, strain TWAR pneumonia. *Annu. Rev. Med.* **43**, 317–323 (1992).
49. Kuo, C. C., Jackson, L. A., Campbell, L. A. & Grayston, J. T. Chlamydia pneumoniae (TWAR). *Clin. Microbiol. Rev.* **8**, 451–461 (1995).
50. Peeling, R. W. & Brunham, R. C. Chlamydiae as pathogens: new species and new issues. *Emerg. Infect. Dis.* **2**, 307–319 (1996).
51. Outbreak of community-acquired pneumonia caused by *Mycoplasma pneumoniae*--Colorado, 2000. *Can. Commun. Dis. Rep. Relevé Mal. Transm. Au Can.* **27**, 104–107 (2001).
52. Klement, E. *et al.* Identification of risk factors for infection in an outbreak of *Mycoplasma pneumoniae* respiratory tract disease. *Clin. Infect. Dis. Off. Publ. Infect. Dis. Soc. Am.* **43**, 1239–1245 (2006).
53. Centers for Disease Control and Prevention. Centers for Disease Control and Prevention, Disease Listing: *Mycoplasma pneumoniae* Web Site.
54. Group A Strep | Surveillance | GAS | CDC. <https://www.cdc.gov/groupastrep/surveillance.html>.
55. Bisno, A. L. Acute Pharyngitis: Etiology and Diagnosis. *Pediatrics* **97**, 949–954 (1996).
56. Ebell, M. H., Smith, M. A., Barry, H. C., Ives, K. & Carey, M. The rational clinical examination. Does this patient have strep throat? *JAMA* **284**, 2912–2918 (2000).
57. Brandt, C. M. & Spellerberg, B. Human infections due to *Streptococcus dysgalactiae* subspecies *equisimilis*. *Clin. Infect. Dis. Off. Publ. Infect. Dis. Soc. Am.* **49**, 766–772 (2009).
58. Baracco, G. J. Infections Caused by Group C and G *Streptococcus* (*Streptococcus dysgalactiae* subsp. *equisimilis* and Others): Epidemiological and Clinical Aspects. *Microbiol. Spectr.* **7**, (2019).
59. Thai, T. N., Dale, A. P. & Ebell, M. H. Signs and symptoms of Group A versus Non-Group A strep throat: A meta-analysis. *Fam. Pract.* **35**, 231–238 (2018).
60. Tiemstra, J. & Miranda, R. L. F. Role of non-group a streptococci in acute pharyngitis. *J. Am. Board Fam. Med. JABFM* **22**, 663–669 (2009).
61. Gerber, M. A. *et al.* Community-wide outbreak of group G streptococcal pharyngitis. *Pediatrics* **87**, 598–603 (1991).
62. Efstratiou, A. Outbreaks of human infection caused by pyogenic streptococci of Lancefield groups C and G. *J. Med. Microbiol.* **29**, 207–219 (1989).
63. Kakuya, F. *et al.* Acute Pharyngitis Associated With *Streptococcus dysgalactiae* Subspecies *equisimilis* in Children. *Pediatr. Infect. Dis. J.* **37**, 537–542 (2018).
64. Biosafety in Microbiological and Biomedical Laboratories (BMBL).

Revision History

Version	Revision Date	Description of Revision(s)
01	August 2023	Initial release
02	May 2024	<p>Additions:</p> <ul style="list-style-type: none"> • 'Procedure' section Step 3: Prepare Patient or QC Sample – new statement to visually inspect the volume of the sample in the fixed volume transfer pipette. • 'Procedure' section Step 5: Start Run - caution statement for selecting the correct sample type. • 'General Precautions' section includes two new precautions regarding selecting the correct sample type (#2) and desiccant statement (#5). • 'Materials Provided' section– ingredients statement added. • Carry Over section <p>Updates:</p> <ul style="list-style-type: none"> • Clinical Performance section – reorganization of data to integrate with new testing data from September 2022 to May 2023 (TS specimens only). • SARS-CoV-2 <i>in silico</i> reactivity data updated based on analysis of available sequences through December 31, 2023. • Table 37 updated to include new cross-reactivities • Other updates for clarity or accuracy. • Customer Support phone number and email • 'Materials Required But Not Provided' section updated with new table

INDEX

Intended Purpose	1
Intended Use	1
Intended User and Use Environment	2
Summary and Explanation of the Test	3
Summary of Detected Organisms	3
Principle of the Procedure	6
Materials Provided	6
Materials Required But Not Provided	7
Warnings and Precautions	8
General Precautions	8
Safety Precautions	8
Laboratory Precautions	10
Precautions Related to Public Health Reporting	11
Precaution Related to REACH Regulation (EC 1907/2006)	11
Reagent Storage, Handling, and Stability	12
Sample Requirements	12
Quality Control	13
Internal Process Controls	13
External Controls	13
Procedure	13
Step 1: Prepare Pouch	13
Step 2: Hydrate Pouch	14
Step 3: Prepare Patient or QC Sample	15
Step 4: Load Patient or QC Sample	16
Step 5: Start Run	16
Interpretation of Results	17
Assay Interpretation	17
Organism Interpretation	18
SPOTFIRE R/ST Panel Test Report	20
Limitations	22
Expected Values	25
Performance Characteristics	31

Clinical Performance	31
Analytical Performance Characteristics	36
Limit of Detection	36
Analytical Reactivity (Inclusivity)	38
Analytical Specificity (Cross-Reactivity & Exclusivity).....	43
Near-LoD/Reproducibility Evaluation	47
Interference	51
Carry-Over.....	53
Appendix A	54
Symbols Glossary	54
Appendix B.....	55
Contact and Legal Information	55
Warranty Information.....	55
Appendix C.....	56
References	56
Index.....	59