3P® SMART Culture Media

Marilyn Velez – bioMérieux, Inc. Lombard, IL Site

PERFORMANCES STUDY – CT 3P IR, TSA NEUT 3P IR, TSA 3P IR

DESCRIPTION

Traditional culture media is integral for any environmental monitoring program.^{1,2} Whether it is used for microbial contamination monitoring or making product release decisions, choosing the right prepared culture media to deliver accurate results is essential.³

bioMérieux's irradiated 3P[®] SMART (Pharmaceutical Proven Performance) culture media is a premium product offering specifically designed for use in pharmaceutical cleanroom environments and have been validated for optimal performances particularly

on environmental strains.

The following presents the performances of the 3P SMART culture media range during development, demonstrating to be the most comprehensive solution for an environmental monitoring (EM) program.



MATERIALS AND METHODS

bioMérieux Product Development

Product development at bioMérieux complies with ISO 13485 guidelines⁴ and strict internal standards.

For each 3P SMART reference (Table 1), a minimum of 3 separate batches are manufactured for the purpose of Verification and Stability testing (shelf-life validation).

Verification testing is executed to ensure that the product can perform as intended.

Stability testing monitored in real-time, evaluates the 3 batches with applied thermal shocks to confirm that the product can consistently perform through expiration. Thermal shocks are applied to mimic storage conditions that could occur during shipping. Batches are then assessed at regular intervals for microbiological and physicochemical performance (i.e., growth promotion, neutralization capability when applicable, pH, and cosmetic aspect).

Table 1. 3P SMART References				
Ref. # Item Description		Plate Dish Format		
420765	Irradiated Count-Tact 3P agar LOCKSURE (CT 3P IR)	Contact plate		
423723	Irradiated Trypcase Soy 3P Agar with Neutralizers (TSA NEUT 3P IR)	90 mm settle plate		
423725	Irradiated Trypcase Soy 3P Agar (TSA 3P IR)	90 mm settle plate		

Microbiological Performance

Table 3 Strains tested during Stability

During verification, growth promotion was evaluated with a total of 110 strains for CT 3P IR & TSA 3P IR and 94 microorganisms for TSA NEUT 3P IR that include pharmacopeia-required strains and cleanroom wild isolates sourced from an internal library collection.

Table 2. Number of strains tested per each microorganism category					
Microorganism Type	CT 3P IR 110 Strains total	TSA NEUT 3P IR 94 Strains total	TSA 3P IR 110 Strains total		
Molds	15	14	15		
Yeasts	9	7	9		
Gram Positive Cocci	22	21	21		
Gram Positive Rods	19	18	19		
Gram Negative Rods	38	33	40		
Anaerobes	7	7	7		

For stability, growth promotion was evaluated with the strains listed in Table 3 in BioBall[®] Multishot-550 format

Ref. #	Microorganism	BioBall Strain	Equivalent ATCC Strain	
56011	Aspergillus brasiliensis	NCPF 2275	16404	
56013	Candida albicans	NCPF 3179	10231	
56012	Bacillus subtilis	NCTC 10400	6633	
56016	Escherichia coli	NCTC 12923	8739	
56017	Pseudomonas aeruginosa	NCTC 12924	9027	
56019	Staphylococcus aureus	NCTC 10788	6538	
	1.5			

List of Microorganisms:

Molds	Gram Positive Cocci	Gram Positive Rods	Gram Negative Rods	K. oxytoca 104844	Yeast
Alternaria spp. 0411781	A. viridans 12950	Aureobacterium spp 0502715	A. baumanii ATCC 19606	K. pneumoniae ATCC 13883	C. albicans ATCC 10231
A. brasiliensis ATCC 16404	E. durans ATCC19432	B. cereus ATCC 13824	A. hydrophilia 101139	L. adecarboxylata 113532	C. glabrata 300502
A. corymbifera ATCC 46773	E. faecalis ATCC 29212	B. cereus 200025	A. sobria 111035	L. adecarboxylata 113764	C. sake 301624
A. fumigatus 0411768	E. faecium 12834	B. coagulans ATCC 9172	A. cryaerophilus 140160	M. mesophylicum 0306751	C. neoformans ATCC 32045
A. fumigatus 1401601	E. gallinarum 12838	B. firmus ATCC 29789	B. vesicularis 101556	P. vulgaris 104972	G. candidum 300575
A. terreus 0411763	K. kristinae 18527	B. licheniformis 200019	B. cepacia ATCC 25416	P. stuartii ATCC 43258	R. mucilaginosa 305157
C. keratinophilum 0411774	K. rhizophila ATCC 9341	B. megaterium 200061	B. gladioli 110263	P. aeruginosa ATCC 9027	S. cerevisiae ATCC 9763
E. floccosum 0411793	K. rosea 8425	B. pumillus 2000111	C.indologenes ATCC 29897	P. aeruginosa ATCC 14207	T. mucoïdes 304040
P. lilacinus ATCC 90461	Micrococcus spp. 10085	B. simplex 202066	C. luteola 0011707	P. fluorescens 106364	Z. rouxii 305593
P. aurantiogriseum ATCC 16025	S. aureus ATCC 6538	B. subtilis ATCC 6633	Citrobacter braakii 100772	P. putida 101109	
P. commune 0909600	S. aureus ATCC 25993	B. nesterukovii 1401619	Citrobacter freundii ATCC 43864	R. aquatilis 109226	Anaerobes
P. iltalicum ATCC 48114	S. capitis 2550	B. brevis 0010703	C. terrigena ATCC 11330	R. picketti 109228	B. fragilis ATCC 25285
R. pusillus 0411761	S. epidermidis 3052	B. epidermidis ATCC 35514	D. acidovorans 106523	R. planticola 104937	B. vulgatus ATCC 8482
S. brevicaulis 0411778	S. epidermidis ATCC 12228	C. haemolyticum 9602770	E. aerogenes ATCC 13048	S. enterica ATCC 19126	C. perfringens ATCC 13124
T. mentagrophytes ATCC 9533	S. haemolyticus 3047	C. striatum ATCC BAA-1293	E. cloacae 105620	S. gallinarum 106305	C. sporogenes ATCC 11437
	S. haemolyticus ATCC 29970	L. mesenteroides ATCC 23386	E. coli ATCC 8739	S. fonticola 106356	P. acnes ATCC 6919
	S. hominis 8045	L. fusiformis 202043	E. coli ATCC 25922	S. liquefaciens ATCC 27592	P. granulosum ATCC 25564
	S. pyogenes ATCC 19615	P. lautus 202122	F. oryzihabitans 0012704	S. flexeneri ATCC 33948	
	S. saprophyticus 3039	V. pantothenticus 202138	H. alvei ATCC 13337	S. paucimobilis 108706	
	S. warnerii 12953			S. maltophilia 109782	

Calibrated inoculum between 10 and 100 CFU from fresh cell suspension or BioBall format was prepared for each microorganism. Plates were inoculated in duplicate and incubated at 20-25°C or 30-35°C according to the growth requirements of the microorganism for an appropriate culture duration.

The recovery rates were calculated from a non-irradiated TSA control with an acceptance criteria between 50-200%.

Neutralization Capability

CT 3P IR and TSA NEUT 3P IR contain neutralizing components to inactivate any residual disinfectants present in the air, operator gloves, or surface to be tested.

These 3P products were challenged to evaluate their ability to neutralize commercially available disinfectants from different classification categories - four (4) for CT 3P IR and three (3) for TSA NEUT 3P IR. The recovery was compared to counts on non-irradiated TSA control (without disinfectant) with an acceptance criteria between 50-200%.

Table 4. Commercially available disinfectants tested per disinfectant type					
Alcohol	Oxidizing Agents	Phenolic Compounds	Quaternary Ammonium Compounds		
70% IPA	Clorox [®] Spor-Klenz [®] *	Vesphene®* Sporicidin®	Decon Quat [®] 100 Texwipe [®] Tex Q™ TX 650 Coverage [®] Spray TB*		
			concluge opila) in		

*Disinfectants tested only on TSA NEUT 3P IR Disinfectants in bold tested during Stability

Neutralization was assessed against two groups of microorganisms: USP <61> required⁵ and common wild isolates during verification (Table 5), and strains tested during stability are shown in Table 6.

Table 5. Microorganisms tested during Verification		Table 6. Microorganisms tested during Stability		
PharmacopeiaSupplemental OrganismsRequired(including wild isolates)		Strain BioBall Multishot Format	Equivalent ATCC Strain	
C. albicans ATCC 10231	S. paucimobilis 108706	Bacillus subtilis NCTC 10400	ATCC 6633	
A. brasiliensis ATCC 16404	C. striatum ATCC BAA-1293	Staphylococcus aureus NCTC 10788	ATCC 6538	
P. aeruginosa ATCC 9027	K. rhizophila ATCC 9341	Candida albicans NCPF 3179	ATCC 10231	
S. aureus ATCC 6538	B. pumilus 200111			
B. subtilis ATCC 6633	E. coli ATCC 8739			

Each disinfectant was prepared based on a pre-determined concentration that inhibits growth on standard TSA (without neutralizers) but are effectively neutralized by the compounds in the 3P media according to the microorganisms' sensitivity to each product. Disinfectant solutions of 200 µL for contact plates and 500 µL for settle plates were applied directly to the surface and allowed to dry prior to inoculating the designated microorganism.

Incubation conditions: up to 48 hours at 30-35°C for bacteria and up to 7 days at 20-25°C for yeast and molds.

Acceptance criteria for microorganism recovery was expected to be between 50-200% compared to counts on non-irradiated TSA control.

Cosmetic Aspect

Condensation and Exudation

3P culture media was evaluated for condensation/exudation inside the plate and sleeve packaging over the shelf-life of the product at both 2-8°C and room temperature (RT) storage using the qualitative and quantitative scale outlined in Table 7.

Table 7. Visible conder	sation and exudation scoring sca	ale	
Condensation Score	Description	Dro Diam	
0 – Dry	No visible fog or water droplets		
1 – Fog and / or Tiny Droplets	Presence of fog or foggy patches with no separate drops of water	<1	
2 – Small Drops	Distinguishable small droplets	1 -	
3 – Full Drops	One or more drops larger than 4mm	>4	
4 – Water	Droplets pooling to form large areas of visible moisture	large co dro (Pooling	

Dehydration (Resistance to Laminar Flow)

Settle plate formats have a 30 mL fill volume optimized for air sampling applications (active & passive) to reinforce resistance to dehydration effects. To mimic passive air sampling, TSA 3P IR and TSA NEUT 3P IR were left to dry for 4.5 hrs under vertical laminar flow and then incubated 4 days at 20-25°C followed by 3 days a 30-35°C.

Dehydration was assessed over the shelf-life of the product using the scale shown in Figure 1.



Figure 1. Dehydration evaluation scale



RESULTS AND DISCUSSION - VERIFICATION

Microbiological Performance



Figure 2. Average percent recovery for all strains tested per microorganism type for each 3P SMART reference

Excellent growth promotion performance for all 3P SMART references confirmed by testing a wide array of microbial strains to enhance detection of any contamination during environmental monitoring.

Neutralization Capability





Neutralization Capability Comparison

Figure 4. Average recovery rate for all strains tested against 3 disinfectant types on CT 3P IR & TSA NEUT 3P IR

CT 3P IR and TSA 3P with Neutralizers demonstrate effective neutralization for all disinfectants tested. Both products are able to neutralize various disinfectants and recover pharmacopoeia and supplemental strains within 50-200% percent recovery.

Isolator Compatibility

Specific 3P media packaging has been designed to be resistant to Vaporized Hydrogen Peroxide (VHP) cycle typically performed to decontaminate all materials being passed into the isolator.

Vaporized Hydrogen Peroxide (VHP)



w/o VPH Decontamination Figure 5. Percent recovery rate observed before (w/o) and after VHP decontamination cycle for all 3P SMART references

Similar growth performance of the media without and with VHP decontamination demonstrate no organisms were inhibited by VHP exposure. All 3P SMART references are validated with VHP cycle for isolator compatibility.

STABILITY – (Shelf-life Validation)



Figure 6. AVG percent recovery for all 3P SMART references through expiration

Neutralization Capability

Microbiological

Performance

Consistent growth

recovery for all 3P

SMART references

through the end of

shelf-life.

Only two disinfectant types were tested through the end of stability for both CT 3P IR and TSA NEUT 3P IR: phenolic and guaternary ammonium compounds.



quaternary ammonium disinfectants on TSA NEUT 3P IR

High neutralization performances for CT 3P IR and TSA NEUT 3P IR to obtain good growth recovery against phenolic and quaternary ammonium disinfectants through the end of shelf-life.

Cosmetic Aspect

Condensation and Exudation





temperature.

Dehydration – Resistance to laminar flow exposure

TSA NEUT 3P IR

TSA 3P with and without neutralizers in settle plate format demonstrate robust media behavior with no major signs of dehydration during passive air sampling throughout its shelf-life duration.

In addition to successfully support growth on an array of microorganisms, neutralization efficacy across different disinfectants, VHP compatible during isolator decontamination, bioMérieux guarantees:

	Contact Plate		90 mm Settle Plate			
	3P	3P SMART	3P	3P SMART	3P	3P SMART
	CT3P Ref. 410250/418049	CT 3P IR Ref. 420765	TSA3PN Ref. 419013/419014	TSA NEUT 3P IR Ref. 423723	TSA3P Ref. 410251/418140	TSA 3P IR Ref. 423725
Appearance	✓	~	✓	~	~	~
рН						
Growth Promotion						
Neutralization Activity					NA	NA
pH Growth Promotion Neutralization Activity	~	~	\checkmark	\checkmark	✓ NA	✓ NA

REFERENCES:

- Manufacturing Practice.

Phenolic

Quaternary

Ammonium



Figure 9. Condensation overview for 3P SMART references inside plate lid and innermost sleeve when stored at 2-8°C & RT

Through the end of the stability study, all 3P SMART plates demonstrate minimal condensation formation inside plate lids and innermost sleeve at either storage



Figure 10. Cosmetic defaults overview for TSA NEUT 3P IR post dehydration process

Figure 11. Cosmetic defaults overview for TSA 3P IR post dehydration process

CONCLUSION

Performances of 3P SMART culture media during development demonstrate excellent biological and physicochemical behavior with consistent proven product performance up to the end of shelf-life.

 Completely controlled manufacturing process ensuring consistent reliability between all batches and monolot capacity.

 Maximum convenience for storage and inventory management thanks to a flexible storage between 2-25°C.

Innovative packaging system that enables a long shelf-life.

• Security of your crucial EM samples during critical transport and incubation steps with bioMérieux's unique LOCKSURE[®] closure system design.

• Support data integrity and traceability solutions with GS1 data matrix barcode.

Equivalent performance to the original 3P[®] with video-ink jet marking:

3P SMART culture media represents the state of the art in culture media technology from both a formulation and packaging perspective allowing bioMérieux to offer the most comprehensive, robust, and reliable solution for environmental monitoring.

1. USP <1116> Microbiological Control and Monitoring of Aseptic Processing.

2. FDA. 2004. Guidance for Industry: Sterile Drug Products produced by Aseptic Processing – Current Good

3. Sandle, Tim. 2014. Assessment of Culture Media in Pharmaceutical Microbiology. American Pharmaceutical Review. Retrieved from https://www.americanpharmaceuticalreview.com/Featured-Articles/163589 Assessment-of-Culture-Media-in-Pharmaceutical-Microbiology/.

4. ISO 13485 Medical Devices – Quality Management Systems - Requirements for Regulatory Purposes. 5. USP <61> Microbial Examination of Nonsterile Drug Products: Microbial Enumeration Tests.