

BIOMÉRIEUX

BACT/ALERT® VIRTUO®

FAST FULLY AUTOMATED BLOOD CULTURE SOLUTION



PIONEERING DIAGNOSTICS

BACT/ALERT® VIRTUO®

INTELLIGENT INNOVATION IN BLOOD CULTURE.
BUILT ON 10 YEARS OF PROVEN PERFORMANCE.

Maximize efficiency while providing clinically impactful results.

Experience matters. Over 60 years of serving laboratories provide the foundation for our unshakable commitment to diagnostic innovation, quality, and reliability.

BACT/ALERT® VIRTUO® brings a decade of real-world performance, trusted by labs worldwide to deliver fast, accurate and fully automated blood culture results. With proven expertise in system integration and customer support, **BACT/ALERT® VIRTUO® adapts to your lab's unique needs, offering peace of mind and operational excellence.**

BACT/ALERT® VIRTUO® can help gain valuable time for high-risk patients and support antimicrobial stewardship, with fast detection of microorganisms responsible for bloodstream infections.

FROM 4 TO 9 HOURS FASTER TIME TO DETECTION IN PATIENTS ON ANTIMICROBIALS^{1,2}

- ✓ Enhanced speed through thermal stability and ultra-sensitive algorithmic detection
- ✓ Significantly shorter time to detection in 72.7% of the tested organisms³
- ✓ Save up to 10 hours in sample reception area, compared to conventional blood culture workflow during the night shift⁴
- ✓ Move on more quickly to identification and antimicrobial susceptibility testing⁵

LEVERAGE RESOURCES WHERE MOST VALUABLE

- ✓ Unique Drop & Go technology with fully automated loading—easy for anyone, anytime, 24/7
- ✓ Reduced steps and hands-on time by 63%¹
- ✓ Up to 6x faster bottle loading and 14x faster negative bottle unloading⁶
- ✓ Automated negative bottle unloading allows staff to focus on added-value samples

EFFICIENT & ADAPTABLE WORKFLOW

- ✓ Unsurpassed 24-hour delayed entry: flexible bottle transport and loading^{7,10}
- ✓ Scalable to meet your needs: from 428 to 1712 bottles
- ✓ Space-efficient: best capacity/footprint ratio¹¹
- ✓ Modular design: fits core labs and decentralized workflows



IN CONTROL, IN REAL TIME

- ✓ Remote alarms and alert notifications on LED signal tower for flexible lab organization
- ✓ Real-time notifications:
 - Blood fill level measurement per bottle
 - Non-conforming blood culture samples after loading, based on personalized thresholds
- ✓ Built for simplicity: real-time access to blood culture KPIs with MAESTRIA™
- ✓ Optimized quality: cell quality control check, calibration, and bottle scanning for traceability

TRUSTED CHOICE FOR ENHANCED PATIENT CARE

- ✓ BACT/ALERT® FAN® PLUS bottles offer superior antimicrobial binding kinetics and better antibiotic neutralization¹²⁻¹⁵
- ✓ Increased recovery by 44%, and up to 71% more with antibiotics^{14,16}
- ✓ Safe handling in transit: colorimetric technology offers stable results and up to 50% fewer false negatives due to delayed entry than fluorometry⁷
- ✓ Designed and validated for sterile body fluids and blood bank sterility control^{13,17}

Empower Your Lab with Comprehensive Blood Culture Solutions.

Our integrated approach – from blood sample collection, to blood culture detection and data management – positions your lab to optimize workflow and improve patient management.



SECURE HOLDERS

Pre-Analytical



BACT/ALERT® VIRTUO®

Analytical



MAESTRIA™ MICROBIOLOGY
MIDDLEWARE




SUBCULTURE DEVICES

Post-Analytical

Transform Data into Actionable Insights.

MA

MAESTRIA™ is the next generation microbiology middleware to help you increase visibility over your laboratory operations, integrate continuous improvement processes, and reduce time spent on non-value-added tasks.



KEY PERFORMANCE DATA

Monitor lab performance KPIs with MAESTRIA™ and target areas for improvement in your blood culture process.

- Contamination associated with blood collection
- Blood volume collected
- Positive vs. negative analysis
- Time-to-detection for blood culture isolates
- Blood culture utilization: number of bottles per set, per patient and/or ward



EASY ACCESS

- Blood culture data analytics and results consolidation
- Remote access, including bottle status remote notification

EVEN MORE
WITH MAESTRIA™

- ✓ Monitor lab activities and KPIs
- ✓ Infection control and AMS statistical reporting
- ✓ Lab efficiency reporting

Leverage Our Support & Expertise to Streamline Your Lab Operations.

Our **lab consultancy** experts can help your lab optimize workflows and enhance productivity, while our **blood culture medical education** experts support better quality management and optimized practices. In addition, our **highly responsive service team** is here to help keep your lab operations running continuously, via our hotline or on-site intervention.

NEED HELP? JUST ASK. CONNECT WITH YOUR BIOMÉRIEUX REPRESENTATIVE TODAY.

1. Halperin A.V. et al. (2022). Impact of automated blood culture systems: crossover clinical trial. *Microbiol Spectr*.
2. Mazulli et al. (2023). Impact of Implementation of BacT/Alert Virtuo on Blood Culture Time to Positivity in Sepsis Patients. *Microbiol Spectr*.
3. Somily AM. et al. (2018). Time-to-detection of bacteria and yeast with the BACTEC FX versus BacT/Alert Virtuo blood culture systems. *Ann Saudi Med*.
4. Pean de Ponfily G. et al. (2021). Impact of 24/7 bottle loading in an automated incubator. *Eur J Clin Microbiol Infect Dis*.
5. Dai Y. et al. (2024). Clinical and economic evaluation of blood culture whole process optimization in critically ill adult patients with positive blood culture. *Int J Antimicrob Agents*.
6. Internal R&D study data: BACT/ALERT VIRTUO vs. BACT/ALERT 3D
7. Akan O.A., Yildiz E. (2006). Effect of delayed entry in two blood culture systems. *Diagn Microbiol Infect Dis*.
8. Wilms M.C. et al. (2009). Effect of preincubation temperature on fastidious organism detection. *J Microbiol Methods*.
9. Adamik et al. (2021). Effect of delayed entry on performance of the BACT/ALERT FAN PLUS bottles in the BACT/ALERT VIRTUO blood culture system. *Eur J Clin Microbiol Infect Dis*.
10. Deslandes et al. (2022). Effect of delayed entry of blood culture bottles in BACTEC automated blood culture system in the context of laboratory consolidation. *Sci Rep*.
11. Internal data. BACT/ALERT VIRTUO vs BD BACTEC FX systems comparison.
12. Hang Y. et al. (2024). Comparison of antimicrobial neutralization across three culture media. *Microbiol Spectr*.
13. BACT/ALERT FAN Plus Package Insert (Ref. 410851/410852).
14. Lovern D. et al. (2016). Antimicrobial binding and growth kinetics in FA Plus vs. BACTEC Aerobic/F Plus media. *Eur J Clin Microbiol Infect Dis*.
15. Katzin B. (2019). A Controlled Study Comparing Anaerobic Blood Culture Bottles Tested Using Automated Blood Culture Systems. *29th ECCMID Amsterdam*.
16. Qin et al. (2024). Comparative evaluation of BacT/ALERT VIRTUO and BACTEC FX400 blood culture systems for the detection of bloodstream infections. *Microbiol Spectr*.
17. She R.C. et al. (2018). Performance of the BacT/ALERT Virtuo for sterile body fluid cultures: multicenter study. *Clin Microbiol Infect*.