

BIOFIRE® SPOTFIRE® REDEFINING CONNECTIVITY



SEAMLESS CONNECTIVITY

To streamline your workflow, the BIOFIRE SPOTFIRE System offers bi-directional communication in HL7 or POCT01-A2 to provide end-to-end data flow between your connected BIOFIRE SPOTFIRE device(s) and your LIS, your Point of Care (POC) middleware and/or any data manager to post to the patient's record.

Take a step further and unleash the potential of your data. BIOFIRE® FIREWORKS™ will transform your operations into valuable insights, while VILINK® will be your partner in maximizing the efficiency of your system.

DISCOVER THE POWER OF CONNECTIVITY

COMMUNICATION PROTOCOL FLEXIBILITY



POC DATA MANAGEMENT SYSTEM

Connect your BIOFIRE SPOTFIRE device(s) using the industry standard POCT01-A2 protocol to facilitate communication with your POC middleware or any data management system.

- Offer consolidated support and supervision for operators at locations across the organization, and uniformize training certification.
- Automatically transmit Quality Control (QC) results to your POC middleware, or any data management system, to help monitor & maintain test reliability across your BIOFIRE SPOTFIRE devices network.

LIS

Connect your BIOFIRE SPOTFIRE device(s) using the HL7* protocol to facilitate communication with your LIS.

*XML format for uni-directional communication is also an option



GET RESULTS SOONER

BIOFIRE SPOTFIRE enables electronic patient test results reporting, reducing manual entry & decreasing result turnaround time. All you need is access to a local network via an ethernet or Wi-Fi connection.



REDUCE REPORTING ERRORS

BIOFIRE SPOTFIRE accurately translates patient test results into a machine-readable format, eliminating the need for manual data entry and reducing the potential for misreported patient results.

TRANSFORM DATA INTO ANALYTICS

BIOFIRE® FIREWORKS™

Connect your BIOFIRE SPOTFIRE device(s) with BIOFIRE FIREWORKS, our centralized data management portal. Gain total insight into system performance, utilization, pathogen surveillance, workflow management and more.

OPTIMIZE INSTRUMENT UPTIME

VILINK®

Connect your BIOFIRE SPOTFIRE device(s) with VILINK®, our intelligent instrument management and support tool for real-time system monitoring, updates management as well as remote access for troubleshooting, repairs and training purposes.

CUSTOMER SUPPORT

Our goal is to provide you with the best support available while minimizing any disruption to your workflow. For BIOFIRE products technical assistance and support, please contact your local bioMérieux representative or distributor.

Product availability varies by country. Please contact your local bioMérieux representative for details.

bioMérieux S.A.
69280 Marcy l'Etoile • France
Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com

Manufactured by:
BioFire Diagnostics, LLC
515 Colorow Drive, Salt Lake City
UT 84108, USA