



E-TEST[®] MEV


(0.004/8 - 64/8 µg/mL)

FDA 510(k) CLEARED

E-TEST Meropenem/Vaborbactam

(0.004/8 - 64/8 µg/mL)

REF. 421560



Single pack | 30 Strips


For determining meropenem/vaborbactam MIC for *E.coli*, *K. pneumoniae*, and *Enterobacter cloacae* complex.

Available from bioMérieux, Inc.

Performance characteristics are provided below.


Organism	Strains	% Essential Agreement	% Category Agreement
	(N)	(EA)	(CA)
<i>Enterobacteriaceae</i>	594	95.8%	99.3%

Notes: EA = % of MIC values within ± 1 dilution of the reference method. Refer to Package Insert for additional performance data.



To place an order, visit

www.biomerieuxDIRECT.com



E-TEST® TLA

(0.002-32 µg/mL)

FDA 510(k) CLEARED

E-TEST Telavancin (0.002-32 µg/mL)

REF. 418731

Single pack | 30 Strips

For determining telavancin MIC for *Staphylococcus aureus* (including methicillin resistant isolates).

Available from bioMérieux, Inc.

Performance characteristics are provided below.

Organism	Strains	% Essential Agreement	% Category Agreement
	(N)	(EA)	(CA)
<i>S. aureus</i> (including methicillin resistant isolates)	376	98.4%	97.9%

Note: EA = % of MIC values within ± 1 dilution of the reference method. Refer to Package Insert for additional performance data.



For more information, contact your local bioMérieux, Inc. representative.

Customer Service: (800) 682-2666 | www.biomerieux-usa.com/ETESTTLA

To place an order, visit

www.biomerieuxDIRECT.com