



# VIDAS<sup>®</sup> 3

Automation that fits your lab



PIONEERING DIAGNOSTICS

# Automate Your Specialty Testing with VIDAS 3

- Sepsis and other specialty assays at your fingertips
- On-demand results and superior accuracy you can trust
- High flexibility and traceability for improved productivity

# VIDAS® 3

## Automation that fits your lab



### Touchscreen

- For guided instrument operation, reporting and troubleshooting

### LED Communication Bars

- Alerts user to status of instrument and individual sections

### 4 Sections of 3 Positions

- Up to 12 tests run at the same time across 4 independent sections for random access, greater flexibility and traceability

### Automated Pipettor

- Ease of use and improved workflow through on-board automated pipetting, calibration and dilution

### Sample Racks

- 3 sample racks offering 27 positions for patient samples, controls and diluents in a thermoregulated unit, with barcoded traceability

### Bidirectional Connectivity

- For order set entry, traceability, result transmission and remote support

### Consumables Storage and Disposal

- Storage for 96 disposable pipettor tips and 32 single-use dilution cups, along with on-board waste drawer for safe disposal after use





## VIDAS® 3 uses VIDAS reagents:

- **Specialty menu:** critical care and infectious diseases assays
- **On-demand testing:** 1 patient, 1 test, 1 result
- **Ready-to-use reagents**
- **One order reference per test:** kits contain all components needed for one test

The VIDAS 3 diagnostic platform enhances the clinical lab with a new generation of automation, traceability and productivity.

## Preparation

Our menu of **high-value, specialized assays** can be quickly loaded and analyzed on the VIDAS 3 system to meet your rapidly evolving needs as they arise.

Designed for **continuous loading**, the VIDAS 3 comes equipped with four compartments—each capable of running three assays, totaling **36 assays per hour** and holding up to 27 samples, diluents and calibrators in a temperature-controlled loading bay.

## Specialization

With an extensive menu of specialty diagnostic assays in a **single-sample, single-test format**, the VIDAS 3 enables you to run specialty esoteric assay autonomously in-house, saving time and money in the long run.

The **easy-to-learn, plug and play** assay kits allow lab technicians to quickly adopt the protocol into their routine. VIDAS 3 even offers ready-to-use VIDAS reagents, giving you everything you need to get started.

## Reliability

The VIDAS 3 smart read system **improves security and QC management** and guarantees traceability of runs with barcode scanning and **bidirectional LIS connectivity**. Combined with **24/7 access**, the VIDAS 3 gives you on-demand, accurate results you can trust.

With a facilitated lab workflow, **automated pipettor action, calibration and on-board dilution**, along with the **unparalleled reliability** and rich heritage of the VIDAS family, VIDAS 3 offers you the utmost peace of mind, ready to embrace a shifting healthcare landscape.

For more information on the VIDAS 3 including our full assay list, visit [www.biomerieux-usa.com/vidas3](http://www.biomerieux-usa.com/vidas3)