

THE WORLD OF VIDAS®  
TRUSTED. ROBUST.  
UNIVERSAL.

**Your lab is evolving.** We understand your needs, whether you are new to automation or upgrading to more advanced technology.

We recognize your need for **simple, built-to-last** solutions with excellent service to keep your lab up and running.

We designed VIDAS® KUBE™ **to offer you advanced technology** with reagents specifically designed for veterinary needs.

By partnering with bioMérieux, you benefit from **over 30 years** of worldwide experience in immunoassays.

VIDAS® OVER 30 YEARS  
OF INNOVATION



04/2025 / BMX 1165750 / 9322986 008/GB/A / This document is not legally binding. bioMérieux reserves the right to modify specifications without notice. / bioMérieux logo, PIONEERING DIAGNOSTICS, VIDAS, VIDAS KUBE, VILINK are used, pending and/or registered trademarks belonging to bioMérieux or one of its subsidiaries or one of its companies. Any other name or trademark is the property of its respective owner / Photos : Getty Images / bioMérieux S.A. RCS Lyon 673 620 399 / Printed in France / Héra / RCS Lyon 8 398 160 242. Do not litter



INSTRUMENT SPECIFICATIONS

|                                      |  |
|--------------------------------------|--|
| <b>VIDAS® KUBE™ References</b>       | 423905 (single unit)<br>47 08 186 (veterinary configuration)   |
| <b>Weight</b>                        | 23.7 kg (single unit)  |
| <b>Dimensions</b>                    | L. 43.4 cm x H 32.1 cm x P 48.5 cm (single unit)   |
| <b>Electrical</b>                    | 100/240 V, 56-60Hz   |
| <b>Power</b>                         | 150 W  |
| <b>Energy efficiency</b>             | Sleep mode and programmable morning start up   |
| <b>Operating temperature</b>         | 18-30°C  |
| <b>Operating humidity</b>            | 20-80% relative humidity non-condensing  |
| <b>Maximum altitude</b>              | 2700 m   |
| <b>Sound level in operating mode</b> | ≤ 70 dBa   |
| <b>Usability criteria</b>            | <ul style="list-style-type: none"><li>• Data purge: Automatic data archiving</li><li>• User safety during testing: Auto-locking of section doors</li><li>• Calibration setting: Modifiable</li></ul> |

PC specifications

An additional PC is required to configure 3 to 6 VIDAS® KUBE™ modules.

|                         |  |
|-------------------------|--|
| <b>OS Configuration</b> | Windows 10   |
| <b>CPU</b>              | Intel Core i3 – 8100<br>3.6 GHz                                |
| <b>Memory</b>           | 27000 tests results<br>Approx. 5 years of data                 |
| <b>RAM</b>              | DDR4 - 16GB<br>2400 MHz  |
| <b>Connectors</b>       | USB port x 6<br>Ethernet port x 1<br>Display port (screen) x 1 |

VIDAS® KUBE™  
Enter a World of Simplicity

ADVANCED DIAGNOSTIC SOLUTIONS  
FOR ANIMAL HEALTH





# VIDAS® KUBE™ ENTER A WORLD OF SIMPLICITY

## VETERINARY IMMUNOASSAY NEEDS ARE EVOLVING

You need to provide accurate results in a timely manner to ensure animal owner and veterinarian satisfaction remains high.


**Step into a world of expertise with VIDAS® KUBE™** a truly flexible, cost-effective automated platform that provides rapid results in complete confidence. Scaling with ease, VIDAS® KUBE™ is a simple and reliable addition to your evolving lab!



## VIDAS® CONCEPT INSIDE

### UNIQUE SINGLE TEST CONCEPT

- Ready-to-use reagents with long shelf-life
- Everything is included in the kit



## EASY IMPLEMENTATION

**VIDAS® KUBE™ is easy to implement in all circumstances:**

- Load-and-go system with ready-to-use reagents
- Run just a single sample or up to 36 samples per hour
- Stackable benchtop design - 1 - 6 modules
- Convenient benchtop design
- Integrated computer with intuitive software

## INCREASED CONFIDENCE IN RESULTS

**VIDAS® KUBE™ helps your lab achieve reliable and accurate results:**

- Published performance data on veterinary samples
- Alerts in case of test issues (sample detection, Strip & SPR compatibility, leakage/clogging issues and more)
- Validated long calibration intervals

## RECOGNIZED HIGH QUALITY SERVICE

**VIDAS® KUBE™ offers long-term peace of mind with:**

- Maintenance supported by a high-quality service team
- Evolutive and adaptable training journey
- Optimized productivity with remote maintenance via VILINK® connectivity

**VIDAS® KUBE™**  
adaptable automation  
will help your lab be  
more accurate and more  
innovative - all in one  
simplified platform



\*Please consult your local representative to find out more