

BIOMÉRIEUX

GENE-UP® ASPERGILLUS PRO

Quality product deserves quality testing.



Microbiology. Simple. Reliable.

PIONEERING DIAGNOSTICS

Aspergillus testing you can count on from an industry leader.

Simple, easy to use, and robust – the GENE-UP® Aspergillus PRO Screening and Upfront Speciation Kit comes from one of the most respected microbiological diagnostics companies in the world.

Uniquely designed for higher accuracy with minimal sample prep, our patented assay helps your lab streamline *Aspergillus* testing, increase overall efficiency, and maintain a higher degree of regulatory compliance.



STREAMLINED WORKFLOW

- Minimal sample prep eliminates the need for robotics
- No secondary enrichment step
- No immunomagnetic separation required



ENRICHMENT COMPATIBILITY ACROSS MATRICES

- 24 or 48 hour single-step enrichment protocols
- Compatible with cannabis and hemp matrices
- One protocol for all matrices



UPFRONT SPECIATION

- 5-channel assay
- Upfront speciation of *Aspergillus flavus*, *niger*, *fumigatus*, and *terreus*
- Superior specificity



EASILY INTERPRETED RESULTS

- No melt curve interpretation
- Complex algorithms eliminated
- Objective, repeatable results



HIGH/LOW THROUGHPUT

- Run 1 sample, or up to 96 at one time
- Full integration with most LIMS systems



OPEN PLATFORM

- Other AOAC approved assays also available, including *Salmonella*, STEC, *Listeria*, and more
- Robust R&D resources to meet the industry's evolving needs




UNPARALLELED SERVICE

- The bioMérieux team has decades of experience to assist you-from field based Cannabis Segment Specialists and Field Service Engineers to remote based Scientific Affairs and Technical Support.

Game-Changing Workflow:

Sample Prep: < 5 min hands-on time * 48-hour *Aspergillus* protocol



Order information:
Item IS1086 GENE-UP
Aspergillus PRO KIT
24 tests

Kit includes:
Assays & Reagents
PDB not included in Kit

