



BIOMÉRIEUX

GENE-UP® PRO STEC/SALMONELLA

Double your throughput with a single test.



Microbiology. Simple. Reliable.

PIONEERING DIAGNOSTICS

Detect STEC and *Salmonella* in Cannabis with a Single Test

Double your throughput with the GENE-UP® PRO STEC/*Salmonella* assay, designed to revolutionize Cannabis microbiology testing. This easy-to-use, scientifically robust solution offers simultaneous PCR detection of STEC and *Salmonella* with a single enrichment, delivering reliable results in just 24 hours.

Optimize your laboratory workflow and achieve your compliance goals with this efficient, single-test approach. Integrate the GENE-UP® PRO STEC/*Salmonella* assay to harmonize your operations, and further enhance your efficiency by running it concurrently with the GENE-UP® *Aspergillus* PRO assay. Experience the ease and reliability of STEC + *Salmonella* + Internal Control, all in one test.

BENEFITS OF A GENE-UP® PRO STEC / SALMONELLA WORKFLOW



ORDERING INFORMATION

ISO683	GENE-UP® PRO STEC / <i>Salmonella</i> Kit	24 Assays
414057	GENE-UP® Lysis Kit	192/Kit
419194	GENE-UP® Biotix Tips	960 tips
ISEV013	<i>Salmonella</i> gDNA Positive Control	8 rxns.
ISEV015	<i>E.coli</i> O158:H7 gDNA Positive Control	8 rxns.

*additional formats available, ask for details

Learn more about our full suite of microbiology solutions for Cannabis & Hemp at [biomerieux.com/us/cannabis](https://www.biomerieux.com/us/cannabis)

HOW IT WORKS

