



IDENTIFICATION & CONFIRMATION

List of validated methods

Reference	Identification/ confirmation Method	Principle	ISO	AOAC RI	AOAC OMA	Health Canada	India (FSSAI), Rapid Analytical Food Testing (RAFT)		
				Certificate #	Certificate #	Certificate #		Expiration	
references according to microorganisms	API®/ID32 range	Phenotypic	Included in ISO 7218 standard for confirmation of ISO standard reference methods* Included in ISO 16140-2 validations for confirmation of all bioMérieux pathogenic detection methods			MFHPB-20 & 30 (for <i>Salmonella</i> and <i>Listeria</i> methods confirmation)	Direct approval	09.01.2023	
20100 - 20160	API® 20E				OMA 978.24 (<i>Salmonella</i>)				
411071 -	VITEK MS and VITEK MS PRIME	Maldi-toff	Included in ISO 16140-2 validations** for confirmation of bioMérieux alternative methods	PTM 011501					
21341/21342	VITEK® 2C GN VITEK® 2C GP	Phenotypic				MFHPB-20 & 30 (for <i>Salmonella</i> and <i>Listeria</i> methods confirmation)	Direct approval	09.01.2023	
21345	VITEK®2 BCL			PTM 080801					
21341	VITEK®2 GN					OMA 2011.17			
21342	VITEK®2 GP					OMA 2012.02			

*e.g., ISO 6579-1: Salmonella - ISO 11290-1&2: Listeria
ISO 21872-1: Vibrio - ISO 10272-1&2: Cronobacter - ISO 10273: Yersinia
ISO 10272-1&2: Campylobacter

**For: VIDAS® SPT - VIDAS® Easy *Salmonella* - VIDAS® *Campylobacter*
CampyFood® agar detection - CampyFood® agar enumeration - ESIA™
One Day

bioMérieux S.A. • 69280 Marcy l'Étoile • France • Tel: +33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux-industry.com