



# 3P<sup>®</sup> RELIABILITY MEETS GS1 DATA INTEGRITY.




Keeping your environmental monitoring data accurate, organized, and accessible not only keeps you in compliance—it increases patient safety and allows faster, more effective resolutions to problems.

## **GS1 BARCODES ENSURE ISO-COMPLIANT DATA INTEGRITY AND TRACEABILITY IN YOUR ENVIRONMENTAL MONITORING DATA.**

Thanks to its international recognition and compliance with ISO standards, GS1 two-dimensional (2D) DataMatrix barcodes allow easy, rapid identification of 3P<sup>®</sup> SMART culture media with a unique identification number for every sample.

- GS1 barcodes ensure compliance with ISO standards
- Well-defined barcode structure for easy data traceability
- Unique permanent ID number for all samples and culture media
- Easily identify and track both samples and media in real-time
- Compatible with commonly used barcode readers and software

## A WELL-DEFINED DATA STRUCTURE MAKES GS1 BARCODES EASY TO IMPLEMENT AND USE.



### bioMérieux well-defined data structure using GS1 standard Application Identifiers

NAME	GS1 APPLICATION IDENTIFIERS	LENGTH
GTIN Global Trade Item Number	01	14
Expiration Date YYMMDD	17	6
Batch or lot number	10	10
Serial number	21	6



**Unique ID  
Number**



**Approved for  
Healthcare Industries**



**DataMatrix  
Legibility**



**Defined DataMatrix  
Structure**



**Compliant with  
International Standards**

To learn more about 3P SMART PLATES, talk to your bioMérieux representative, or visit:

**[biomerieux-industry.com](http://biomerieux-industry.com)**