



MEDIA

List of validated methods

Reference	Target	Method	ISO 16140-2 (2016)			AOAC RI	Health Canada	India (FSSAI), Rapid Analytical Food Testing (RAFT)	
			Reference method	Certificate #	Expiration	Certificate #	Certificate #	Expiration	
AEB520020 / AEB620027 / AEB184135 / AEB184135/10 /AEB150022	β-glucuronidase-positive <i>Escherichia coli</i>	Rebecca™	ISO 16649-2	AFNOR AES 10/06-01/08	2028 01			Direct approval	renewal in progress
AEB520020 / AEB620027 / AEB184135 / AEB184135/10 /AEB150022	<i>Enterobacteriaceae</i>	Rebecca™	ISO 21528-2	AFNOR AES 10/07-01/08	2028 01				
42017	Coliforms (37°C)	chromID® coli	ISO 4832	AFNOR BIO12/20-12/06 (coliforms 37°C)	2026 12				
42017	β-glucuronidase-positive <i>Escherichia coli</i> (37°C)	chromID® coli	ISO 16649-2	AFNOR BIO12/19-12/06	2026 12				
42017	β-glucuronidase-positive <i>Escherichia coli</i> (44°C)	chromID® coli	ISO 16649-2	AFNOR BIO 12/5-01/99	2027 01				
AEB520069 / AEB620067	Detection <i>Salmonella</i>	SMS®	ISO 6579-1	AFNOR AES 10/04-05/04	2028 05				
418246 / 418247	Detection <i>Salmonella</i>	SALMA® One Day	ISO 6579-1	AFNOR BIO 12/41-03/17	2025 03				renewal in progress
AEB520080 / AEB520079 / AEB120082 / AEB150072 / AEB684420/ 423151	Detection <i>Listeria monocytogenes</i> and <i>Listeria</i> spp.	ALOA® One Day	ISO 11290-1	AFNOR AES 10/03-09/00	2028 09		MFHPB-30 (for <i>Listeria</i> methods confirmation)	Direct approval	
AEB520080 / AEB520079 / AEB120082 / AEB150072 / AEB684420/ 423151	Enumeration <i>Listeria monocytogenes</i> and <i>Listeria</i> spp.	ALOA® Count	ISO 11290-2	AFNOR AES 10/05-09/06	2026 09		MFHPB-30 (for <i>Listeria</i> methods confirmation)		
43471	Enumeration <i>Campylobacter</i> spp.	CampyFood® agar	ISO 10272-2	MICROVAL 2009LR28*	2026 03				
43471	Detection <i>Campylobacter</i> spp.	CampyFood® agar	ISO 10272-1	AFNOR BIO 12/30-05/10*	2026 05			renewal in progress	
AEB520100 / AEB120102 / 423148 / 423849 / 423868 / 423871	Enumeration <i>Bacillus cereus</i>	BACARA® method (BACARA® and BACARA® 2)	ISO 7932	AFNOR AES 10/10-07/10	2026 07				
AEB520010 / AEB150002	Detection <i>Cronobacter</i> spp.	ESIATM One day	ISO 22964	AFNOR BIO 12/37-11/14 *	2026 11				

bioMérieux media may also be standardized in ISO reference methods, FDA BAM methods etc... (e.g., BAM Media M10a for ALOA®, BAM chapter 14 for BACARA®)

*The scope includes the use of the VITEK® MS and VITEK® MS PRIME instruments for confirmation

bioMérieux S.A. • 69280 Marcy l'Étoile • France • Tel: +33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com