

BIOMÉRIEUX

GENE-UP[®] ENVIROPRO[™]

Simultaneously Detect
Salmonella and *Listeria* species



Your Trusted Partner in Augmented Diagnostics

PIONEERING DIAGNOSTICS

Introducing GENE-UP® ENVIROPRO™

Simultaneous detection of *Salmonella* & *Listeria* Species from environmental swabs.


A multiplexed configuration that simultaneously detects *Salmonella* species and *Listeria* species from environmental samples.

By combining these assays into one tube, it provides both cost-savings and a streamlined workflow to the end-user. The best-in-class workflow eliminates cumbersome DNA purification steps and reduces hands-on time.


Optimized Environmental sampling: One swab, one enrichment, one sample prep, two results!

Ultimately, GENE-UP ENVIROPRO is a game changer, improving overall sampling efficiencies and providing rapid data for quality decisions.

Benefits of the GENE-UP® ENVIROPRO™ *Salmonella* & *Listeria* Environmental Workflow



How it works



Ordering Information

IS1071	GENE-UP ENVIROPRO Kit	24 tests/kit
ISO374	Dehydrated Universal Enrichment Media	2kg
414057	GENE-UP Lysis Kit	192 tests/kit
419194	GENE-UP Tips	960 tips
IS0935	VSD Sample Tubes	24 tubes

Contact your local bioMérieux representative for availability