

BIOMÉRIEUX

TEMPO[®] FOOD QUALITY AUTOMATION

Because You Count!



Your Trusted Partner in Augmented Diagnostics

PIONEERING DIAGNOSTICS




Automated Quality Indicator Solution


ONLY **3** STEPS

as simple as 1, 2, 3!

1 Automated sample preparation



2 Incubation
in your standard incubators



3 Automated reading and recording



TEMPO® is a fully automated enumeration system that tests quality indicators in food products and environmental samples within your laboratory.



What does it mean for you?



People Empowerment

TEMPO® empowers your team with easy learning, prevents errors, and reduces extra hours and weekend work.

- Automation
- Reading time: <5mins - 20 cards
- No confirmation



Confidence In Quality And Results

TEMPO® offers an easier, faster process with no serial dilution and no cross-contamination.

- Simple and intuitive operation: only 3 steps
- Officially validated: AFNOR, AOAC, MICROVAL & Health Canada
- Sealed Cards



Traceability

TEMPO® ensures traceability, connectivity, and audit readiness.

- LIMS Connectivity
- Integrated Control plan & barcode
- VILINK® Remote Access



TEMPO[®] parameters

	TEMPO [®] parameters	NEXT DAY SOLUTION	ISO 16140-2	AOAC	MLG*
411113	AC Aerobic mesophilic total flora**	✓	✓	✓	✓
80002	STA Staphylococcus (coag+)	✓	✓	✓	
80003	EB Enterobacteriaceae	✓	✓	✓	✓
80004	EC Escherishia coli	✓	✓	✓	✓
80006	TC Coliforms	✓	✓		
80044	CC Coliforms	✓		✓	✓
80106	BC Bacillus cereus	✓	✓	✓	
80001	YM Yeast & Mold	72-76 h		✓	
80071	LAB Lactic acid bacteria	40-48 h			
421509	CAM Campylobacter	48 h	✓	✓	

*USDA Microbiology Laboratory Guidebook

** TEMPO[®] Now allowing AC thermophilic micro-organisms enumeration - 55°C



Do you need an automated quality indicator solution?

Contact your local bioMérieux representative about TEMPO[®] and learn more on our website.